
M1450- 2023-04-24 - Rev. 3.0

ENGLISH	3 - 5
DEUTSCH	6 - 8
NEDERLANDS	9 - 11
DANSK	12 - 14
NORSK	15 - 17
SVENSKA	18 - 20
SUOMI	21 - 23
ESPAÑOL	24 - 26
FRANÇAIS	27 - 29
ITALIANO	30 - 32
PORTUGUÊS	33 - 35
PORTUGUÊS DO BRASIL	36 - 38
РОССИЯ	39 - 41
POLSKI	42 - 44
ČESKY	45 - 47
SLOVENSKY	48 - 50
中國	51 - 53
日本	54 - 56
ΕΛΛΗΝΙΚΗ	57 - 59
ÍSLENSKA	60 - 62


Roadability

Driving UP kerbs (max. 100 mm)

Forwards


This technique is recommended for the experienced users.

- Make sure the anti-tip device is deactivated
- Drive forward towards the kerb
- Balance on the rear wheels and lift the front wheels from the ground. Bend forward and pull hard on the hand rims.

Forwards, with assistance


- Make sure the anti-tip device is deactivated
- Drive forward towards the kerb
- Tilt the chair backwards to lift the front wheels above the kerb.
- Lift on the handles to drive the rear wheels above the kerb.

Backwards

This technique is only to be used at low edges and is depended on the height of the footplate.

- Make sure the anti-tip device is deactivated
- Drive backwards towards the kerb
- Bend forward and pull hard on the hand rims.

Backwards, with assistance


- Make sure the anti-tip device is deactivated
- Drive backwards towards the kerb.
- Tilt the chair backwards to lift the front wheels.
- Lift on the handles to pull the chair up. Pull the chair away from the kerb and put down the front wheels.

Driving DOWN kerbs (max. 10 cm)

Forwards


This technique is recommended for the experienced users.

- Make sure the anti-tip device is deactivated
- Drive forward towards the kerb
- Pull hard on the hand rims and drive forward. All wheels will touch the ground at the same time.

Forwards, with assistance


- Make sure the anti-tip device is deactivated
- Drive forward towards the kerb.
- Tilt the chair backwards to lift the front wheels.
- Drive carefully down the kerb and put down the front wheels.

Backwards


This technique must never be used if the height difference is more than 10 cm and only by experienced users.

- Make sure the anti-tip device is deactivated
- Drive backwards towards the kerb
- Bend forward and drive carefully down the kerb.
- High tilting risk at this technique

Backwards, with assistance

- Make sure the anti-tip device is deactivated
- Drive backwards towards the kerb.
- Drive carefully down the kerb and pull the chair away from the kerb.
- Put down the front wheels.

Driving UP stairs

- Always ask for help
- Never enter an escalator even with help from others

Backwards, with assistance

- Make sure the anti-tip device is deactivated and the handles is tightened securely
- Drive backwards toward the first step
- Tilt the chair and pull up the chair one step at a time. Make sure to keep balance on the rear wheels
- Pull the chair away from the last step and put down the front wheels

Driving DOWN stairs

- Always ask for help
- Never enter an escalator even with help from others

Forward, with assistance

- Make sure the anti-tip device is deactivated
- Drive towards the first step and tilt the chair
- Drive down one step at a time. Make sure to keep balance on the rear wheels
- Put down the front wheels after the last step
- If two helpers are present, one can assist by lifting at the front
- The helpers need to make sure to lift the correct way - with the strength from the legs and not from the back

Driving DOWN/UP a hill


Follow the below to avoid tilting when driving down/up a hill.


Driving down:

Sit upright in the chair. Correct the speed with the hand rims and not with the brakes.

Driving up:

- Bend forward to correct the point of gravity.
- Avoid to turn around in the middle of the hill.
- Always sit as upright as possible
- Better to ask for help than to take a risk on your own.

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R82 A/S