Transportation in motor vehicles

ISO 7176-19:2008
ISO 16840-4:2009

2021-11-rev.5.0

www.etac.com

ENGLISH	4 - 7
DEUTSCH	8 - 11
NEDERLANDS	12 - 15
DANSK	16 - 19
NORSK	20 - 23
SVENSKA	24 - 27
SUOMI	28 - 31
ESPAÑOL	32 - 35
FRANÇAIS	36 - 39
ITALIANO	40 - 43
PORTUGUÊS	44 - 47
PORTUGUÊS DO BRASIL	48 - 51
РОССИЯ	52 - 55
POLSKI	56 - 59
ČESKY	60 - 63
中國	64 - 67
日本	68 - 71
ΕΛΛΗΝΙΚΗ	72 - 75

ISO 7176-19:2008

- Cricket
- Stingray
- High low:x outdoor
- High low:xo
- Combi Frame:x
- Kudu
- Multi Frame:x
- Multi Frame

EN

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ENGLISH

Transportation in motor vehicles

The instructions on how to prepare the seat and chair, must be carried through before transportation.

The user should transfer to the vehicle seat and use the vehicle-installed restraint system whenever possible, and the unoccupied wheelchair should be stored in a cargo area or secured in the vehicle during travel.

The wheelchair must be placed in a forward-facing position, when used as a seat in a motor vehicle.

The front wheels need to be turned under the frame so as to be in line with the frame before fixation.

Use an approved 4-point strap-type tie-down system according to ISO 10542-2. to secure the product. Use either a hook or a strap in the fitting, see figure 1.

The location of all wheelchair securement points are marked with a label, see figure 2.

The pelvic restraint should be worn low across the front of the pelvis, so that the angle of the pelvic restraint is within the preferred zone of 30° to 75° to the horizontal, similar to that shown in figure 3, a steeper (greater) angle within the preferred zone is desirable.

The back must be placed in vertical position and the seat plane must be horizontal.


Figure 1.

Figure 2

On yellow background


Figure 3.

ENGLISH

Use an approved 3-point belt according to ISO 10542-1 in the product.

The belt restraints should be adjusted as tightly as possible, consistent with user comfort. Further the belt webbing should not be twisted when in use.

Belt restraints should not be held away from the body by wheelchair components or parts, such as the wheelchair armrests or wheels. See figure 4.

Shoulder belt restraints should fit over the shoulders, similar to the illustration provided in figure 5.

An updated table covering configurations for seats and frames, max. load/user weight in transportation can be found in the overview “frame_and_seat_combination.pdf” at

www.etac.com/support/support-center/support-documents/

Check the User Guide and or Mounting Guides for the products mentioned on page 3 for information about:

- Minimum weight limit
- Securement points
- Max load for seat or frame, also stated on each product

Figure 4: Belt restraints should not be held away from the body by wheelchair components or parts, such as the wheelchair armrests or wheels.


Figure 5

03 - Warnings!

The product is approved according to ISO 7176-19:2008 furthermore the product complies to ISO 16840-4:2009 if the product is a seat..

The wheelchair has been dynamically tested in a forward facing orientation with the ATD restrained by both pelvic and shoulder belts (e.g. a shoulder belt as part of a three-point belt restraint).

Pelvic and shoulder belt restraints should be used to reduce the possibility of head and chest impacts with vehicle components.

Loose parts and trays must be removed from the wheelchair and secured separately in the vehicle. Other auxiliary wheelchair equipment should be either secured to the wheelchair or removed from the wheelchair and secured in the vehicle during travel, so that it does not break free and cause injury to vehicle occupants in the event of a collision.

The product should be inspected by the dealer before reuse following involvement in any type of vehicle collision.

Alterations or substitutions should not be made to the wheelchair securement points or to structural and

frame parts or components without consulting the wheelchair manufacturer. The approval is not valid on 'custom made' chairs.

Care should be taken when applying the occupant restraint to position the seatbelt buckle so that the release button will not be contacted by wheelchair components during a crash.

Postural supports should not be relied on for occupant restraint in a moving vehicle, unless they are labelled as being in accordance with the requirements specified in ISO 7176-19:2008.

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© 2017 R82 A/S All rights reserved.
The R82 logo and the icons are registered trademarks of R82 A/S
R82 A/S