

Molift Air 200 / 350

molift[®]
by Etac

User Manual - Bruksanvisning - Brukermanual - Brugervejledning - Käyttö-
pas - Bedienungsanleitung - Gebruikershandleiding - Manuel d'utilisation
- Manuale per l'utente - Manual de usuario - Příručka pro uživatele

BM16199 Rev. L 2024-02-14

Contents

User manual - English	- 1
Gebrauchsanweisung - Deutsch	- 13
Brukermanual - Norsk	- 27
Bruksanvisning - Svenska	- 39
Käyttöohje - Suomi	- 51
Manual utente - Italiano	- 63
Manual d`utilisation - Francais	- 75
Manual de usuario - Español	- 87
Handleiding - Nederlands	- 99
Brugsvejledning - Dansk	- 111
Příručka pro uživatele - Čeština	- 123

Important

This User Manual contains important safety instructions and information regarding the use of the sling and accessories. In this manual the user is the person being lifted. The carer is the person operating the hoist.


This symbol indicates important information related to safety. Follow these instructions carefully.

Read User Manual before use!
It is important to fully understand the content of the user manual before attempting to use the equipment.

Visit www.etac.com for download of documentation to ensure you have the latest version.

Molift Air 200/350

Hoist components:

1. Trolley
2. Hoist
3. Side covers
4. Emergency stop/ Emergency lowering
5. Lifting belt
6. Quick release pin
7. Hand control
8. Sling bar

About Molift Air

Molift Air is a strong and smooth ceiling lift that enables users and residences to be transferred in a comfortable and safe way. Suitable for sitting and lateral transfers as well as standing and gait training situations together with the comprehensive accessory program from Molift.

Molift Air is a very light motor which in combination with the quick release coupling system gives an easy handling when mounting and servicing.

General

CE declaration

The product, and its related accessories, described in this instruction for use (IFU), is in compliance with the regulation (EU) 2017/745 of 5. April 2017 – as a medical device, risk class I.

The product has been tested and approved by a third party according to standards IEC 60601-1, IEC 60601-1-2 and EN/ISO 10535:2021.

Any serious incident that occurs in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/ or patient is established.

Conditions for Use

Lift and transfer of a user will always pose a certain risk and only trained personnel are allowed to use the equipment and accessories covered by this user manual. The rail system must be installed by certified personnel in accordance with applicable installation instructions. The hoist is not intended to be operated by the user being lifted.

Modifications and use of components made by other manufacturers.

We recommend only using Molift components and spare parts. Declaration of conformity is not valid and Etac is not responsible for warranty if any modifications are made to the product. Etac shall not be liable for faults or accidents that can occur when using components made by other manufacturers. It can be unsafe to use accessories, detachable parts and materials that is not listed in this user manual.

Only certified personnel are allowed to open hoist or accessories to perform service or repAir. Risk of injury from rotating parts and electric shock.

Warranty notice

Two-year warranty against defects in workmanship and materials of our products. One-year warranty for batteries. Please refer to www.etac.com for terms and conditions.

Definitions

The person described as “the user” in this Instruction for Use is the person who is sitting in the product. “The carer” is the person who manoeuvres the product.

Intended purpose

To be used for lifting and transferring users, who require support to be lifted or transferred within a rail system; e.g. to/from a bed, the floor, a wheelchAir, a shower commode chAir, a toilet or similar. Also for gait training in combination with ambulating vest.

Intended target group

Used for users who require support to be lifted and transferred; e.g. persons with limited ability to move and who require body support during transfers. Not an exhaustive list.

Intended environment

For indoor use within a rail system in acute care, long term care, home care.

Indications

The device can be used for children, adults and elderly people with a need for support to be lifted and transferred, or when gait training, within a railsystem.

The target group for the device is based on individual health and mobility function, and not a specific diagnose or age.


Contra indications

There are no known contra indications regarding the hoist, but make sure the sling fits the user.

Product identification


Product label

The Product labels barcode contain article number, serial number and production date:

Warning labels and Symbols

Symbols used on the product, explained in more detail:

-  Class II double insulated
-  CE marked
-  Refer to user manual
-  Indoor Use only
-  Max user weight
255 kg / 560 lbs
-  Do not dispose in general waste
-  Emergency stop/
Emergency lowering
-  Service Light
-  Battery Light
-  Hex tool symbol for
Manual Emergency lowering
-  Certification Body mark
c N us
0601-1
-  Humidity limitation
-  Atmospheric pressure limitation
-  Temperature limitation
-  Manufacturer
-  Date of manufacture
YYYY-MM-DD
-  REF Catalogue number
-  SN Serial number
-  MD Medical Device

EN
DE
NO
SV
FI
IT
FR
ES
NL
DA
CZ

Technical data

Safe Working Load (SWL)

Molift Air 200: SWL = 205 kg (450 lbs)
 Optional SWL range (setting) 160-205 kg
 Molift Air 350: SWL = 350 kg (771 lbs)
 Optional SWL range (setting) 160-350 kg

Weight of unit:

Chassis incl. battery, Excl. sling bar:
 8,56 kg (18,87 lbs)
 4-point Medium Sling bar: 2,2 kg/4.8 lbs

Battery:

Art: 2920204 25.6 V LiFePO4 2.4 Ah
 Art: 2920007 12 V SLA 2.9 Ah

Battery charger:

Charger for handcontrol, Mascot - type 3743
 Art. no.: 2540100
 100-240 V AC, 50-60 Hz max 0.5A
 Output 0.9A
 IRC Power Supply, Mascot - type 3721
 Art. no.: 2610119
 Input: 100-240 VAC, 50-60 Hz max 1.0A
 Output: 24VDC 1.75A/42W
 IRC Charger, Mascot - type 3844
 Art. no.: 2640200
 In: 20-30VDC 1.3A
 Output 0.5A

Lifting speed

60mm/second (2,36 inches/second)
 with 75kg (165,35 lbs) load

Propulsion motor speed

0,2m/second (7,87 inches/second)

Protection class

Hand control IP24
 Lift motor IP24

Sound pressure

52 dBA

Sound power (L_{WA})

63 dB

Operating forces button

Buttons on handset: 3.4 N

Material:

Aluminum, plastic, steel

Expected Lifetime:

The hoist has an expected lifetime of 30 000 cycles SWL or 10 years.

Dimensions:

L x W x H (Length, Width, Height without trolley)

36 x 19 x 19,5 cm

(14,2 x 7,5 x 7,7 inches)

A Air 200 Hoisting Range:


300 cm (118 Inches)

Air 350 Hoisting Range:

260 cm (102 Inches)

B (Minimum distance from ceiling to sling bar coupling point) 32,7 cm (12.7/8 Inches)

For extended information about technical data, how to preform service and repAir, see technical manual.


Installation

The hoist is marked with Safe Working Load (SWL), this should not exceed rail systems max load capacity.

 The rail system can only be installed by certified personnel in accordance with applicable installation instructions.

 Do not start using hoist before completing control according to checklist after installation.

Trolley installation


Mounting trolley in rail system. Remove End stop and insert Trolley into rail. Mount end stop, and make sure it is securely fastened. See BM3401 for installation details. Etac offers trolley for mounting Molift's lifts in rail systems from other manufacturers. Contact your local Etac representative for more information on trollys compatible with other manufacturers' rail systems.

 Make sure end stops are mounted in all ends of the Rail system before using the hoist!


Hoist installation

 Be carefull during elevated work situations. We recomend two or more installers, to prevent loosing balance and injury during installation!

Mounting of hoist on trolley. No tools required. The hoist has two connection points.

Connect one point first. Push button all the way in, and place hoist in connection point on trolley. Release button and make sure it is clearly showing green before connecting the next point.

Make sure both connection point buttons completely return and clearly show green after installation. Emergency lowering cord and tube can be adjusted to the correct length (height) by cutting tube and cord.

Installation

Use this checklist to verify that the hoist is properly installed and can operate correctly and safely before use.

- End stops on rail are installed after mounting of trolley as described in BM3401
- Make sure that the hoist is properly fastened to the rail and that the lift does not have any loose parts.
- Perform one lift with load (60-80 kg)
- Make sure battery is fully charged.

How to use Molift Air

General Safety Precautions

Only use accessories and slings that are adjusted to fit the user, type of disability, size, weight and type of transfer.

Working pause ratio/ Duty Cycle.

Molift Air 200 / 350 should not be run constantly for more than 2 minutes (with maximum load), and rest for minimum 18 minutes. Duty Cycle 10%. (Intermittence according to standard ISO-EN 10535)

 Do not position hoist in a way that makes it difficult to disconnect mains cable. To isolate the system, pull out the main cable from the electrical wall socket.

 If maximum load (SWL) differs between hoist sling bar and body support unit, then the lowest maximum load shall always be used

 Molift hoists shall only be used to lift users. Never use the hoist to lift or move objects of any kind.

Before use / Daily check



Inspection to be performed by carer daily or before use:

- Inspect hoist and rail system has no visible damage, defects or deformations
- Make sure sling bar connection and all detachable parts are properly connected and secured
- Test emergency stop button and emergency lowering
- Make sure battery indicator or service light is not illuminated
- Make sure lifting strap does not have visible damage or frays.
- Test run hoist's operation functions and make sure hoist does not make any abnormal sounds

If there are any faults or defects, the lift needs to be taken out of operation and marked "out of order". No parts on the Molift Air should be serviced or maintained while in use with a user.

Hand control


The Hand control has 2 buttons for lifting and lowering, or 4 buttons if the hoist is equipped with propulsion. The Hand control has an indicator light that will illuminate when battery level is low and the hoist requires charging

Push with finger or use tube on cord to push button back in to reset emergency stop.

Electrical Emergency Lowering

Pull and hold to start lowering. Hold until user is lowered and can be released from sling bar.

Do not pull Hand control to move the hoist along the rail

The lift will not work if hand control is connected to a charger. Disconnect hand control from charger before operating the lift.


Only for emergency use when hand control does not work!

Manual Emergency lowering

Manual emergency lowering is only to be performed by service personnel!


Allen key or bit for power drill can be found inside the lifter. Remove cover to locate.

Push down and press cover outwards with a flat screwdriver or coin in the two slots to open cover.

Locate sticker over the hole for manual emergency lowering. Punch hole or remove sticker and insert tool. Rotate counterclockwise to lower lifting belt. Using a power drill is recommended.

Mounting sling bar


1. Align sling bar in connection point.
2. Push button on locking pin and insert all the way through.
3. Make sure locking pin is properly fastened
4. Push and hold down the button on the locking pin, and pull to remove locking pin.


Emergency stop / emergency lowering

Emergency Stop

Activation and reset of emergency stop:

Pull to activate emergency stop. The button will come out, and hoist will stop.

Service is NOT required if:

- 1) the Air manual emergency down is used for a known reason, and
- 2) the belt is undamaged, and
- 3) the unit is operating normally after the manual emergency lowering procedure.

The new brake override will contribute to this by making emergency manual lowering less stressful to lift components.

If the unit have had a malfunction with unknown reason, and manual emergency lowering system has been used, Service is required.

Electronics

Transport and Storage

For long time storage it is recommended that the emergency stop button is activated (pulled out). The hoist can be stored and transported under temperatures between -25 - 70 °C.

Operating

The hoist is designed for use at standard room temperatures (+5 to +40°C).

Air Pressure: 70 - 106 kPa

Relative Humidity: 15 - 90 %

Following storage or transport at other temperatures leave the hoist in a room with a suitable temperature until it reaches a safe operating temperature.

Medical electrical equipment requires special precautions regarding electromagnetic compatibility (EMC). Portable or mobile radio communication equipment may affect the medical electrical equipment, and should be kept minimum 25 cm (10 inches) from the hoist's electronics.

Batteries

Molift Air 200/350 is available with 2 x 12V SLA battery or a 25.6V LiFePO4 battery.

The SLA battery has a life expectancy of approx. 500 charge cycles.

The LiFePO4 battery has a life expectancy of approx. 2000 charge cycles.

Recharge battery upon receipt, and then every six months, to avoid battery failure.

Batteries must be disposed of as special waste according to local rules and regulations. Do not dispose in general waste. Replacement of lithium batteries or fuel cells by inadequately trained personnel could result in a hazard

Battery and service indicator on hoist

	Wrench symbol / Service indicator
	Battery symbol / Battery indicator

The electrical system has a power-save function which will turn off the electrical system five minutes after last operation. All indicators will turn off. The system is activated when pushing one of the operating buttons on the hand control.

Battery indicator on hoist

When the multi level battery indicator is yellow the hoist will have sufficient power available for at least one full lifting cycle with max load. When battery is flashing yellow it is only possible to lower the sling bar.

Battery indicator	Battery level	
	Green	80% or more
	Green	80 - 60%
	Green	60 - 40%
	Green	40 - 20%
	Yellow	20% - 0%
	Yellow flashing	0%

Service indicator / Periodic inspection

The hoist's electronics record the weight lifted and number of lifts. After a certain period of operation a signal is given to indicate that service is required.

Service indicator	Mode
No light	No service needed
Yellow	Periodic inspection required
Green flashing	Order service
Red	Perform service

Service scope

When the hoist has performed a number of lifts, (5,000 lifts with the maximum allowed weight (SWL)), the service indicator will light up red.

A service must then be performed by authorized personnel, including the inspection and replacement of worn parts.

Contact Etac or your local dealer or visit etac.com to find an approved service partner.

Charging

Molift Air can be charged through the hand control with a battery charger, or through the rail with In Rail Charging.

Indicator on charger and power supply indicator


<i>Charger indicator</i>	<i>Mode</i>
Green flash	Not connected to hoist
Yellow	Charging
Yellow flash	Top of charge
Green	Fully charged

<i>Power supply indicator</i>	<i>Mode</i>
Blue	Connected to power outlet

Hand control charging


The battery charger is plugged into a power outlet. Wall-mounted docking for hand control is available. Make shure not to damage the cable. The charger can be connected to the power outlet at all times. Use only Molift chargers. The Hand control has an indicator light that will illuminate when battery level is low and the hoist requires charging.

Charger for Molift Air hand control: Art. no.: 2540100


IRC - In Rail Charging

The hoist is charged through a built in charger in lift motor for continuous charging through the entire rail system length. Molift Air IRC has a built in charger, connected to the battery. It must be installed with IRC Trolley in a rail system with a conductive tape. The IRC trolley is connected to the battery charger inside the Molift Air IRC. The IRC trolley is always in contact with the the conductive tape inside the rail. The IRC power supply is connected to the conductive tape through IRC connection at the end of the rail, and supplies power to the battery charger inside Molift Air IRC.

IRC fits both Single rail and Traverse systems.

IRC Power Supply

Art. no.: 2510119

Charger 3844 DC/DC (Part of Molift Air IRC)

Art. no.: 2640101

Transfer

Plan the lifting operation in advance to ensure that it is as safe and smooth as possible. Remember to work ergonomically. Assess the risks and take notes. The carer is responsible for the safety of the user.

Molift Air can be used with both 2-point and 4-point sling bar for different types of transfers; sitting, recumbent or ambulating.

Using slings

It is important that the sling has been tested with the individual user and for the intended lifting situation.

 Read User Manual for the sling prior to use.


 Do not to use damaged or badly worn slings.

Slings made by other manufacturers


We recommend using Etac or Molift slings. The hoist can be combined with slings from other manufacturers with similar suspension loops. Visit www.etac.com for compatibility statements between lifts and slings from other manufacturers. Please note that the natural or legal person using CE-marked compatible devices from different suppliers is responsible for the activity. Etac cannot therefore be held responsible for any incident.

Lifting and lowering


When moving the user, stand to the side of the user you are lifting. Make sure that arms and legs do not obstruct the seat, bed, etc. Keep eye contact with the user to help them feel safe.

The 4-point sling bar must always be positioned across the user, to prevent the user from unintentionally sliding out of the sling.

The sling bars hook design prevents sling to be inadvertently detached when mounted correct.


1. Check that the sling is correctly fitted around the user and that the strap loops are correctly fitted to the sling bar hooks.
2. Stretch the sling straps without lifting the user. Ensure that all four loops of the sling are securely fastened to avoid the user slipping or falling.
3. Lift user, and perform transfer.

 Do not use the emergency stop to turn OFF the hoist, except in an emergency. The hoist will switch OFF (stand by) automatically after 5 minutes of inactivity.

Angle sensor

The hoist has a directional safety feature that prevents lifting if angle on lifting strap is too steep.

If the hoist doesn't respond to Hand control "up" button, it might be because lifting strap angle is too steep, or lifting belt is twisted and sensor is activated.

Accessories

Charger cords and plugs.

Plug adapter EU, art.no 2510365

Plug adapter UK, art.no 2510366

Plug adapter US, art.no 2510367

Plug adapter AUS, art.no 2510368

Cord EU, art.no 1140105

Cord UK, art.no 1140106

Cord US, art.no 1140104


Cord AUS, art.no 1140107

Wall bracket for charger, art. no. 2610123

Sling bars

Molift Air can be used with 2 or 4 point sling bar.

Aluminum sling bar

2-point Small, aluminum, 350 mm,


Art. no.: 2530115

2-point Medium, aluminum, 450 mm,

Art. no.: 2530120

2-point Large, aluminum, 550 mm,

Art. no.: 2530125

4-point Small, aluminum, 350 mm,

Art. no.: 2530210

4-point Medium, aluminum, 450 mm,


Art. no.: 2530211

4-point Large, aluminum, 550 mm,

Art. no.: 2530212

4-point X-Large, aluminum, 650 mm,

Art. no.: 2530213

8-point sling bar for fabric stretcher, steel

Art. no.: 2140003

Scale

Molift Scale Set (without Sling bar)

Art. no.: 1840000

Can be combined with all Molift Air compatible sling bars. Scale is approved Class III. The scale can also calculate BMI (Body Mass Index) with a single touch. Read the manual that comes with the weight scale prior to use.

Slings

Etac supplies a wide selection of slings for different types of transfers. The Molift RgoSling sling series is developed to be combined with a 4-point sling bar with loop connections. The slings also works with a 2-point sling bar. See the combination list in the slings user manual for the correct sling and sling bar combination. The Molift RgoSling sling series is available in sizes XXS – XXL, in polyester and polyester mesh. Complete list of slings and sling bars can be found on www.etac.com

Molift Rgosling Mediumback Padded (XS-XXL)

Molift Rgosling Highback Padded (XS-XXL)

Molift Rgosling Mediumback Net (XXS-XXL)

Molift Rgosling Highback Net (XXS-XXL)

Molift Rgosling Toilet Lowback (XXS-XXL)

Molift Rgosling Toilet Highback (XXS-XXL)

Molift Rgosling Ampu Mediumback (XXS-XXL)

Molift Rgosling Ampu Highback (XXS-XXL)

Molift Rgosling Ambulating Vest

Molift Rgosling Comfort Highback (S-L)

Molift RgoSling Fabric Stretcher

Molift RgoSling Repositioning

Tillbehör:


Molift Rgosling Extension Loops

Art.nr.: 1721600


Molift Rgosling Ambulating Vest Groin strap (XXS-XL)

Transition coupling

Some rail systems have transition couplings ① installed. These are used for transferring the user from one rail system to another e.g. between two different rooms.


Components


- ① Transition coupling unit
- ② Transition coupling unit with pulley
- ③ Pulley with knob
- ④ Gate (two pcs per unit)

States


Connected

Disconnected


To disconnect

Before disconnecting the transition coupling, ensure trolley is min. 0.5 m away from transition coupling unit.

Pull pulley down until stop.

You hear a sound, which confirms disconnection.


The transition coupling is now disconnected.



Transfer - disconnected transition coupling

Do not perform any transfer of user before you have checked the following:


Both gates ① are visible (transition coupling unit seen from below).

Locking bolt ② must **not** be visible when coupling is disconnected.



Push traverse rail in the desired direction.


To connect

Before connecting the transition coupling units, ensure trolley is min. 0.5 m away from transition coupling unit.

Push the traverse rail towards the stationary rail, so the transition coupling units connect, and you hear a sound.

Transfer - connected transition coupling

Do not perform any transfer of user before you have checked the transition coupling units are completely locked together. Pushing the rails apart by hand must be impossible.

The user can now be transferred from one rail to another.

Cleaning / Disinfection

Clean on a regular basis. Clean surfaces with a damp cloth using an appropriate pH-neutral detergent. Do not use solvents or strong liquids, this may damage surfaces on the hoist. For disinfection when needed; use isopropyl alcohol. Avoid abrasive cleaning products. Check emergency stop and emergency lowering after cleaning. The lift should not be exposed to running water. When installing in environments with a risk of moisture, it is important that the rail system and the motor are regularly washed and dried. A complete guide to reducing the risk of corrosion is available on request.

 **Make sure not to damage or remove labels when cleaning.**

 **The hoist and hand control have IP24 protection class which indicates: >12.5mm protection against access to hazardous parts, and that water splashed against the enclosure from any direction shall have no harmful effects.**

 **Carer and user is exposed to the following materials: Polypropilene and polyester (slings*), PC/ABS (hand control and hoist cover), powder coated aluminium, powder coated zink (sling bars).
* = Applied part**

 **In rare cases user or carer exposed to polypropilene and polyester may suffer from allergic reactions such as skin irritation, mild itching and skin tenderness.**

Reconditioning

Follow cleaning and installation procedure, complete periodic inspection and use checklist after installation to recondition the hoist.

Recycling

Refer to "Recycling instructions" for how to properly dispose of product. This can be found on www.etac.com.

Spare parts

A list of spare parts is available on request.

Periodic inspection scope

Periodic inspection is a visual examination (particularly of the hoist's load bearing structure and lifting mechanism with attachments, brakes, controls, safety devices and person-support devices) according to Periodic Inspection Report for Molift Air 200 and 350.
See www.etac.com.

 **If there is any problems with the hoist that could jeopardize someones safety, the hoist shall immediately be taken out of service and marked "out of order". Do not use the hoist until it is repAired. Contact you local Molift/Etac dealer/distributer to report any problems.**

Periodic Inspection shall be performed at least once a year or more frequently if required by local requirements. The inspection must be performed by service personnel authorized by Etac. Contact Etac at molift@etac.com for training and authorization or recommendation of an approved service partner.

When performing a periodic inspection, the inspector shall fill out the inspection report for Molift Air. The reports should be retained by the person(s) responsible for servicing the hoist. If the inspection reveals defects and damages, the hoist shall be taken out of operation and marked "out of order" until it is repAired.

Troubleshooting

Symptom	Possible Cause/Action
The hoist does not respond to hand control action	Emergency Stop is activated. Deactivate by pushing button back in.
	The hoist's electronics is overheated. Wait for it to cool down.
	Hand control is not plugged in properly. Hand control or plug or cord can be broken an should be replaced.
The hoist does not respond to Hand control up button	Lifting strap angle sensor is activated. Adjust lifting strap or move hoist to reduce lifting strap angle.
When using Electrical Emergency lowering the lifting band goes upwards	Lifting band has gone completely out, and has been winded the wrong way. Use manual emergency lowering.

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Molift Rgosling Mediumback Padded (XS-XXL)
 Molift Rgosling Highback Padded (XS-XXL)
 Molift Rgosling Mediumback Net (XXS-XXL)
 Molift Rgosling Highback Net (XXS-XXL)
 Molift Rgosling Toilet Lowback (XXS-XXL)
 Molift Rgosling Toilet Highback (XXS-XXL)
 Molift Rgosling Amпу Mediumback (XXS-XXL)
 Molift Rgosling Amпу Highback (XXS-XXL)
 Molift Rgosling Ambulating Vest
 Molift Rgosling Comfort Highback (S-L)
 Molift RgoSling Fabric Stretcher
 Molift RgoSling Repositioning

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Molift Rgosling Mediumback Padded (XS-XXL)
 Molift Rgosling Highback Padded (XS-XXL)
 Molift Rgosling Mediumback Net (XXS-XXL)
 Molift Rgosling Highback Net (XXS-XXL)
 Molift Rgosling Toilet Lowback (XXS-XXL)
 Molift Rgosling Toilet Highback (XXS-XXL)
 Molift Rgosling Ampu Mediumback (XXS-XXL)
 Molift Rgosling Ampu Highback (XXS-XXL)
 Molift Rgosling Ambulating Vest
 Molift Rgosling Comfort Highback (S-L)
 Molift RgoSling Fabric Stretcher
 Molift RgoSling Repositioning

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- Molift Rgosling Mediumback Padded (XS-XXL)
- Molift Rgosling Highback Padded (XS-XXL)
- Molift Rgosling Mediumback Net (XXS-XXL)
- Molift Rgosling Highback Net (XXS-XXL)
- Molift Rgosling Toilet Lowback (XXS-XXL)
- Molift Rgosling Toilet Highback (XXS-XXL)
- Molift Rgosling Ampu Mediumback (XXS-XXL)
- Molift Rgosling Ampu Highback (XXS-XXL)
- Molift Rgosling Ambulating Vest
- Molift Rgosling Comfort Highback (S-L)
- Molift RgoSling Fabric Stretcher
- Molift RgoSling Repositioning

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Etac A/S