
Molift UnoSling Limblift

molift[®]
by Etac


User manual

BM30199 Rev. A 2022-11-01

MD

CE

Index

Labels and symbols.....	3
User manual – English.....	4
Bruksanvisning – Svenska.....	6
Brukermanual – Norsk.....	8
Brugervejledning – Dansk.....	10
Käyttöohje – Suomi.....	12
Gebrauchsanweisung – Deutsch.....	14
Handleiding – Nederlands.....	16
Manual d'utilisation – Français.....	18
Manual utente – Italiano.....	20
Manual de usuario – Español.....	22

Labels and symbols/ Etiketter och symboler/ Etiketter og symboler/ Mærkater og symboler/ Tarrat ja symbolit/ Beschriftungen und Symbole/ Labels en symbolen/ Étiquettes et symboles/ Etichette e simboli/ Símbolos y etiquetas:

Molift UnoSling LimbLift
 A patient specific, disposable sling
 Art. no: xxxxxxx
 DK HMI: xxxxxx

Size: XX

Material: Polypropylene, Polyester

Barcode according to GS1-128
 (01)xxxxxxxxxxxx(11)xxxxxx(21)xxxxxxxx
 (01):EAN number (11)YYMMDD(21):Serial number

Name: _____

UP/UPP OUT/UT

YYYY-MM-DD

1

Date of first use


Year	2022	2023	2024	2025	2026	2027
Month	01	02	03	04	05	06
	07	08	09	10	11	12

Etac A/S
 Parallevej 3
 DK-8751 Gedved

MD

CE UK

DO NOT WASH
 DO NOT BLEACH
 DO NOT TUMBLE DRY
 DO NOT IRON
 DO NOT DRY CLEAN

Symbols:

 EN - Manufacturer
 SV - Tilverkare
 NO - Produsent
 DA - Producent
 FI - Valmistaja
 DE - Hersteller
 NL - Fabrikant
 FR - Fabricant
 IT - Produttore
 ES - Fabricante

 EN - This Way up, this side is outside
 SV - Denna sida upp, denna sida utåt
 NO - Denne siden opp, denne siden ut
 DA - Denne side op, denne side er ydersiden
 FI - Tämä puoli ylöspäin, tämä puoli ulospäin
 DE - Hier oben, diese Seite ist außen
 NL - Deze kant boven, dit is de buitenkant
 FR - Vers le haut, côté extérieur
 IT - Lato alto, lato esterno
 ES - Este lado hacia arriba, este lado fuera

 EN - CE marked
 SV - CE-märkt
 NO - CE-merket
 DA - CE-mærket
 FI - CE-merkintä
 DE - CE-konform
 NL - CE-markering
 FR - Marquage CE
 IT - Marchio CE
 ES - Marchio CE

 EN - max user weight
 SV - Max. brukarvikt (säker arbetsbelastning)
 NO - Maks. brukervekt (SWL)
 DA - Maks. brugervægt (sikker arbejdsbelastning)
 FI - Käyttäjän enimmäispaino (SWL)
 DE - Max. Benutzergewicht (Tragfähigkeit)
 NL - Max. gewicht gebruiker (veilige belasting):
 FR - Poids max. de l'utilisateur :
 IT - Peso massimo utente (carico di sicurezza)
 ES - Peso máximo del usuario (CTS)

 EN - Refer to user manual
 SV - Se manual
 NO - Se bruksanvisningen
 DA - Se brugervejledning
 FI - Katso käyttöohje
 DE - Lesen Sie die Bedienungsanleitung
 NL - Zie de bedienings-handleiding
 FR - Consultez le manuel d'utilisation
 IT - Fare riferimento al manuale dell'utente
 ES - Consulte el manual del usuario

 EN - Medical device
 SV - Medicinteknisk produkt
 NO - Medisinsk utstyr
 DA - Medicinsk udstyr
 FI - Lääkinnällinen laite
 DE - Medizinprodukt
 NL - Medisch hulpmiddel
 FR - Dispositif médical
 IT - Dispositivo medico
 ES - Producto sanitario

 EN - Date of manufacture. YYYY-MM-DD
 (year/month/date)
 SV - Tillverkningsdatum. ÅÅÅÅ-MM-DD
 (år/månad/dag)
 NO - Produksjonsdato. ÅÅÅÅ-MM-DD
 (år/måned/dag)
 DA - Produktionsdato. ÅÅÅÅ-MM-DD
 (år/måned/dag)
 FI - Valmistuspäivä. VVVV-KK-DD
 (vuosi/kuukausi/Päivä)
 DE - Herstellungsdatum. JJJJ-MM-DD
 (Jahr/Monat/Tag)
 NL - Productie-datum JJJJ-MM-DD
 (jaar/maand/dag)
 FR - Date de fabrication. AAAA-MM-DD
 (année/mois/jour)
 IT - Data di produzione. AAAA-MM-DD
 (anno/mese/giorno)
 ES - Fecha de fabricación. AAAA-MM-DD
 (año/mes/día)

 EN - Single patient - Multiple use
 SV - En patient - flera användningar
 NO - Kun til bruk på én pasient - flegangsbruk
 DA - Enkelt patient - flegangsbrug
 FI - Potilaskohtainen - voidaan käyttää useita kertoja
 DE - Einzelner Patient - mehrfach anwendbar
 NL - Eén patiënt - meervoudig gebruik
 FR - Un seul patient - à usage multiple
 IT - Singolo paziente - uso multiplo
 ES - Un solo paciente - uso múltiple

 EN - UKCA marked
 SV - UKCA-märkt
 NO - UKCA-merket
 DA - UKCA-mærket
 FI - UKCA-merkintä
 DE - UKCA-konform
 NL - UKCA-markering
 FR - Marquage UKCA
 IT - Marchio UKCA
 ES - Marchio UKCA

Wash symbols:

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Molift UnoSling Limblift

The Molift UnoSling LimbLift is a consumable for single patient use.

Molift UnoSling LimbLift is a lifting band designed for different uses such as lifting up and holding a body part and providing support to the user. It may be used in situations such as treatment or examination of pressure ulcers or foot care.

Molift UnoSling LimbLift is available in two sizes, Standard and Wide.

Definition

The person described as “the user” in this Instruction for Use is the person who is lying in the product. “The carer” is the person who applies the product.

Intended purpose

To be used with a hoist for lifting a body part on one specific user for a short time; e.g. a leg or arm during examination, repositioning, dressing, undressing or similar.

Intended target group

Used for users who require support to lift a body part for a short time; e.g. persons with limited ability to move and who require body support during positioning or dressing. Not an exhaustive list.

Intended environment

For indoor use in acute care, long term care, home care.

Indications

The device can be used for children, adults and elderly people when there is a need for a patient specific product.

The target group for the device is based on individual health and mobility function, and not a specific diagnose or age.

Contra indications

Make sure the sling fits the user. The sling is non-washable and should be disposed when dirty or when patient leaves the facility.

The Molift UnoSling LimbLift sling is a consumable for single patient use, it has a limited service life.

Write the user’s name and the date of first use on the product label before using the product.


The Molift UnoSling LimbLift is manufactured from non-woven material. Ensure that no pointed or sharp objects come into contact with the sling. Avoid pulling hard on the sling when applying it to or removing it from the user.

The Molift UnoSling LimbLift should be taken out of use:

- when it has become dirty or damp, or if there is a risk that the sling has been exposed to contamination;
- if it has been damaged;
- when the user no longer requires it.

The sling should be disposed of in accordance with applicable procedures for the health care facility.

Standard size illustration

Wide size illustration

Things to remember


- Read through the manual before starting to use the sling.
- It is important to understand the contents of the user manual before using the equipment. Read the manual for both hoist and sling.
- Molift slings shall only be used with persons. Never use the sling to lift or move objects of any kind.
- Make a decision on whether one or more assistants are required. Plan the lifting operation in advance to ensure that it is as safe and smooth as possible. Remember to work ergonomically.
- Lifting and transferring a person’s body part always involves a certain risk. Only informed personnel should use the equipment described in this user manual. If the maximum load (SWL) differs between the sling bar and the sling, the lowest maximum load should always be used.
- Carry out a risk assessment and document it. As the assistant, you are responsible for the safety of the user.
- Never leave a user unattended in a lifting situation.
- Remove the sling carefully to ensure that there is no sudden drop or fall.

The following inspection should be carried out daily or before use


- The user must have the bodily functions required for the use of the sling.
- Make sure the sling is suitable for the hoist/sling bar to be used.
- Check that the sling has no visible damage or faults.

Sling checkpoint

Checkpoint	Troubleshooting
1 The sling is placed under user's limb as required with no folds.	If the sling is placed with folds it risks discomfort and pressure on the user's body part.
2 Do not lift higher than the user's mobility of the joints can sustain.	Lifting the user higher than what the joints can sustain poses the user to a risk of injury.
3 The lifting straps are stretched simultaneously when the lift is started. See 'Lifting checkpoints' below.	If the straps are stretched unequally, the sling may be pulled out of position, and there may be a risk of the user's body part sliding out.


Lifting checkpoints

- 1 Check that the sling is correctly fitted around the user's limb and that the lifting loops are correctly fitted to the sling bar hooks.
- 2 Start lifting. Ensure that the lifting straps are stretched tense, and pause before lifting. Check that both lifting loops are attached correctly to the sling bar before lifting the user.
- 3 Carry out the planned LimbLift.


Sling application

- 1 Use the application pocket on the sling for easier fitting.
- 2 Push the sling under the user so it is possible to pull out on the opposite side.
- 3 Pull the sling and center it under user's body part(s).
- 4 Make sure the lifting straps are fastened with equal length on right and left side.

Removing the sling

- 1 Lower the user's body part. Release the lifting band from the sling bar.
- 2 Remove the sling by folding the lifting band under the sling itself and pull out.


The Molift UnoSling LimbLift, as described in this manual, is CE labeled and complies with the Council of Europe's directive on medical devices MDR (EU) 2017/745, class 1, and has been tested and approved in accordance with the EN ISO 10535:2021 standard.

- This symbol indicates especially important safety information. Follow these instructions carefully.
- Do not wash. If the sling is washed, this symbol disappears and another symbol appears. The sling should then be disposed of.
- Do not use. If this symbol appears the sling has been washed and should be disposed of.
- The sling should not be stored in direct sunlight.

Safe Working Load (SWL): 200 kg (440 lbs)
Material: Polypropylene, polyester

2-point sling bar	Size of sling	
	Standard	Wide
Small 340 mm	✓	✓
Medium 440 mm	✓	✓
Large 540 mm	✓	✓

Item no.	Sling model	Size
3100020	Molift UnoSling LimbLift Single-patient sling, consumer item	Standard
3100030	Molift UnoSling LimbLift Single-patient sling, consumer item	Wide

Size in mm	A	B	C
Standard	380	950	150
Wide	380	950	250

Visit www.etac.com to download documentation to ensure that you have the latest version.

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Etac Ltd, Unit 60