
# User manual

Molift RgoSling StandUp

Molift RgoSling StandUp with support

**molift**<sup>®</sup>  
by Etac

BM40399 - 2024-02-22 - Rev. 4.0

[www.etac.com](http://www.etac.com)

	<b>Symbols</b>	<b>4</b>
	<b>Device overview</b>	<b>8</b>
	<b>General</b>	<b>10</b>
	English	10
	Svenska	14
	Norsk	18
	Dansk	22
	Suomi	26
	Deutsch	30
	Nederlands	34
	Français	38
	Italiano	42
	Español	46
	<b>Handling</b>	<b>50</b>
	Sitting application	51
	Sitting removal	53
	Warnings	54
	Information	55
	<b>Combination list</b>	<b>58</b>
	<b>Sling measurements</b>	<b>58</b>
	<b>Size guide</b>	<b>58</b>

**EN:** Warning: This symbol appears in the user manual along with a number referring to the instructions below. It draws attention to situations where the product or the user or carer's safety may be in danger.

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**EN:** Information: This symbol appears in the user manual along with a number referring to supplementary information which can be found in the information section.

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




## General

Thank you for choosing an Etac device.

This manual must be read thoroughly before you start using the device as it contains important safety information and gives instructions on how to use the device correctly. Visit [www.etac.com](http://www.etac.com) for all relevant product documentation, for example the user manual, periodic inspection document and recycling instructions. You can select your language via the “International” and “Local websites” link.

You can also consult the QR code on the product label for further information and instructions.

### Copyright

This document is the property of Etac A/S. The content shall not be modified in whole or part without prior approval of Etac A/S.

### Disclaimer

Etac continuously makes improvements to our devices and information. We therefore reserve the right to modify our devices and instructions for use without prior notice. Consult [www.etac.com](http://www.etac.com) to find the latest version of the documentation.

## CE Declaration of Conformity

The device mentioned in this manual and the different models hereof comply with the Medical Device Regulation (EU) 2017/745.

The device is classified as a Class I medical device.

### Product standards

The device has been tested and complies with the requirements applicable for this type of device in the below mentioned standards:

- ISO 10535:2021 Assistive products – Hoists for the transfer of persons – Requirements and test methods.

### Notice

Any serious incident that has occurred in relation to the device should be reported to your distributor/representative and the national competent authority within your region in a timely manner. The distributor/representative will forward the information to the manufacturer.

### Pre-sales information

Please visit [www.etac.com](http://www.etac.com).

## Device description

Etac supplies a wide selection of slings for different types of transfers. The Molift RgoSling series is designed to provide high comfort for users by evenly distributing weight and pressure. Molift RgoSlings are available in different sizes, models and materials. See [www.etac.com](http://www.etac.com) for a complete overview of possible variants.

Molift RgoSling StandUp can be used for transferring a user from a sitting to a standing position and vice versa with a sit-to-stand hoist with a rope lock lifting arm.

Molift RgoSling StandUp with support is intended for users with extra need for support when standing up. This model comes with a padded support band for the buttocks to provide additional support.

### Benefits

The sling is padded to distribute the weight, which provides a comfortable support for the user through the entire transfer.

The sling has a waist belt with non-slip lining inside which helps keep it in place.

Using the sling provides a high level of safety and security for the user.

### Materials

Fabric: Polyester

Padding: Polyethylene

## Introduction

The person referred to as “the user” is the person who is sitting in the product. “The carer” is the person who applies the product.

### Intended purpose

The sling is an assistive device intended for alleviation of or compensation for a functional impairment due to an injury or disability. The device is designed for an individual lacking the ability to stand up and transfer themselves over shorter distances to another sitting position to/from a bed, a wheelchair, a chair, a toilet or similar, due to reduced mobility or physical strength.

### Intended user

The device can be used for children, adults and elderly people with a need for support to stand up and be transferred over shorter distances between sitting positions.

### Intended environment

For indoor use in acute care, long term care and home care.

### Indications

For users who are in need of support in order to stand up and be transferred to another sitting position to/from a bed, a wheelchair, a chair and a toilet or similar.; e.g. people suffering from pain, reduced range of motion in their joints, stiffness, and muscle weakness. Not an exhaustive list.

The user group for the device is based on individual health and mobility function, and not on a specific diagnosis or age.

### Contraindications

There are no known contraindications.

### Precautions

Make sure the sling fits the user and the sling bar. The user should be able to follow instructions, maintain a good core stability and load at least one leg.

### Training requirements

Only carers who have read and understood the user manual are allowed to use the device.

## Safety information

The safety of both user and carer must never be jeopardized – if you suspect that the device has a malfunction you must put it in quarantine marked “out of order” and call your local distributor/representative.

Safety signs, safety measures and warnings symbols are given in the context relevant for the actual handling and/or operating situation – it is mandatory the given instruction is followed to avoid dangerous/hazardous situations.

Signs, symbols and instructions placed on this device must never be covered or removed and must remain present and clearly legible throughout the entire lifetime of the device. Immediately replace or repair illegible or damaged signs, symbols, and instructions. Contact your distributor/representative for instructions.

### Patient transfer

Transfer of a user will always pose a certain risk. Therefore, a risk assessment must always be carried out before the transfer is performed.

Never leave a user unattended when hoisted.

Molift slings shall only be used to hoist persons. Never use the sling to hoist objects of any kind.

### Risk assessment

Do a risk assessment of the sling, the individual user and the intended transfer, and document your analysis. You as the carer are responsible for the safety of the user.

Make a decision on whether one or more carers are required.

Plan the transfer in advance to ensure that it is performed as safe and smooth as possible.

Remember to work ergonomically according to local recommendations.

Before using the sling, mark the date of first use on the periodic inspection label.

## Slings

Always check – before each use:

- that the user has the required ability for the sling.
- that the model, size and fabric are suitable for the user and the transfer situation.
- that the sling is suitable for the hoist and the sling bar.
- that the sling does not have visible damage or frays to seams, fabric, straps and loops.
- that the loop connections are correctly attached to the sling bar hooks before hoisting the user.
- that the wheels on the wheelchair or bed are locked during the transfer. If using a mobile hoist the wheels on the hoist must not be locked.

If the inspection reveals any safety risks, the sling must immediately be taken out of service and be disposed of.

Never use a faulty or damaged sling as it can break and cause personal injury. Destroy and discard damaged and old slings.

## Materials

All used materials which are to be in direct contact with skin are biologically assessed and evaluated according to relevant standards. However, users/carers who are allergic to any of the materials are recommend not to use the device.

## Accessories

User manuals supplied with accessories must be read and understood before use. This includes any limitations in use and combinations. Information about accessories can be found on [www.etac.com](http://www.etac.com).

## Combinations

This device can be combined with Molift sling bars and hoists. See the combination list on page 58.

Statement on sling and hoist compatibility for other manufactures of hoists can be found on [www.etac.com](http://www.etac.com).

When combining the device with other devices, the combination must not change the intended purpose of the devices or modify the devices in such a way that compliance with the applicable requirements may be affected. The person/company responsible for the combination of the devices must ensure that the applicable requirements are fulfilled.

If safe working load (SWL) differs between hoist, sling bar and sling, the lowest load should always be used.

## Troubleshooting

Problem	Action
There is more than 10 cm between the left and right side padding of the waist belt.	Try a larger size.
The waist belt is too loose after tightening it.	Try a smaller size.
The sling is sliding up.	Tighten the waist belt properly. Make sure the sling is positioned low enough on the user's back when applying it. Make sure the user leans backwards when hoisting. Otherwise, evaluate whether the user is able to follow the instructions or requires a passive hoist. Make sure the user is not wearing slippery clothes.
The user is not able to take load on the legs.	Evaluate whether the user has the physical ability for this type of sling. Otherwise, the user requires a passive hoist.
The ropes are not loaded at the same time.	Make sure the ropes are attached symmetrically. Make sure the centre of the sling is at the midline of the user.

## Service, cleaning and maintenance

Maintenance and periodic inspections of the device are wholly the responsibility of the owner of the device. Failure to maintain a device in accordance with the instructions may compromise the safety of the user and/or carer and may affect the expected lifetime of the device.

### Expected lifetime

The device has an expected lifetime of 1 to 5 years under normal use. The lifetime of the device may vary depending on usage frequency, loads and washing/drying frequency and conditions.

Washing and drying at high temperatures wear out the material faster.

### Periodic inspection

Periodic inspection must be carried out at least every 6 months. More regular inspection may be required if the sling is used or washed heavily. See the periodic inspection label on the sling and the periodic inspection document at [www.etac.com](http://www.etac.com).

### Cleaning

60-85°  
140-185F

Max. 60°  
/140F

## Transportation, storage and disposal

### Transportation and storage

When the device is not in use, it can be stored under the following conditions:

### Disposal

The device contains no hazardous substances that require special disposal. Dispose of the device according to the national provisions that apply.

### Recycling

See the recycling instructions on [www.etac.com](http://www.etac.com), and follow the national provisions that apply.

## Technical data

### Size guide

Measurements in the table on page 58 should be used as a reference. The correct sling size depends on the user's weight, function and body shape. Consult the size guide, and measure the user's relevant measurements in a sitting position. Try the sling with the most correct measurements to make sure it fits.

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**EN:** Information

01. Slide the sling down behind the user's back until it reaches the user's waist.
02. Instruct the user to lean back slightly when lifting, and if possible hold onto the handles.
03. It is recommended to be two assistants when applying RgoSling with support. Tighten the support so it supports when lifting. It must not take up the whole load.
04. Make sure the lifting ropes are fastened symmetrically on left and right side.
05. Before performing the transfer – make sure that:
  - Sling is placed at waist height.
  - Waist belt is locked and tightened.
  - Both ropes are attached to the lifting arm.
  - Waist belt is retightened when the sling is loaded.
  - When using the Molift RgoSling StandUp w/ Support, the support must support the user's buttocks.
06. When using RgoSling Standup with support, the support band may be removed when the user is in standing position.

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**Etac Ltd.**