

A copy of this document is available for download at www.etac.com.

Scan QR code on product label for accessing product home page. The user manual can be found in the document section.

Article no. (see ID label):

.....

Serial no. (see ID label):

.....

**Production Year:** 

.....

Owner:

.....

Sling shall be inspected separately and is not included in this inspection.

## Service light indicator

Note colour of service light indicator:

Status	Service light	Mode
	(No light)	Power saving (Stand by)
	Green	Ready for use
	Yellow	Order service
	Red	Perform service
	Red (and sound when activated)	Perform service immediately

### Situation of use

Home Hospital Nursing home Other



The periodic inspection is to be performed by a person who is suitable and properly qualified and well acquainted with the design, use and care of the hoist



The hoist's software Molift Service tool provides option to enter date of PI and name of inspector along with additional information. Molift recommends adding all observations / remarks concerning the hoist in the Molift Service tool for complete history log of the hoist

# Visual examination



Visual examination of load bearing structure to make sure there is no damage, cracks, frays or deformation. All checkpoints must be checked to approve hoist for further use

OK Not OK

Lifting column

Lifting arm

Product label

Plastic covers

Hand control

Wheels / castors

Brakes on wheels / castors

Cables

Battery

Other electronics

Sling bar

Leg spreading mechanism

No corrosion

Weldings

Instruction for use / user manual is present / available


# Functional examination

	OK Not OK
Test FUNCTION and check for wear and damage. All checkpoints must be checked to approve hoist for further use	Battery. With fully charged battery inserted, after initial cycle, bat- tery light is off, and service light indicator is green
OK Not OK	Accessories:
Wheels / castors	
Wheel / castor brakes	
With legs in inner position (par- allel): Check that the distance between legs can not decrease when pushed together by hand	
Lifting column is completely fixed when locking handle is engaged	Performed by Full name:
Emergency stop button is fully operational and interrupts power when engaged. Reset switch	Date/Place:
Tool for manual emergency lowering	
Electrical emergency lowering	Signature:
Manual emergency lowering	
Hand control (buttons, strain relief on hand control cable, hook, and battery indicator)	Approved without faults
Sling bar coupling (swivel joint moves freely and without noise)	Next inspection (YYYY / MM):/
If all points so far are "OK", hoist shall be tested with SWL. Perform one lifting cycle with weight equal to SWL (see prod- uct label for value kg/lbs).	Hoist is marked with "Out of order" and sent for repair Hoist is not eligible for repair and taken out of service
Motor and spindle or must not slip when loaded (causing lifting arm to lower), and movement must be smooth and even	If periodic inspection reveals any defect, wear or other damage that jeopardises the safety of the patient the hoist may not be used until the deficiency has been eliminated
Leg spreading mechanism, outer and inner position with SWL	The owner is notified
Perform new visual control; Damage, play and deformations as described above. Any dam- aged parts must be repaired or replaced and test is repeated once more until the hoist performs correctly Charger is operational	



#### Due to legal requirements this document or a copy must be filed with the owner's service log book

If approved without faults, apply inspection label and mark with month and year for **next** inspection. When ordering label, use item no. 1100306



## Notes and observations

.....