
Molift Smart 150




Periodic inspection checklist (EN)

In accordance with ISO:10535

PI09201 Rev. F 2022-06-16





A copy of this document is available for download at www.etac.com.

Scan QR code on product label for accessing product home page. The user manual can be found in the document section.

Article no. (see ID label):

Serial no. (see ID label):

.....

.....

Production Year:

.....

Owner:

.....

Sling shall be inspected separately and is not included in this inspection.

Service light indicator

Note colour of service light indicator:

Status	Service light	Mode
	(No light)	Power saving (Stand by)
Green		Ready for use
	Yellow	Order service
	Red	Perform service
	Red (and sound when activated)	Perform service immediately

Situation of use

Home Hospital Nursing home Other



The periodic inspection is to be performed by a person who is suitable and properly qualified and well acquainted with the design, use and care of the hoist



The hoist's software Molift Service tool provides option to enter date of PI and name of inspector along with additional information. Molift recommends adding all observations / remarks concerning the hoist in the Molift Service tool for complete history log of the hoist

Visual examination



Visual examination of load bearing structure to make sure there is no damage, cracks, frays or deformation. All checkpoints must be checked to approve hoist for further use

OK Not OK

	Lifting column
	Lifting arm
	Legs
·	Product label
	Plastic covers
	Hand control
	Actuator
	Wheels / castors
	Brakes on wheels / castors
	Charger
	Cables
	Battery
	Other electronics
	Sling bar, incl. padding
	Leg spreading mechanism
	No corrosion
	Weldings
	Instruction for use / user manual is present / available

Functional examination



Test FUNCTION and inspect for wear and damage. All checkpoints must be checked to approve hoist for further use

service light is green

OK	Not OK		
		Wheels / castors	
		Rear wheel / castor brakes	Performed by
		With legs in inner position (par- allel): Check that the distance between legs can not decrease when pushed together by hand	Full name:
		Lifting column is completely fixed when locking handle is engaged	Date/Place:
		Hand control (buttons, strain relief on hand control cable, hook, and battery indicator)	Signature:
		Electrical emergency lowering	
		Manual emergency lowering	Approved without faults
		Emergency stop button is fully operational and interrupts power when engaged. Reset switch	Next inspection (YYYY / MM):/
	Sling bar coupling (swivel joint moves freely and without noise) If all points so far are " <i>OK</i> ", hoist shall be tested with SWL.	for repair Hoist is not eligible for repair and taken out of service	
		Perform one lifting cycle with weight equal to SWL (see prod- uct label for value kg/lbs). Actuator must not slip when loaded (causing lifting arm to lower)	If periodic inspection reveals any defect, wear or other damage that jeopardises the safety of the patient the hoist may not be used until the deficiency has been eliminated The owner is notified
		Leg spreading mechanism, outer and inner position with SWL	
		Charger is operational	
		Battery. With fully charged battery inserted, after initial cycle, battery light is off, and wrench/	

Accessories:


.....

.....



Due to legal requirements this document or a copy must be filed with the owner's service log book

If approved without faults, apply inspection label and mark with month and year for **next** inspection. When ordering label, use item no. 1100306



Notes and observations

.....