

Molift Nomad


molift[®]
by Etac

Periodic inspection checklist (EN)


In accordance with ISO:10535

PI15001 Rev. C 2022-06-17

Components



Product label



SWL value

ID label



Serial no.

Article no.

Production date



Article no.

Production date
(YYYY-WW)

Serial no.

A copy of this document is available for download at www.etac.com.

Scan QR code on product label for accessing product home page. The instruction for use can be found in the document section.

Article no. (see ID label):

.....

Serial no. (see ID label):

.....

Production Year:

.....

Owner:

.....

Sling shall be inspected separately and is not included in this inspection.


Service light indicator


Note colour of service light indicator:


Status	Service light	Mode
	(No light)	Power saving (Stand by)
	Green	Ready for use
	Yellow	Order service
	Red	Perform service
	Red (and sound when activated)	Perform service immediately

Situation of use


- Home
- Hospital
- Nursing home
- Other

 **The periodic inspection is to be performed by a person who is suitable and properly qualified and well acquainted with the design, use and care of the hoist**

 **The hoist's software Molift Service tool provides option to enter date of PI and name of inspector along with additional information. Molift recommends adding all observations / remarks concerning the hoist in the Molift Service tool for complete history log of the hoist**

 **This inspection document covers hoist only. Rail system must be inspected according to PI19201**

Visual examination

 **Visual examination of load bearing structure to make sure there is no damage, cracks, frays or deformation. All checkpoints must be checked to approve hoist for further use**

OK Not OK

	Lifting belt with hook. Must be replaced every 5 years. See lable on belt.
	Safety hook
	Safety hook for extension arm (if applicable)
	Labels, both sides
	Plastic covers
	Hand control and cable
	Lifting hooks
	Charger
	Instruction for use / user manual is present / available

Functional examination



Test FUNCTION and inspect for wear and damage.
All checkpoints must be checked to approve hoist for further use

<i>OK</i>	<i>Not OK</i>
	Emergency stop button is fully operational and interrupts power when engaged. Reset switch. Service light indicator must turn green
	Hand control (buttons, wire sockets, and battery indicator)
	Hand control functionality
	Up/down buttons on hoist
	End stop
	Lifting belt with safety hook. Must be replaced every 5 years. See lable on belt.
	Belt checked for damage, and tear and wear
	Safety hook checked for damage, and tear and wear
	If applicable: Safety hook for extension arm checked for damage, spring function, loose parts, and tear and wear
	Swivel joint at trolley moves freely when turned 360°
	Checked that the hoist functions correctly and runs without load all the way up and down without abnormal noise
	If all points so far are "OK", hoist shall be tested with SWL. NOTE! Ensure SWL for the complete hoist system exceeds or equals hoist SWL before any load testing
	Perform one lifting cycle with weight equal to SWL with min. range of 500 mm in the top of lifting range. Brake must not slip when loaded. Movement must be smooth and even, and without unusual noise

<i>OK</i>	<i>Not OK</i>
	Move hoist loaded with SWL through entire rail system. Check deflection and suspension is acceptable
	Perform new visual control; Damage and deformations as described above. Any damaged parts must be repaired or replaced and test is repeated once more until the hoist performs correctly
	Charger is operational
	Battery. Battery light is OFF and service light indicator is green

Accessories:

.....

.....

.....

Performed by

Notes and observations

Full name:

Date/Place:

Signature:

Approved without faults


Next inspection (YYYY / MM):/.....

Hoist is marked with "Out of order" and sent for repair

Hoist is not eligible for repair and taken out of service


If periodic inspection reveals any defect, wear or other damage that jeopardises the safety of the patient the hoist may not be used until the deficiency has been eliminated.

The owner is notified



Due to legal requirements this document or a copy must be filed with the owner's service log book

If approved without faults, apply inspection label and mark with month and year for **next** inspection. When ordering label, use item no. 1100306



Mark label with month and year for next inspection