



Immedia 2Move

Transfer board
for horizontal transfer

 **etac**[®]

Immedia 2Move

Transfer board for horizontal transfer

Immedia 2Move

The transfer board Immedia 2Move is used to bridge the gap between two horizontal surfaces and alters the height difference if needed. The transfer board is used for horizontal transfer between beds, bed to shower trolley, bed to X-ray table, stretcher to bed, etc.

2Move is easy to transport and store due to its low weight of only 1,5 kg and since it is foldable.

The handles facilitate carrying and make it easy to hang up the transfer board.

The foam material is comfortable for the user to be transferred on and the surface is easy to clean and disinfect according to general standards.

2Move is always used in combination with a low friction cover (accessory) which slides around the board, minimising friction and facilitating the transfer.

Intended Use

- Bridges the gap when transferring between two supine positions
- Transfer between beds
- Transfer between bed and shower trolley
- Transfer between bed and stretcher, X-ray table, operating table
- Can be used in X-ray and scanners



Immedia 2 Move Large foldable

Options



Small non-foldable

Immedia 2Move Small can be used when the user is partly dependent with either upper or lower body weakness.



Large foldable

Immedia 2Move Large foldable can be used for all horizontal transfers and is easy to transport and store since it is foldable.



Large non-foldable

2Move large non-foldable is similar to the Large foldable, but is rigid - can not be folded.



X-Large foldable

X-Large foldable is a wider option to the large transfer board and can be used for either wider transfers or for large or heavy users.

Tutorial - how to transfer a user from a bed to a stretcher



1. Place the Immedia 2Move halfway under the users body. Place the user's head and feet on the Immedia 2Move.



2. Place the two horizontal surfaces close to each other. One caregiver tilts and pushes the user towards the Immedia 2Move. A second caregiver pulls the user towards the stretcher.



3. When the transfer is complete, remove the Immedia 2Move.

Accessories



Immedia 2Move nylon cover

Immedia 2 Move cover is a low friction cover which slides around the board. It minimises friction and facilitates the transfer. The cover is available in three sizes. Material: Waterproof PU-nylon



Immedia 2Move Disposable cover

The cover is available in three sizes. Small: 100 pcs/roll Large: 40 pcs/roll Extra large 40 pcs/roll Material: PE-LD



Immedia 2Move Display

The wall display assures the accessibility of the transfer board and has a dispenser for the disposable cover.

Immedia 2Move

		Measure	Item no.
Small	Non foldable	500 x 950 mm	IM8050095
Small set	Non foldable	incl IM8050095 + IM4047	IM4319
Large	Non foldable	500 x 1900 mm	IM8150190
Large set	Non foldable	incl IM8150190 + IM4048	IM4320
Large	Foldable	500 x 1900 mm	IM8050190
Large set	Foldable	incl IM8050190 + IM4048	IM4318
Extra large	Foldable	700 x 1900 mm	IM8070190
Extra large set	Foldable	incl IM8070190 + IM4049	IM4321

Accessories

			Item no.
Cover - Small	Nylon	515 x 790 mm	IM4047
Cover - Large	Nylon	515 x 1610 mm	IM4048
Cover - Extra large	Nylon	715 x 1610 mm	IM4049
Cover - Small	Disposable	600 x 1000 mm	IM4035
Cover - Large	Disposable	600 x 2100 mm	IM4034
Cover - Extra large	Disposable	800 x 2100 mm	IM4036
Wall display			IM4041

Material

Board: Plastazote

Weight

1,5 kg.

Maximum user weight

150-300 kg*

*depend on the distance between the horizontal beds

Washing instructions transfer board and display



Clean the product with a solvent-free cleaning agent with a pH level between 5-9, or with a 70% disinfectant solution.

Washing instructions nylon cover



Wash as instructions. Do not use fabric softener or bleach



Etac is a world-leading developer and provider of ergonomic assistive devices and patient handling equipment. Our heart lies in the solutions that optimise quality of life for the individual, their family and caregivers.

For the latest news and continuously updated product information – visit: www.etac.com

Etac AB
Etac Immedia
Egeskovvej 12
DK-8700 Horsens
info@etacimmedia.com
www.etac.com


Creating Possibilities