
Transportation in motor vehicles

R82[®] 
by Etac

ISO 7176-19:2022

M1480 - 2025-09-08 - Rev. 1.0

Content	
ENGLISH	5 - 11
DEUTSCH	12 - 18
NEDERLANDS	19 - 25
DANSK	26 - 32
NORSK	33 - 39
SVENSKA	40 - 46
SUOMI	47 - 53
ESPAÑOL	54 - 60
FRANÇAIS	61 - 67
ITALIANO	68 - 74
PORTUGUÊS	75 - 81
PORTUGUÊS DO BRASIL	82 - 88
РОССИЯ	89 - 95
POLSKI	96 - 102
ČESKY	103 - 109
中國	110 - 116
日本語	117 - 123
ΕΛΛΗΝΙΚΑ	124 - 130
ÍSLENSKA	131 - 137
 Dahl docking system mounting on Kudu	138 - 177

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Transportation in motor vehicles

The instructions on how to prepare the seat and chair, must be carried out before transportation. The user should transfer to the vehicle seat and use the vehicle-installed restraint system whenever possible, and the unoccupied wheelchair should be stored in a cargo area or secured in the vehicle during travel.

Occupant restraint

Use an approved vehicle anchored 3-point belt certified according to ISO 10542-1 in the product. The pelvic restraint should be worn low across the front of the pelvis, so that the angle of the pelvic restraint is within the preferred zone of 30° to 75° to the horizontal, similar to that shown in figure 1. A steeper (greater) pelvic belt angle between 45° to 75° to the horizontal is preferred.


Figure 1

The back must be placed in vertical position and the seat plane must be horizontal. The belt restraints should be adjusted as tightly as possible, consistent with user comfort. Further, the belt webbing should not be twisted when in use. Belt restraints should not be held away from the body by wheelchair components or parts, such as the wheelchair armrests or wheels. See figure 2.


Figure 2

Shoulder belt restraints should fit over the shoulders, similar to the illustration provided in figure 3.


Figure 3

An updated table covering configurations for seats and frames, max. load/user weight in transportation can be found in the overview “frame_and_seat_combinations.pdf” at: www.etac.com/support/support-center/support-documents/

Check the user manual or assembly instructions for the product, for information about:

- the minimum weight limit
- the securement points
- maximum load for the seat or frame (can also be found on the product label on the product).

The wheelchair’s size and turning radius will affect its ease of access to and maneuverability in motor vehicles.

Small wheelchairs and/or wheelchairs with a shorter turning radius will generally provide greater ease of access and maneuverability.

This product has no provision for wheelchair mounted belt restraints.

Compatible securements:

The wheelchair’s securement points comply with the requirements of ISO 7176-19:2022 Annex B, subclause B2 Geometric specifications.

All approved fittings and securements that are compatible with this specification can be used for securement of the wheelchair during transportation.

Approved 4-point strap-type tie-down system

The wheelchair/product must be placed in a forward-facing position, when used as a seat in a motor vehicle.

The front wheels need to be turned under the frame so as to be in line with the frame before fixation.

Use an approved 4-point strap-type tie-down system according to ISO 10542-1 or ISO 7176-19,

to secure the product. Use either a hook or a strap in the fitting. See figure 4.


Figure 4

The location of all wheelchair securement points are marked with a label. See figure 5. It is only allowed to use the marked securement points.

Figure 5

Approved docking system (Dahl MADS™)

Use a Dahl MADS™ adaption kit for Kudu (REF 503036) in combination with either docking systems Docking MK II (REF. 501750) or Dahl VarioDock™ (REF. 503600) from Dahl Engineering. For ordering, instructions and user manuals for installation, service and use of the Adaption kit and docking system, please refer to the manufacturer home page: www.dahlengineering.dk or directly to the overview of approved wheelchairs and manuals: www.dahlengineering.dk/produkter/dahl-docking-systemer/godkendte-koerestole-og-manualer/


Figure 6 - Dahl MADS™ adaption kit for Kudu

Installation of the Dahl Docking stations in the vehicle

Only professional companies in the business of converting or building wheelchair accessible vehicles can order the docking systems from Dahl Engineering. A qualified and experienced technician must carry out the installation.

Dahl Engineering can provide vehicle specific installation instructions for a large range of vehicles, which must be respected by the installer.

Changes or replacements must not be made to the anchorage points/car fastenings for docking systems on the wheelchair or 4-point strap tie down systems, or to constructional elements or parts of the frame, without approval from the manufacturer.

Please contact Dahl Engineering for further information about approved vehicles and installation positions.

Dahl Engineering contact details are available at: www.dahlengineering.dk

WARNING

Incorrect installation of the MADS adaptation kit on the wheelchair could cause serious injury or death in a crash.

Always read R82's user manual and mounting instructions for your Kudu and Dahl's manuals for the mounting and use of the Dahl Docking System before use.

Do not drive the vehicle until wheelchair and user are correctly secured.

Identification of the Dahl Docking system

If your wheelchair is equipped with the locking plate for the Dahl Docking systems, the following label will be present on your wheelchair.

See figure 7.


Figure 7

The height noted on the label indicates the value that the height adjustable VarioDock™ control panel display must be set to.

For the MK II, this is the height that the docking unit must be installed at.

Failure to set/install the docking unit to the correct height might impede the locking plate from sliding smoothly into the docking unit or prevent the Docking station and wheelchair from connecting altogether.

For the Docking Station to function properly the Locking Plate must be mounted horizontally. When the entire assembly process is completed measure this height H with the user sitting in the wheelchair.

Make sure the tire pressure is correct beforehand!

Securing the wheelchair in the MK II and VarioDock™ docking stations:

Maneuver the wheelchair slowly and in a straight line over the docking station. The locking plate positioned under the wheelchair helps to guide the wheelchair into place in the docking station. When the locking plate is fully engaged in the docking station, a spring-action locking pin automatically secures the lock plate.

The Dahl docking stations are equipped with a control switch that indicates whether the locking plate is correctly secured in the docking station. As soon as the locking plate comes into contact with the locking pin, a warning tone will sound (a high-pitched howl), and the red diode/lamp/LED, in the control panel will light up until the lock plate is either fully engaged or the wheelchair is removed from the docking station.

As an indication that the wheelchair locking plate is fully entered into the docking unit and properly secured, the warning tone will cease, the red lamp/LED in the control panel will go out and the green lamp/LED will light up.

See pages 174 to 177.

Warning!

Do not move the vehicle:

- While the wheelchair is being maneuvered into position in the docking station;
- If the wheelchair and user are not correctly secured;
- If the warning tone sounds and/or the red warning lamp (LED) in the control panel flashes or is lit.

Always check if the lock plate is properly engaged in the docking station by trying to reverse the wheelchair out of the docking station before moving the vehicle.

It must not be possible to reverse out of the docking station without pressing the red release button in the control panel.

Buckle the seat belt before driving.

General occupant restraint instructions

Use a Dahl 3-point occupant restraint system to secure the occupant;

- Both pelvic and upper torso restraint belts must be used to restrain the occupant to reduce the possibility of head and chest impacts with the vehicle components;
- Any wheelchair anchored occupant restraint i.e. 3-point belt, harness or postural supports (lap straps, lap belts) should not be used or relied on for occupant restraint in a moving vehicle, regardless if they are labeled ISO 7176-19, ISO 10542-1, SAE J2249 or any other standard. Use a vehicle anchored and certified occupant restraint system instead.




Use a suitable positioned headrest during transport.

Restraint belts must be adjusted as tightly as possible consistent with user comfort. Restraint belt webbing must not be twisted when in use. Care should be taken when applying the occupant restraint to position the seatbelt buckle so that the release button will not be contacted by wheelchair components while driving or during a crash.

Positioning the occupant restraint when using it with the Dahl Docking systems

The pelvic restraint belt must be worn low across the front of the pelvis so that the angle of the pelvic belt is within the optional or preferred zone of 30° to 75° to the horizontal.

A steeper (greater) angle within the preferred zone, 45° to 75° is desirable i.e. closer to, but never exceeding 75° degrees. See figure 8.

Dimensions in millimeters

Figure 8

Proper belt restraint fit

Occupant belt restraints should make full contact with the shoulder, chest, and pelvis.

Pelvic belts should be positioned low on the pelvis near the thigh-abdominal junction. See figure 9

and 10.

The upper torso restraint belt must fit over the midpoint of the shoulder and across the chest as illustrated in figure 9 and 10.

The occupant restraint can be installed/routed on the outside of the wheels if the wheelchair has small wheels and the wheels do not prevent the seat belt from making full contact with the hip of the user.

When using a Kudu with wheels up to and including 16", the belt restraint should be routed on the outside of the wheels. See figure 9.

When using a Kudu with wheels larger than 16", the belt should be routed on the inside of the wheels. See figure 10.


Figure 9


Figure 10

Illustration of improper belt restraint fit

Occupant belt restraints must not be held away from the body by wheelchair components such as armrests or wheels. See figure 11.


Figure 11

Manual release in case of electric fault of the docking system:

A manual emergency release lever is located at the front edge of the docking unit.

In order to use the manual release:

- Move the wheelchair forward to remove the pressure on the lock pin and the release lever should be pushed sideways and held in order to release while the wheelchair is being moved away;
- A cable-activated manual operating lever can also be installed (available for MK II only). The red release arm is pushed to one side and should be held there whilst the wheelchair moves away.

If the described manual release procedures fail, a red emergency release tool comes with each docking unit.

In order to release the wheelchair using the emergency release tool:

1. Move the wheelchair forward to remove the pressure on the lock pin.
2. Place the emergency release tool in the gap between the locking plate and the docking station.
3. Push the release tool and wheelchair forward until the locking pin has been forced down - after which the wheelchair can reverse out of the docking station. See figure 12.

Figure 12

03 - Warnings

The product is approved according to ISO 7176-19:2022.

The wheelchair has been dynamically tested in a forward facing orientation with the anthropomorphic test device restrained a three-point belt restraint.

Pelvic- and shoulder-belt restraints shall be used to reduce the possibility of head and chest impacts with vehicle components.

Wheelchair-mounted trays not specifically designed for crash safety should be removed from the wheelchair and secured separately in the vehicle or be secured to the wheelchair but positioned away from the occupant with energy-absorbing padding placed between the tray and the occupant.

Other auxiliary wheelchair equipment should be either secured to the wheelchair or removed from the wheelchair and secured in the vehicle during travel.

Postural supports should not be relied on for occupant restraint in a moving vehicle. Use a vehicle anchored and certified occupant restraint system instead.

The product should be inspected by the dealer before reuse following involvement in any type of vehicle collision.

Alterations or substitutions should not be made to the wheelchair securement points or to structural and frame parts or components. If it is considered necessary to make these kinds of alterations, the wheelchair manufacturer shall be consulted.

The approval is not valid on 'custom made' chairs.

Care should be taken when applying the occupant restraint to position the seatbelt buckle so that the release button will not be contacted by wheelchair components during a crash.

The product shall only be used as indicated in the instructions for use.

The product is not approved for users below a user mass of 22kg.

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www.etac.com/support/support-center/support-documents/

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Note 29:

EN: When the screw heads, which are mounted in the ends of the red locking plungers 501832, are positioned correctly, the spring-loaded red locking mechanisms 501832 must plunge into the sockets of the locking brackets 501833. Hereafter the console has been secured. To release the console carry out the same steps in reversed order, after the red spring-loaded locking plungers 501832 have been unlocked by sliding them inwards towards each other.

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R82 A/S