


Nova Standard low / low net


Manual

A

	XS	S	M	L	XL	XXL
A	43	43	43	44	53	53
B	36	47	53	61	83	93
C	53	53	53	54	54	54
D	76	87	94	108	126	136
E	16	16	18	19	19	19

(DA) Indhold

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General information

Thank you for choosing a sling from Etac.

Etac's slings is tested by accredited and independent testing institutes and comply with the requirements of the Medical Devices Directive Class 1 products (MDD 93/42/EEC).

The slings also comply with the requirements in standard EN ISO 10535.

In this manual, the user is the person being lifted. The assistant is the person operating the sling.

This symbol appears in the manual along with text. It indicates a need for caution when the safety of the user or assistant may be at risk.

Etac continuously makes improvements to our products.

We therefore reserve the right to make changes to our products without prior notice.

The dimensions shown in the drawings or other material are for guidance only.

We cannot be held responsible for any printing errors or omissions.

Warranty: Two years for material and manufacturing defects with correct use.

For further information on Etac's range of slings and lifts, see www.etac.com.

Properties

The sling has a low back and divided leg support with a medium cut.

The sling has a centred band at the back and pockets at the bottom to make it easier to position correctly.

The sling can be used for lifting to

and from a sitting position.

Approved for users up to 350 kg/
770 lbs.

Remember

- In order to avoid injury during handling and use, please read through the lift and sling manual
- It is important to choose the right size. 
- This sling is not recommended for people with an amputated leg(s).
- It is important to follow the instructions on how to perform a transfer.
- Check the sling regularly, particularly after washing. Pay particular attention to wear and damage to seams, material, straps and loops.
- Do not use damaged lifting accessories.
- Never use a defective sling!
- The sling should not be stored in direct sunlight.
- If required, consider using one or more co-assistants from time to time.
- Plan the lifting operation prior to the lift to ensure it takes place as safely and as smoothly as possible.
- It is important to check that the straps are correctly fastened to the slingbar once the straps are fully stretched, but before the user is lifted from the support.
- Never lift a user higher from the ground than is necessary to perform the lift.
- Never leave a user unattended in a lifting situation.
- If the lifting frame which is used has a catch to prevent accidental unhooking, the catch must be in working condition. Lifting frames without a catch or where the catch does not work must not be used.

Practical operation

fig. B

Make sure that the wheels of a wheelchair, bed, bunk, etc. are locked during operation.

1. Fit the sling around the user. Make sure the user's arms are on the outside.
2. Suspend the sling from the hooks on the slingbar.
3. Tension the straps of the sling without lifting the user. Check that the sling is correctly suspended from the hooks and is correctly fitted around the user.
4. Lift the user, taking care not to lift him/her any higher than is necessary.

Sources of errors/User sitting incorrectly

- Is the lifting sling the right size?
The sling is too small if the leg supports do not reach at least halfway down the inside of the thigh and the slingbar is too close to the face.
- Is the sling positioned correctly?
The sling may not have been pulled far enough down the back when being fitted round the user..
- Are all the straps lifting simultaneously? If the sling straps lift unevenly, the sling may slide up the user's back, causing the user to slip down into the sling hole.
- Ensure the leg supports are correctly positioned on the thigh.
- Is the sling connected correctly according to individual prescription?
- Do separately connected leg supports increase the risk of the user slipping out?
- If the user is wearing clothing made of slippery material, does he/she slide down/out of the sling more easily?

EN

Cleaning

Washable up to 90 °C/194 F

Tumble-dry up to 45 °C/113 F.

NB!

Do not bleach, autoclave or dry-clean.

The sling must be checked on a regular basis, preferably before each use but especially carefully after washing.

Check there is no wear and damage to seams, material, straps and loops.

NB! Never repair a sling yourself. A worn sling must never be used.

EN

Accessories

Extension strap

For use in situations where a particular seated position, sitting angle or lifting height is required.

Available in lengths 14 cm/6 Inches or 18 cm/7 Inches, sold separately in bags.

Headrest

Helps to keep head in place in low harnesses and provides extra support in high harnesses. Attached in the lift harness.

Marking straps

Marking straps for marking tested loop positions on the sling for lifting.

Two supplied with every sling in plastic bags. Can also be purchased as an accessory.

Extra padding

Intended for users that are extra sensitive to pressure or require extra comfort. Manufactured in soft and pleasant sheepskin imitation polyester. Sold in pairs.

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Etac Sverige AB
Etac A/S