User manual

LeanOnMe Wedge
LeanOnMe Triangle

immedia
by Etac

IMM1070_IFU - 2025-10-08 – Rev. 1.0

www.etac.com

Symbols

4

Device overview

8

General

10

English	10
Svenska	14
Norsk	18
Dansk	22
Suomi	26
Deutsch	30
Nederlands	34
Français	38
Italiano	42
Español	46

Application suggestions

50

LeanOnMe Wedge and LeanOnMe Triangle	51
--------------------------------------	----

Warnings	53
----------	----

Information	55
-------------	----

EN: Warning: This symbol appears in the user manual along with a number referring to the instructions below. It draws attention to situations where the product or the user or carer's safety may be in danger.

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EN: Information: This symbol appears in the user manual along with a number referring to supplementary information which can be found in the information section.

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Example: LeanOnMe Wedge with Non-Slip Cover

Product label

5

Example: LeanOnMe Triangle with Soft Touch Cover

Product label

5

1

EN: Handle
SV: Handtag
NO: Håndtak
DA: Håndtag
FI: Kava
DE: Handgriff
NL: Handgreep
FR: Poignée
IT: Impugnatura
ES: Asa

2

EN: Label inside the cover
SV: Etikett på insidan av överdraget
NO: Etikett på innsiden av trekket
DA: Mærkat på indersiden af betrækket
FI: Tarra kannen sisäpuolella
DE: Etikett auf der Innenseite des Bezugs
NL: Etiket aan de binnenkant van de hoes
FR: Étiquette à l'intérieur de la housse
IT: Etichetta all'interno del rivestimento
ES: Etiqueta en el interior de la funda

3

EN: Zipper
SV: Blixtlås
NO: Glidelås
DA: Lynlås
FI: Vetoketju
DE: Reißverschluss
NL: Ritsluiting
FR: Fermeture à glissière
IT: Cerniera
ES: Cremallera

4

EN: Non-Slip on the bottom
SV: Halksäker undersida
NO: Sklisikker underside
DA: Non-Slip i bunden
FI: Liukumaton pohja
DE: Rutschfeste Unterseite
NL: Antislip onderaan
FR: Antidérapant sur le dessous
IT: Antiscivolo sulla parte inferiore
ES: Parte inferior antideslizante

5

EN: QR-code
SV: QR-kod
NO: QR-kode
DA: QR-kode
FI: QR-koodi
DE: QR-Code
NL: QR-code
FR: Code QR
IT: Codice QR
ES: Código QR

General

Thank you for choosing an Etac device.

This manual must be read thoroughly before you start using the device as it contains important safety information and gives instructions on how to use the device correctly. Visit www.etac.com for all relevant product documentation, for example the user manual, maintenance and recycling instructions. You can select your language via the “International” and “Local websites” link.

You can also consult the QR code on the device for further information and instructions.

Copyright

This document is the property of Etac Immedia A/S. The content shall not be modified in whole or part without prior approval of Etac Immedia A/S.

Disclaimer

Etac continuously makes improvements to our devices and information. We therefore reserve the right to modify our devices and instructions for use without prior notice. Consult www.etac.com to find the latest version of the documentation.

CE Declaration of Conformity

The device mentioned in this manual and the different models hereof comply with the Medical Device Regulation (EU) 2017/745.

The device is classified as a Class I medical device.

Product standards

The device has been tested and complies with the requirements applicable for this type of device in the below mentioned standards:

- ISO 21856:2022 Assistive products — General requirements and test methods.

- DS/EN ISO 20342-1:2019 Assistive products for tissue integrity when lying down - Part 1: General Requirements.

- The devices comply with EN 1021-1:2014 and EN ISO 12952-1:2010 resistance to ignition by cigarette (but not flame).

Warranty

Please refer to www.etac.com/support.

If the customer makes adjustments, repairs or uses combinations not predetermined by Etac, the Etac CE certification and Etac warranty will not be valid. If in doubt, please contact Etac.

Notice

Any serious incident that has occurred in relation to the device should be reported to your distributor/representative and the national competent authority within your region in a timely manner. The distributor/representative will forward the information to the manufacturer.

Pre-sales information

Please visit www.etac.com.

Device description

Etac supplies a wide selection of Positioning Cushions for different types of positioning situations.

LeanOnMe Wedges are designed to provide a stable support in supine, lateral or sitting position, and are available in three sizes with either Soft-Touch or Non-Slip cover.

LeanOnMe Triangles are designed to provide support for smaller parts of the body and fill gaps in both supine or sitting position, and are available in two sizes with either Soft-Touch or Non-Slip cover.

See www.etac.com for a complete overview of possible variants.

Clinical benefits

The clinical benefits of the devices are:

- Can reduce risk of pressure injury
- Can relieve pain perception
- May reduce risk of joint stiffness
- Can relieve edema in upper and lower extremities by elevation
- Can facilitate respiration by supported positioning
- May enhance well-being and relaxation by supported positioning
- Can reduce physical strain on carer

Materials

Filling: Polyurethane foam.

SoftTouch Cover:

Outside: 75% Lyocell (Tencel™ + Smartcel™ sensitive), 25% Modal (Tencel™). Inside: 100% polyester filling and polyurethane membrane.

Bottom: 100% polyester with polyurethane membrane.

Non-Slip Cover:

100% polyester with polyurethane membrane.

Introduction

Definitions

The person referred to as “the user” is the person who is positioned with the device. “The carer” is the person who manages the device.

Intended purpose

The Positioning Cushions are assistive devices intended for alleviation of, or compensation for, a functional impairment due to an injury or disability. The devices are designed to provide support and relieve in multiple positions for an individual lacking the ability to position themselves due to reduced mobility or physical strength.

Intended user

The intended users of the devices are based on the individual's functional ability and not on a specific diagnosis, health condition or age. Ensure a clinical risk assessment is conducted. Secondary users of the devices are carers providing assistance and clinicians/technicians who set up the device.

Intended environment

The devices are intended for indoor use.

Intended application

The device is intended for ‘Long term’ use, meaning intended for continuous use for more than 30 days. The device is intended for use with intact skin only. The device is intended for refurbishment and re-use.

Expected service life

The expected service life of positioning cushions are 2 years. The lifetime of the device may vary depending on usage frequency, loads and washing/drying frequency and conditions.

Indications

Indication for use is disability, including, but not limited to, users who are in need of support to be safely and comfortably positioned in lying or sitting positions, e.g. individuals at risk of developing pressure ulcers, including those diagnosed with stroke, cancer, diabetes mellitus, urinary tract infections, pneumonia, Huntington's disease, critically ill patients, Obstructive Sleep Apnea (OSA), particularly positional OSA (POSA). Not an exhaustive list. Furthermore, these devices are prescribed for relieving a number of symptoms such as back pain, forward head posture, hypertonia, contractures, edema and for use in bariatric care. Not an exhaustive list.

Contraindications

There are no known contraindications.

Precautions

Important that the cushion do not block the airways.

Training requirements

Only carers and/or users who have read and understood the user manual are allowed to use the device.

Safety information

The safety of both user and carer must never be jeopardized – if you suspect that the device has a malfunction you must put it in quarantine marked “out of order” and call your local distributor/representative.

Safety signs, safety measures and warnings symbols are given in the context relevant for the actual handling and/or operating situation – it is mandatory the given instruction is followed to avoid dangerous/hazardous situations.

Signs, symbols and instructions placed on this device must never be covered or removed and must remain present and clearly legible throughout the entire lifetime of the device. Immediately replace or repair illegible or damaged signs, symbols, and instructions. Contact your distributor/representative for instructions.

Risk assessment

Do a risk assessment of the LeanOnMe device, the individual user and the positioning, and document your analysis. You as the carer are responsible for the safety of the user.

Plan the positioning in advance to ensure that it is performed as safe and smooth as possible.

Remember to work ergonomically according to local recommendations.

Encourage the user to assist where possible.

The device

Warnings:

- Never use a product without cover.
- Always use the right cover.
- Always position the zipper away from the body.
- Check regularly to ensure that the device(s) are properly placed.
- Never leave the user alone on the edge of the bed.
- Never leave the product on the floor.
- Clean and/or disinfect the device between users to avoid cross contamination.
- Always transport the device away from the floor.
- Never use a defective product. If the product shows signs of wear-and-tear, it must be scrapped.

Safety measures

Important that the cushion do not block the airways

Materials

All used materials which are to be in direct contact with skin are biologically assessed and evaluated according to relevant standards. However, users/carers who are allergic to any of the materials are recommend not to use the device.

Spare parts

Article numbers and spare part can be found on www.etac.com

Combinations

When combining the device with other devices or accessories, the combination must not change the intended purpose of the devices or modify the devices in such a way that compliance with the applicable requirements may be affected. The person/company responsible for the combination of the devices must ensure that the applicable requirements are fulfilled.

Mounting/assembly

Apply or remove the Soft Touch/Non-Slip cover:

Open the zipper on the cover and insert or remove the foam. Close the zipper after applying the cover.

Service, cleaning and maintenance

Maintenance and service

Maintenance and inspections of the device are wholly the responsibility of the owner of the device. Failure to maintain a device in accordance with the instructions may compromise the safety of the user and/or carer and may affect the expected lifetime of the device.

Make sure:

- that the materials, seams, stitching, zippers, handles are intact/faultless.
- never to repair a product yourself.

Cleaning

80 °
176F

It is recommended to wash the cover before use. Remove the foam before washing.

Do not over dry.

Do not drain dry.

The Non-Slip Cover can be wiped off with detergent with pH-level between 5-9 or 70% disinfection solution or disinfectant solution or commercial and/or household cleaning agents.

Refurbishment and reuse

The LeanOnMe Wedge and LeanOnMe Triangle and covers are suitable for reuse after inspection and cleaning.

Periodic inspection

The device should be checked regularly, preferably each time it is used, and especially after washing.

Transportation, storage and disposal

Transportation and storage

When the device is not in use, it can be stored under the following conditions:

Keep packaging away from children

Keep away from fire

The device can be stored and transported under temperatures between -30 to +40 °C.

Disposal

The device contains no hazardous substances that require special disposal. Check if there are options available for recycling in an environmentally friendly way and follow the national provisions that apply. Dispose of the device according to the national provisions that apply.

Recycling

See the recycling instructions on www.etac.com, and follow the national provisions that apply.

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EN: Information

01. Application suggestion, the user's special needs and risk assessment define the final solution.
02. LeanOnMe Wedge and LeanOnMe Triangle can beneficially be used together with LeanOnMe Cushions.
03. Firm and stable positioning is required, either in a lateral or supine position.
04. The LeanOnMe Triangle supports smaller parts of the body in supine, lateral, or seated positions. They can be used to fill gaps where additional support is needed.
05. In lateral positioning, the wedges can be placed under the knees or at the foot of the bed to support the feet.
06. When transitioning from a supine to a seated position, the wedges provide effective support for the back, or as lateral support when leaning to the side or standing up.

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**Etac Ltd.**