

Product	Task								Uniqueness								SWL	Levels of Mobility				
	Turning in bed	Repositioning lying	Transfer between two horizontal surfaces	Sitting Positioning	Transfer between two sitting positions	In and out of a car	Support stand up/sit down	Support standing transfer	Single handed care	For turning users with a hoist	Early rehabilitation	Single Patient Use (SPU)	Minimize friction under pressure points	Stays under user	For incontinence users	Application of sling	Dressing	Safe Working Load (KG)	A - Independent	B - Dependent in some situations	C - Dependent in many situations	D - Dependent in most situations
<b>Positioning Horizontal - Systems and solutions for horizontal positioning, turning and moving higher up I bed</b>																						
SatinSheet BaseSheet (without drawsheet)	.	.						.		.		.	.				200	.	.			
SatinSheet BaseSheet with 2Direction/4Direction DrawSheet	.	.						.	.	.		.	.				200	.	.	.	.	.
SatinSheet BaseSheet with In2Sheet 4Direction	.	.						.	.	.		.	.	.			200	.	.	.	.	.
4WayGlide System	.	.						.	.	.		.	.				200/300	.	.	.	.	.
TwinSheet4Glide System	.	.						.	.	.		.	.	.			200	.	.	.	.	.
MultiGlide	.	.		.				.		.	.	.	.		.	.	300	.	.	.	.	.
In2Sheet	.	.						.		.		.	.	.			200	.	.	.		
GlideCushion	.	.		.				.		.		.	.				200	.	.	.	.	.
TurningMattress	.	.	.					.		.		.	.				200	.	.	.	.	.
<b>Positioning Sitting - Solutions for positioning when sitting</b>																						
AutoTurn				.		.						.	.				150	.	.	.		
EasyTurn				.		.				.		.	.				150	.	.	.		
OneWayGlide				.				.				.	.				150	.	.	.	.	.
<b>Transfer Sitting - Solutions for seated transfers between bed, wheelchair, chair , toilet, etc.</b>																						
PediTurn hard and soft				.								.					150	.	.	.	.	
E-Boards and 3B-Board, Butterfly, Dolphin				.		.		.		.		.					150/250	.	.	.	.	
E-board Oval			.	.				.		.		.			.		150	.	.	.	.	.
Swan	.	.		.		.		.		.		.					150	.	.	.	.	.
<b>Transfer Horizontal - Solutions for horizontal transfers between beds, bed to shower trolley, etc.</b>																						
2Move			.								.	.					150-300			.	.	.
TransferMattress	.	.	.									.					200			.	.	.
<b>Support - Solutions for support, mobilisation and transfer</b>																						
SupportBelt				.		.	.	.		.							150		.	.		
Sling /SmartMove		.		.	.	.		.		.							150		.	.	.	.
BedString		.		.				.		.							150		.	.	.	.
<b>Special Solutions - Solutions for support, mobilisation and transfer</b>																						
SlingOn												.			.		300			.	.	.
MiniBoard			.														300			.	.	.
OneManSling		.						.		.							150		.	.	.	.
LeanOnMe/Wedge		.		.				.		.			.				200	.	.	.	.	.
MultiGlide Glove	.	.		.				.		.		.			.	.	200		.	.	.	.
<b>Sit-to-Stand</b>																						
Etac Turner Pro				.	.		.	.		.							200		.	.		
Molift Raiser Pro				.	.		.	.		.							170		.	.		