


REF 82488-1
82488-2

etac[®]

Mjukdel lårstöd, höger/vänster

Soft pad thigh support, right/left

Oberschenkelpolster, rechts/links



75068.2 2022-08-15

Teknisk data / Technical data / Technische Daten

	<p>Användningstemperatur Temperature for usage Verwendungstemperatur</p>	<p>Lagring 6 månader/ max 25° /mörker Storage 6 months/ max 25° / darkness Lagerung 6 Monate/ Max 25°/ im Dunkel</p>
--	--	--

Monteringsanvisning / Assembly instructions / Montageanweisung

<p>1.</p>	<p>3.</p>	<p>Kardborre klistras i rumstemperatur Glue Velcro in room temperature Klettband im Raumtemperatur zusammenkleben</p>
<p>2.</p>	<p>alkali ph min 6</p> <p>Rengör yta före applicering Clean the surface before application Die Oberfläche reinigen vor dem Auftragen</p>	
<p>4.</p>		
<p>Nyklistrad kardborre belastas först efter 2 timmar Glue joint to be charged first after 2 hours Leimfuge erst nach 2 Stunden beladen</p>		