



Molift RgoSlings

Excellent ergonomics for
comfortable hoisting

 **etac**[®]

Molift RgoSlings

Molift has a long reputation for manufacturing high quality hoists and slings. RgoSling is a portfolio of slings developed to meet the demands for the everyday needs in healthcare.

Maximum Comfort

All our slings are designed to distribute weight and pressure evenly for maximum comfort. As Molift is committed to comfort and safety we offer a wide range of various models and sizes. We recommend individual testing. We have a colour code indicating sling size to facilitate a quick and easy selection.

Anatomical shape

Molift RgoSling has an anatomical shape and a fitted cut with fixed padding both in the back and leg parts. The anatomic design allows for the natural curvature of the spine and stops the user from sliding. In combination with the balanced weight distribution, a safe and comfortable hoist is achieved.


Unique sliding loop function

The unique sliding loop function in the shoulder straps, balances the weight distribution in every stage of the hoist, making it very comfortable to hoist from a lying to a sitting position and vice versa. Thereby an external angle adjustment of the sling bar is redundant.

Scandinavian Design

Molift slings are designed in Norway, applying the same high standards as it does to Molift hoists. All the slings have been approved according to quality standard EN ISO 10535. We recommend the use of a 4-point sling bar for comfortable hoisting, but the slings may also be used with a 2-point sling bar.





Sliding loops

Sliding loops on head straps. Making a smooth transfer from lying to sitting position. Ensures a comfortable and safe position for the patient.

Colour coding

Colour coding of hoisting loops helps you find the right fit.

Comfortable padding

Padding in back and leg support distributes the weight over a larger area, making the sling more comfortable for the user.

Labelling

Clear, practical information on a highly visible label. Washing instructions and technical information are found in a protective pocket at the back.

Guide straps

Guide straps assist with positioning during the lifting process.

Application pocket

Pockets in the back and on the leg supports make application easier.

Standard models: HighBack and MediumBack

Our standard models are designed for use in the usual transfer situations and are suitable for most users. The standard models are available in two back heights: HighBack and MediumBack.

Models with specific features

We know that users have many different requirements, so we've designed additional models with specific features to meet their needs, including the RgoSling Toilet and RgoSling Amputee.



Molift RgoSling HighBack















Standard leg support opening

- Sliding loops on head straps
- Support for the entire back and the head
- Hoisting to and from a lying or a sitting position
- All-round sling that fits most users
- Hoisting from floor

Molift RgoSling MediumBack

Standard leg support opening

- Sliding loops on head straps
- Suitable for a user with headcontrol
- All-round sling that fits most users
- Adjust the sitting angle by the leg straps

Size	HighBack padded	HighBack net	Size	MediumBack padded	MediumBack net
XXS 	1720200	1720400	XXS 	1720100	1720300
XS 	1720210	1720410	XS 	1720110	1720310
S 	1720220	1720420	S 	1720120	1720320
M 	1720230	1720430	M 	1720130	1720330
L 	1720240	1720440	L 	1720140	1720340
XL 	1720250	1720450	XL 	1720150	1720350
XXL 	1720260	1720460	XXL 	1720160	1720360

Safe working load (SWL): 300 kg
HighBack padded XL - XXL 350 kg

Safe working load (SWL): 300 kg
MediumBack padded XL - XXL 350 kg



Molift RgoSling Toilet








Wide leg support opening

- Fixed loops
- Easy access sling, facilitate dressing and undressing
- For users with good core stability and muscle tonus
- Comfortable padding
- Flexible leg support
- Belt included - for extra security








Molift RgoSling Amputee

Narrow leg support opening

- Sliding loops on head straps
- For users with an amputated leg or a user with reduced core stability or muscle tonus
- Ideal solution for user with a high amputee
- Comfortable padding
- Belt included - for extra security

Size	LowBack	HighBack
XXS 	1720500	1720600
XS 	1720510	1720610
S 	1720520	1720620
M 	1720530	1720630
L 	1720540	1720640
XL 	1720550	1721450
XXL 	1720560	1720660

Safe working load (SWL): 300 kg

Size	MediumBack	HighBack
XXS 	1721100	1721200
XS 	1721110	1721210
S 	1721120	1721220
M 	1721130	1721230
L 	1721140	1721240
XL 	1721150	1721250
XXL 	1721160	1721260

Safe working load (SWL): 300 kg



Molift RgoSling Shadow













Standard leg support opening

- Sliding loops on head straps
- Suitable for users in moulded seats
- No straps, seams, pockets or labes in contact with the user
- Polyester net with thin webbing: no straps, seams, pockets or labels in the contact surface to the user


Molift RgoSling Comfort HighBack

Full body support sling without split leg

- Full body and head support
- Can be left in the wheelchair underneath the user
- Polyester net which provides good ventilation
- Red seam indicating the center line of the back, to facilitate correct application

Size	HighBack	Size	HighBack
XXS 	1722900	XS 	1721410
XS 	1722910	S 	1721420
S 	1722920	M 	1721430
M 	1722930	L 	1721440
L 	1722940	XL 	1721450
XL 	1722950	XXL 	1721460

Safe working load (SWL): 300 kg



Molift RgoSling Ambulating Vest

For gait training with assistance

- 4 straps with fixed loops
- Self tightening strap in the waist
- Standing or balance exercise
- Double buckles for extra security
- Standard polyester fabric with padding

Molift RgoSling Fabric Stretcher



A soft fabric hoisting device for horizontal transfers

- Molift RgoSling Fabric Stretcher is suitable when hoisting in a normal sling would be inappropriate due to medical conditions
- Application easily performed
- Standard polyester fabric with padding of the head part
- Used with an 8-point sling bar

Size	Item no	Size	Width x Length	Item no
XXS	1721300	M	450 x 1750 mm	1721520
XS	1721310	L	650 x 1780 mm	1721530
S	1721320	XL	900 x 1830 mm	1721540
M	1721330			
L	1721340			
XL	1721350			
XXS-XS	Groin straps, single pack (order two for a pair)			1721700
S-M	Groin straps, single pack (order two for a pair)			1721720
L-XL	Groin straps, single pack (order two for a pair)			1721740

Safe working load (SWL): 300 kg





Molift RgoSling MediumBack Plus

Full body support for users weighing from 160 kg

- Provides full body support
- Straps with fixed loops and multiple height options
- Comfortable padded leg-support with dual straps
- Colour coded straps
- Split leg opening
- Can be used with a separate headrest

Molift RgoSling Repositioning sheet


An in-bed lift and positioning sheet

- Turning and transferring patients higher in the bed
- Lifting from bed to another supine surface, stretcher etc
- Stays in bed - for comfortable and dignified care
- Facilitates care situations - saves time for caregivers
- Colour coded straps facilitate communication when attaching the straps equally on both sides

Size	MediumBack Plus	Size	Quantity	Repositioning sheet
L ■	1721940	L 1340 x W 1300 mm	1 pc	1721570
XL ■	1721950	L 1340 x W 1300 mm	10 pcs	1721570-C
Safe working load (SWL): 500 kg		Safe working load (SWL): 500 kg		


Facts

Materials



Polyester

The standard polyester fabric is durable and hard wearing. The material is easy to maintain and has a dirt repellent surface.




Polyester net

Hard wearing and durable soft net fabric, breathes well and dries quickly. Suitable for shower and bathing situations.



Non-woven

Polypropylene non-woven. Patient specific slings when the sling is for short term use.



Durable polyester

Extra strong durable polyester fabric for products that has a weight limit up to 500 kg. The material is easy to maintain and dries quickly.

Detailed information

Size	Recommended user weight
XXS	12 - 17kg
XS	17 - 25 kg
S	25 - 50 kg
M	45 - 95 kg
L	90 - 160 kg
XL	160 - 240 kg
XXL	230 - 300 kg

Washing instructions



Washing temperature 60-85 °C / 140-185 °F
 Washing at higher temperatures wears out the material faster.



Accessories




Molift RgoSling Extension Loops

Straps for easier adjustment
RgoSling: Used to extend the hoisting straps, flexible for 3 different length.
N:B: Order 2 for a pair.



Immedia SlingOn

When the user can't be turned onto their side Immedia SlingOn used with a folding technique makes it possible to apply the sling without first turning the user over to the side.
Washing in 80 °C or 176 °F.



Molift Strap Marker

For marking a selected hoisting loop for head or leg support hoisting straps. Mark the hoisting loops selected during individual testing of the sling.

Product	Extension	item no.	Product (sold in pairs)	Size	item no.
Molift RgoSling Extension loops (flexible)	+60, 180 or 260 mm	1721600	Immedia SlingOn - Small	800 x 1000 mm	IM437/S
			Immedia SlingOn - Medium	800 x 1700 mm	IM438
			Immedia SlingOn - Large	1200 x 1500 mm	IM437
Product (sold in pairs)		item no.			
Molift Strap Marker		1723600			





Etac is a world-leading developer and provider of ergonomic assistive devices and patient handling equipment. Our heart lies in the solutions that optimise quality of life for the individual, their family and caregivers.

For the latest news and continuously updated product information – visit: www.etac.com