



molift[®]

Molift UnoSling

Safe and comfortable disposable slings for single-patient use

Molift UnoSling

Molift UnoSling is a portfolio of single-patient use slings (SPU), developed to meet the growing need in healthcare institutions, where infection risk prevention is of utmost importance.

Benefits of SPU products

SPU products are essential for preventing the transmission of infection. They can also enable cost savings and a more efficient operation for the facility by avoiding logistics and laundry costs. It also leads to easier storage of products that are always ready for use.

High infection control requirements

As Molift UnoSlings stays with the patient throughout their stay at the facility and are disposed when they become dirty or damaged or when the patient leaves the facility they are an effective solution for preventing the transmission of infection.

Strong and secure material

Molift UnoSlings are manufactured in a strong non-woven material which is non-washable. The sling is discarded if it becomes dirty or damaged, or if the user no longer needs it.

Models with specific features

Although being disposable products, Molift UnoSlings are still designed with high quality and patient comfort and safety in mind. Several models are available to meet the different requirements of users, UnoSling HighBack, UnoSling Toilet, UnoSling Repositioning Sheet, UnoSling Ambulating Vest and UnoSling LimbLift.





Non-wash label

The symbol will change appearance if the product is accidentally washed.



2-point or 4-point


Suitable for use with either 2-point or 4-point sling bar.

Disposable for single-patient use

To prevent the transmission of infections, the sling is disposed of when soiled, or the patient no longer needs it.

Label with user's name

All products have a large, clear label with practical information. The size is shown clearly and it is easy to add the user's name and date of first use, which ensures safe use.




Materials



Non-woven

Polypropylene non-woven. Applies to all Molift UnoSlings except Molift UnoSling Repositioning Sheet which is made of polyester.



Polyester

The standard polyester fabric is durable and hard wearing. The material is easy to maintain and has a dirt repellent surface.

Molift UnoSling HighBack

Suitable for most users in common hoisting situations.


The Molift UnoSling HighBack is suitable for most users in common hoisting situations, such as between bed and wheelchair. The sling remains with the user during the entire length of their stay in the hospital.

Molift UnoSling HighBack is suitable for use with either 2-point or 4-point sling bar.

- Fits most patients in common hoisting situations
- Hoisting to and from lying or sitting position
- Support for the entire back and head
- Suitable with 2- or 4-point sling bar



Suitable with 2- or 4-point sling bar

Molift UnoSling HighBack				Measurements		
						
Size	Pieces/pack	Item number	Size	(A) Back height	(B) Width	
S	20	3060020	S	80 cm	35 cm	
M	20	3060030	M	93 cm	40 cm	
L	20	3060040	L	100 cm	47 cm	
XL	20	3060050	XL	105 cm	54 cm	
XXL	10	3060050	XXL	110 cm	60 cm	

Maximum user weight 350 kg

Molift UnoSling Toilet

For users with good core stability and muscle tone

The Molift UnoSling Toilet is a low back sling with a wide leg opening to facilitate dressing and undressing when using the toilet. The sling is easy to apply and remove and has a belt for added security during hoisting.

The inside of the sling has a non-slip lining to help the sling stay in place.

The sling has a green seam in the back, as a guide to help that the sling is placed in a centered position.

- For patients with good core stability and muscle tone
- Facilitates dressing and undressing
- Easy to apply and remove
- Safety belt for extra security during hoisting
- Suitable with 2- or 4-point sling bar



Safety belt for extra security during hoisting

Molift UnoSling Toilet				Measurements		
Size		Pieces/pack	Item number	Size	(A) Back height	(C) Recommended waist size
S	■	5	3090020	S	33 cm	76-100 cm
M	■	5	3090030	M	33 cm	85-110 cm
L	■	5	3090040	L	33 cm	94-120 cm
XL	■	5	3090050	XL	34 cm	115-140 cm

Maximum user weight 350 kg



Molift UnoSling Ambulating Vest

Sling for safe rehabilitation activities, such as gait training, walking and balance training

The Molift UnoSling Ambulating Vest is designed for patients in need of safe and secure rehabilitation such as gait training, standing as well as balance exercise. It comes with double buckles and a non-slip material inside the belt. It also includes optional groin straps for preventing the vest from sliding up and out of position.


The sling has a green seam in the back, as a guide to help that the sling is placed in a centered position.

- Ensures a safe rehabilitation and early mobilisation
- Double buckles and non-slip lining for extra security
- Groin straps to prevent the vest from sliding up
- Suitable with a 2-point sling bar




Groin straps for extra security

Molift UnoSling Ambulating Vest



Measurements



Size	Pieces/pack	Item number
S ■	5	3070020
M ■	5	3070030
L ■	5	3070040
XL ■	5	3070050

Maximum user weight 350 kg

Size	(A) Back height	(C) Recommended waist size
S	20 cm	72-92 cm
M	20 cm	90-110 cm
L	21 cm	105-125 cm
XL	21 cm	123-145 cm



Molift UnoSling StandUp

Disposable sling suitable for active sit-to-stand transfers

The Molift UnoSling StandUp is suitable for users that require the support of a sit-to-stand device to transfer from one seated position to another, such as to or from a bed, a wheelchair or a toilet. The sling distributes weight over the back and provides comfortable support for the user through the entire sit-to-stand movement. It is also a useful product in early rehabilitation activities with an active sit-to-stand hoist.

The inside of the sling has a non-slip lining to help the sling stay in place.


Molift UnoSling StandUp is designed for use with sit-to-stand hoists with a 2-point loop-style sling bar.

- Useful for active sit-to-stand transfers
- The belt and non-slip lining provides extra safety
- Distributes the weight over the back, for comfortable support




Suitable with a 2-point loop-style sling bar

Molift UnoSling StandUp



Measurements



Size	Pieces/pack	Item number
S	5	3080020
M	5	3080030
L	5	3080040
XL	5	3080050

Maximum user weight 350 kg

Size	(A) Back height	(C) Recommended waist size
S	34,5 cm	70-95 cm
M	35 cm	90-110 cm
L	35 cm	100-120 cm
XL	35 cm	110-140 cm

Molift UnoSling LimbLift

Single patient use sling for lift and support of limbs

For lift and support of limbs

The Molift UnoSling LimbLift is a lifting band designed for different uses such as lifting and holding a body part. It may be used in situations such as foot care or examination and treatment of pressure ulcers. It is available in two widths.

The sling remains with the user during the entire length of their stay in the hospital.

- Used during treatment, washing or examination
- For a leg, arm or other parts of the body
- Reduces the strain on caregivers
- Pocket on one side to facilitate application




Molift UnoSling LimbLift



Size	Pieces/pack	Item number
Standard	10	3100020
Wide	10	3100030

Maximum user weight 200 kg

Measurements



Size	(A) Length	(B) Width
Standard	95 cm	15 cm
Wide	95 cm	25 cm

Molift UnoSling Repositioning Sheet

Suitable for repositioning in-bed and transfers between supine surfaces

The Molift UnoSling Repositioning Sheet is suitable for a multitude of care situations to make transfers and repositioning tasks faster, easier and safer for both caregiver and user. For instance turning bed-bound patients, lifting patients higher in the bed and for bariatric care.


The sling remains with the user during the entire length of their stay in the hospital.

- Can be used for a variety of care situations for bedridden patients
- For turning and moving higher in bed
- Can remain in bed under the patient
- Suitable for bariatric care




Clear labeling for up/down and top/bottom

Molift UnoSling Repositioning Sheet



Measurements



Size	Pieces/pack	Item number
One size Maximum user weight 500 kg	10	3200020

Size	(A) Length	(B) Width
One size	215 cm	135 cm



Etac is a world-leading developer of innovative assistive devices and patient handling equipment. Since 1973, we have been committed to improving quality of life for the individual, family members and caregivers.

For the products' intended purpose and latest information, please visit www.etac.com

Etac AB

 +45 79 68 58 33

 patienthandling@etac.com

 www.etac.com

 **etac**[®]
Creating Possibilities