

Etac Crissy prescription form

1A15745.6 2021-11-04

1 Model						
Swing away			Active			Item. no.
37.5 cm		45 cm	37.5 cm		45 cm	
40 cm		47.5 cm	40 cm		47.5 cm	
42.5 cm		50 cm	42.5 cm		50 cm	

2 Colour, frame	
White -74	
Black -01	
Ruby red -75	
Champagne -76	
Menthol -77	

3 Colour, rear wheel attachment	
Black -01	
Ruby red -75	
Champagne -76	
Menthol -77	

4 Back support system	Back support height		Back support angle		Cover	
	Standard	High	Fixed	Adjustable	Plush	Airmesh
Back support 2B						

5 Push handle/Bow handle	
Push handles, fixed	
Push handles, adjustable in height	
Bow handle, adjustable in height	

7 Attendant brake	
Without attendant brake	
With attendant brake	

6 User brakes			
Model	H	V	
Standard angled red			
Standard angled black			
With extension arm (foldable)			
With extension arm (fixed)			
Scissor brake			
Single-hand brake right or left mounted			
Reset brake alternatives			

8 Leg support		
Model	R	L
Standard, 70° angle		
Standard, 70° angle with locking device		
Narrow, 90° angle		
Narrow, 90° angle with locking device		

Elevating, foldable calf support		
----------------------------------	--	--

9 Other leg support		
Amputee support, plush cover		
Amputee support, hygiene cover		
Amputee weight		
Amputee support+weight, plush cover		
Amputee support+weight, hygiene cover		

10 Foot support		
Model	R	L
Depth 15 cm, medium attachment tube		
Depth 15 cm, long attachment tube		
Depth 13 cm, medium attachment tube		
Depth 13 cm, long attachment tube		
One-piece foot support, foldable and adjustable in angle		

Seat depth, seat height and balance position

1 Seat depth	
S (41 cm*)	
M (43 cm*)	
L (45 cm*)	

*Functional seat depth 3-6 cm longer depending on the setting of the back support straps/back support cover

2 Camber angle		3 Rear seat height and rear wheel			
0°	2°	Rear wheel	22"	24"	24" Light weight wheel (Carbon)
		F	38,5 cm	41 cm	41 cm
		E	40,5 cm	43 cm	43 cm
		D	42,5 cm	45 cm	45 cm
		C	44,5 cm	47 cm	47 cm
		B	46,5 cm	49 cm	49 cm
		A	48,5 cm	51 cm	51 cm

4 Balance position

F
E
D
C
B
A

5 4 3 2 1

7 Handrim	
Aluminium	
Stainless steel	
Cellular rubber	
Plastic coated	
Titanium	

8 Tyre rear wheel	
Schwalbe Rightrun	
Schwalbe Marathon plus	
Pneumatic 1 3/8"	
Solid	
Tyre with solid inner tube	

9 Seat extender	
Without seat extender	
With seat extender	

5 Front seat height, front fork and castor

Castor	Front fork short	Front fork Medium	6 Front fork attachment	
			Standard	Widened
5"	(40.5-47.5) cm	(45.5-51.5) cm		
6"	(44-49) cm	(44.5-53) cm		
8"	-	(49-55) cm		

For correct setting of the front fork attachment: Max. backward inclination 6 cm higher front than rear (9°), max. forward inclination 2 cm lower front than rear (3°). Measures in brackets are theoretical. Hover over the editable text fields to see possible front seat height.

Accessories

Arm supports/Mudguards		R	L	Item. no.
Solid moulded plate, height adjustable	Length 25 height 23-32 cm			25172-1 25172-2
	Length 25 <small>finger screw, height 23-32 cm</small>			26379-1 26379-2
	Length 38 height 23-32 cm			25173-1 25173-2
	Length 38 <small>finger screw, height 23-32 cm</small>			26380-1 26380-2
	Length 25 height 20-25 cm			28578-1 28578-2
	Length 38 height 20-25 cm			28579-1 28579-2
Mudguard, carbon fibre				
Without arm support/mudguard				
Arm support cover				
Arm support cover dark grey plush, detachable	Foam	25 cm		62282-60
		38 cm		62283-60
	Gel	25 cm		26240-60
		38 cm		26241-60
Arm support cushion Hemi				
Without arm support cover/cushion				

		R	L	Item. no.
Tray, to be mounted on long arm supports				
Tray Hemi	For short arm support			
	For long arm support			
Non-slip cover for tray Hemi				81633
Tray Hemi communication right/left				27352-1 27352-2
Without tray				
Seat cushion, plush cover				
Seat cushion, plush cover				
Comfort wedge				
Comfort wedge				
Seat cushion Tidy, hygiene cover				
Seat cushion Tidy, hygiene cover				
Utan sittdyna				
Calf support cover, for elevating leg support				
Calf support cover, for elevating leg support	Plush			62762-60
	Hygiene			62762-63
Without calf support cover				
Calf strap				
Calf strap				
Padding for calf strap				
Padding for calf strap				
20 cm				26845
30 cm				26846
40 cm				26847
Calf support wide				
Calf support wide				
Without calf strap				
Heel strap				
Heel strap				24458
Without heel strap				
Positioning belt				
Positioning belt			st	25668

Head support	Dark grey plush		28534-60
	Grey Hygiene		28534-65
Without head support			
Head support attachment for bow handle			
Head support attachment for bow handle			25795
Without head support attachment			
Back support brace, foldable			
Back support brace, foldable			
Without back support brace			

Accessories

			Item. no.
Back support wedges	H:2, W:9, L:13 cm	pair	28030
	H:2, W:11, L:18.5 cm	pair	28139
	H:2.5, W:17.5, L:20 cm	pair	28140
Back support wedge plush	Short 15 cm	pcs	63012-60
	Long 30 cm	pcs	63013-60
R L			
Side cushion			63171-60-1 63171-60-2
Without side cushion			
Widening pads, back support			28637
Anti-tipper			28600
Anti-tipper Active			28520-1 28520-2
Without anti-tipper			

		R	L	Item. no.
Cane holder universal				27235
Without cane holder				
Spoke guard	Grey			
	Yellow			
Without spoke guard				
Tool kit				25174
Without tool kit				

Reference

User

Phone

Date

Comments