

# Seat height and balance position

1

| Camber angle |    |
|--------------|----|
| 0°           | 2° |



2

## Rear seat height and rear wheel

| Rear wheel | 20"   | 22"     | 24"   | 24"<br>Lightweight wheel<br>(carbon) | 26"    |
|------------|-------|---------|-------|--------------------------------------|--------|
| H          | 32 cm | 34,5 cm | 37 cm | 37 cm                                |        |
| G          | 34 cm | 36,5 cm | 39 cm | 39 cm                                |        |
| F          | 36 cm | 38,5 cm | 41 cm | 41 cm                                | 43 cm* |
| E          | 38 cm | 40,5 cm | 43 cm | 43 cm                                | 45 cm* |
| D          | 40 cm | 42,5 cm | 45 cm | 45 cm                                | 47 cm  |
| C          | 42 cm | 44,5 cm | 47 cm | 47 cm                                | 49 cm  |
| B          | 44 cm | 46,5 cm | 49 cm | 49 cm                                | 51 cm  |
| A          | 46 cm | 48,5 cm | 51 cm | 51 cm                                | 53 cm  |

\* Not with user brake in this position

3 Balance position



4

## Front seat height and castor

| Castor | Front fork extra short | Front fork short | Front fork Medium | Front fork Long | Castor housing             |                        |
|--------|------------------------|------------------|-------------------|-----------------|----------------------------|------------------------|
| 4"     | cm                     | -                | -                 | -               | Low widened castor housing |                        |
| 5"     | cm                     | cm               | -                 | -               |                            |                        |
| 6½"    | -                      | cm               | -                 | -               |                            |                        |
| 5"     | -                      | cm               | cm                | -               | Standard castor housing    | Widened castor housing |
| 6"     | -                      | cm               | cm                | -               |                            |                        |
| 6½"    | -                      | cm               | cm                | cm              |                            |                        |
| 7"     | -                      | -                | cm                | cm              |                            |                        |
| 8"     | -                      | -                | cm                | cm              |                            |                        |

4 mm fork thickness

5

## Hand rim

|                 |  |
|-----------------|--|
| Aluminium       |  |
| Stainless       |  |
| Cellular rubber |  |
| Vinyl coated    |  |
| Titanium        |  |

6

## Tyre rear wheel

|                            |  |
|----------------------------|--|
| Pneumatic high pressure 1" |  |
| Pneumatic 1 3/8"           |  |
| Solid                      |  |

7

## Single hand drive

|                           |  |
|---------------------------|--|
| Without single hand drive |  |
| Mounted on right side     |  |
| Mounted on left side      |  |

8

## Tyre castor 8"

|                           |  |
|---------------------------|--|
| 8" (200X35) semipneumatic |  |
| 8" (200X50) semipneumatic |  |
| 8" (200X30) pneumatic     |  |

For correct setting of the front fork attachment: Max. backward inclination 6 cm higher front than rear (9°), max. forward inclination 2 cm lower front than rear (3°).

Measures in brackets are theoretical. Hover over the editable text fields to see possible front seat height.