

# Convoid Coaster

## Canadian Retail Price List / Order Form

Account #: \_\_\_\_\_ Purchase Order #: \_\_\_\_\_ Date: \_\_\_\_\_

Contact: \_\_\_\_\_ ATP/RTS: \_\_\_\_\_ Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email:	

### Coaster Standard Features:

- Three fixed tilt setting options 10°, 20°, or 30° (Adjustable with Tools)
- Four Fixed Recline Setting Options (Adjustable with Tools)
- Tension Adjustable Back
- 2-Point Positioning Belt
- Seat Depth Growth of 2"
- Standard Convoid Upholstery
- Zippered Pockets in Seat & Back Upholstery
- Removable, Height Adjustable Swing-Away Footplates
- One Piece Height Adjustable Push Handle (Outside Mount)

### Standard Models:

- Coaster CTR12 - C800121
- Coaster CTR14 - C800141
- Coaster CTR16 - C800161
- Coaster CTR18 - C800181

- 80° Front Rigging
- Quick Release Rolko Rear Wheels
- Single Foot Operated Wheel Lock
- Clamshell Compact Fold
- Closure Strap
- Saddle Bags
- 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tire
- Powder Coated Frame
- Made in the USA

### Scout Standard Features:

- Hand Operated Brake
- Flip-Up Armrests
- 8" x 2" Front Ribbed, 12.5" x 2.5" Rear Pneumatic Knobby Tires
- Height & Tension Adjustable Back Positioner
- Height & Tension Adj Back Positioner Hardware

**\*Total for Scout Standard Features is \$1,047**

**All prices in CAD**

### Coaster Scout Models

- |  | Price          |
|--|----------------|
| <input type="checkbox"/> Coaster Scout CTR12 - C800123 | <b>\$3,622</b> |
| <input type="checkbox"/> Coaster Scout CTR14 - C800143 | <b>\$3,622</b> |
| <input type="checkbox"/> Coaster Scout CTR16 - C800163 | <b>\$3,622</b> |
| <input type="checkbox"/> Coaster Scout CTR18 - C800183 | <b>\$3,622</b> |

- |  | Price          |
|--|----------------|
| <input type="checkbox"/> Coaster CTR12 - C800121 | <b>\$2,574</b> |
| <input type="checkbox"/> Coaster CTR14 - C800141 | <b>\$2,574</b> |
| <input type="checkbox"/> Coaster CTR16 - C800161 | <b>\$2,574</b> |
| <input type="checkbox"/> Coaster CTR18 - C800181 | <b>\$2,574</b> |

### Select Tilt and Recline Setting

Factory Default Setting if No Selection is Made: 10° Fixed Tilt and 105° Recline

- Fixed Tilt set at 10°  
Fixed Recline set at:  100 |  105 |  110 |  115
- Fixed Tilt set at 20°  
Fixed Recline set at:  90 |  95 |  100 |  105
- Fixed Tilt set at 30°  
Fixed Recline set at:  80 |  85 |  90 |  95

### Tilt and Recline Options



### Upholstery Options

#### Colors:

- Sassy Purple
- Cherry Red
- Princess Pink
- Royal Blue
- Forest Green
- Panther Black

- C903851
- C903855
- C903759
- C903850
- C903857
- C903854



### Upholstery Color Options

#### Foam Kit:

- Positioning Wedges (Foam Kit for Zippered Upholstery Pockets)

#### Lumbar Support:

- Height and Tension Adjustable Back Positioner (Hardware Required)
- Height and Tension Adjustable Back Positioner Hardware

- |  | Price        |
|--|--------------|
| <input type="checkbox"/> Positioning Wedges (Foam Kit for Zippered Upholstery Pockets)     | <b>\$131</b> |
| <input type="checkbox"/> Height and Tension Adjustable Back Positioner (Hardware Required) | <b>\$165</b> |
| <input type="checkbox"/> Height and Tension Adjustable Back Positioner Hardware            | <b>\$333</b> |


CTR12	CTR14	CTR16	CTR18
C905567	C905568	C905569	C905570
C905536	C905536	C905536	C905536
C905538	C905618	C905539	C905619

**Seat Depth Extension (Standard 2")**

- 3" Seat Depth Extension Tubes

**Price**  
\$131

C905540	C905540	C905540	C905540
---------	---------	---------	---------



**Seat Depth Ranges**

Seat Width	Standard 2"	Option +3"
12"	10-12	12-15
14"	11-13	13-16
16"	12-14	14-17
18"	14.5-16.5	16.5-19.5

**Frame Options**

- Height Adjustable Flip-Up Armrests
- Caster Lock for 2" Wide Front Tires<sup>1</sup>
- Rear Anti-Tip Tubes (Cannot use with Curb Tipper)
- Curb Tipper (Cannot use with Rear Anti-Tip Tubes)

**Price**  
\$319  
\$131  
\$161  
\$137

<sup>1</sup>Manufacturer Installation Required.

CTR12	CTR14	CTR16	CTR18
C905611	C905611	C905611	C905611
C905526	C905526	C905526	C905526
C905527	C905527	C905527	C905527
C510737	C510737	C510737	C510737

**Transit Required Items<sup>^</sup>**

- Transit Anchor Brackets<sup>1^</sup>
- Headrest Extension<sup>^</sup>
- H-Harness with Padded Covers (3-Point Pelvic Belt Required)<sup>^</sup>
- 3-Point Pelvic Belt<sup>^</sup>

**Price**  
\$280  
\$181  
\$137  
\$68

**Foot Positioners - Select One Size Below:**

- Small 7" - 10"<sup>^</sup>
- Medium 8" - 11" (Default if no size is selected)<sup>^</sup>
- Large 9" - 12"<sup>^</sup>

\$219  
\$219  
\$219

<sup>1</sup>Manufacturer Installation Required.

<sup>^</sup>Transit required items are included in base price when selecting Coaster Transit Model

\*Crash Tested - WC19 compliant

CTR12	CTR14	CTR16	CTR18
C905531	C905531	C905531	C905531
C902863	C900484	C900138	C905571
C903988	C903988	C903988	C903988
C904017	C904017	C904017	C904017
C903415	C903415	C903415	C903415
C903416	C903416	C903416	C903416
C903417	C903417	C903417	C903417

**Transit Option**

- Q'Straint Transit Lap Belt

**Price**  
\$259

CTR12	CTR14	CTR16	CTR18
C903094	C903094	C903094	C903094

**Scout Option**

**Scout Standard Features:**

- Hand Operated Brake
- Flip-Up Armrests
- 8" x 2" Front Ribbed, 12.5" x 2.5" Rear Pneumatic Knobby Tire
- Height & Tension Adjustable Back Positioner
- Height & Tension Adjustable Back Positioner Hardware

**Price**  
\$76  
\$319  
\$154  
\$165  
\$333

- 8" x 2" Front Ribbed, 12.5" x 2" Rear Solid Knobby Tires

\$418

CTR12	CTR14	CTR16	CTR18
C905545	C905545	C905545	C905545

**Tire Options**

- 7.5" x 2" Front Ribbed, 12.5" x 2" Rear Solid Knobby Tires
- 7.5" x 2" Front, 12.5" x 2" Rear Pneumatic Tires

**Price**  
\$418  
\$418

CTR12	CTR14	CTR16	CTR18
C905535	C905535	C905535	C905535
C905561	C905561	C905561	C905561

**Head Supports**

- Padded Headwing
- Occi Headwings
- Headrest Extension

**Price**  
\$340  
\$319  
\$181

CTR12	CTR14	CTR16	CTR18
C905572	C905573	C905574	C905575
C905270	C905271	C905271	C905272
C902863	C900484	C900138	C905571

**Trunk Positioning**

- H-Harness with Padded Covers (3-Point Pelvic Belt Required)
- Full Torso Support Vest

**Price**  
\$137  
\$369

CTR12	CTR14	CTR16	CTR18
C903988	C903988	C903988	C903988
C905598	C905599	C905600	C905601

**Pelvic Positioning**

- 3-Point Pelvic Belt

**Price**  
\$68

CTR12	CTR14	CTR16	CTR18
C904017	C904017	C904017	C904017

**Foot & Leg Positioning**

- Medial Thigh Support, Abductor
- Lateral Thigh Support, Adductor
- Calf Panel\*
- Front Rigging Options: (80° Fixed Front Rigging Standard)**
- Elevating Leg Rest (90° - 180°)
- Padded Elevating Legrest Cover

**Price**  
\$227  
\$227  
\$86

**Footplate Options:**

- Individual Angle Adjustable Footplates (Elevating Leg Rest Required)
  - Elevating Leg Rest (90° - 180°)
- Angle Adjustable Footplates
- One Piece Angle Adjustable Flip-Up Footplate (-15° - +15°)
- Fully Padded Footbox (One Piece Footplate Required)
- Footplate Securement Strap

\$350  
\$131  
\$223  
\$350  
\$223  
\$155  
\$454  
\$43

**Foot Positioners - Select One Size Below:**

- Small 7" - 10"
- Medium 8" - 11"
- Large 9" - 12"

\$219  
\$219  
\$219

\*Calf Panel cannot be used when footplate extension tube is rotated, and in the highest position.

CTR12	CTR14	CTR16	CTR18
C905576	C905577	C905578	C905579
C905580	C905581	C905582	C905583
C904580	C904580	C904580	C904580
C905541	C905542	C905543	C905544
C904619	C904619	C904619	C904619
C905624	C905625	C905626	C905627
C905541	C905542	C905543	C905544
C903445	C903268	C903268	C903269
C905557	C905558	C905559	
C904954	C904955	C904956	
C903584	C903584	C903584	C903584
C903415	C903415	C903415	C903415
C903416	C903416	C903416	C903416
C903417	C903417	C903417	C903417

<u>Medical Support Accessories</u>	<u>Price</u>
<input type="checkbox"/> Medical Necessity Storage Basket	<b>\$412</b>
<input type="checkbox"/> IV Pole (Height Adjustable and Collapsible)^	<b>\$267</b>
<input type="checkbox"/> Adj. Swing Away Hardware (Required for IV Pole)	<b>\$192</b>

CTR12	CTR14	CTR16	CTR18
C905587	C905588	C905588	C905589
C904566	C904566	C904566	C904566
C905612	C905612	C905612	C905612

<u>Accessories</u>	<u>Price</u>
<input type="checkbox"/> Headrest Cover (Canopy)	<b>\$319</b>
<input type="checkbox"/> Extended Headrest (Canopy)	<b>\$350</b>
<input type="checkbox"/> Under Seat Storage Basket	<b>\$160</b>
<input type="checkbox"/> Utility Bag	<b>\$59</b>
<input type="checkbox"/> Rain Cover (Canopy Required)	<b>\$211</b>
<input type="checkbox"/> Mosquito Net (Canopy Required)	<b>\$211</b>

CTR12	CTR14	CTR16	CTR18
C905017	C905018	C905019	C905020
C905602	C905603	C905604	C905605
C905584	C905585	C905585	C905586
C904531	C904531	C904531	C904531
C905590	C905591	C905592	C905593
C905594	C905595	C905596	C905597

The Coaster transit option includes four manufacturer installed wheelchair transportation anchors. The Coaster transit model has been crash tested and conforms to RESNA WC4:2012 section 19 standards. Anthropomorphic Test Device weights are CTR12: 66 lbs./30kg; CTR14: 100 lbs./45.5kg; CTR16 and CTR18: 170 lbs./77kg.

The Convoid Coaster was tested equipped with a head support and wheelchair-anchored pelvic lap belt that comply with RESNA WC4:2012 section 18. The performance and safety of accessories and/or other seating options not tested by Convoid is unknown, and no warranty or claim about their performance is being, or can be, made.

**Standard Shipping Rates: Ground \$71.95 • For expedited shipping rates, please call 1-844-876-6245 or email [convoidsales.us@etac.com](mailto:convoidsales.us@etac.com)**

Prices in CAD & subject to change, ex works Torrance, California.

Effective: 4/24/2024

Convoid  
 2830 California Street  
 Torrance, CA 90503  
 1 844-US-MOBILITY (844-876-6245)  
[www.etac.us.com](http://www.etac.us.com)

OF2602-ENG-REV.02

