

# Convoid Coaster Scout - HCPCS E1236

## Retail Price List / Order Form

Account #: \_\_\_\_\_ Purchase Order #: \_\_\_\_\_ Date: \_\_\_\_\_  
 Contact: \_\_\_\_\_ ATP/RTS: \_\_\_\_\_ Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email:	

**Coaster Scout Standard Features:**

- Three fixed tilt setting options 10°, 20°, or 30° (Adjustable with Tools)
- Four Fixed Recline Setting Options (Adjustable with Tools)
- Height and Tension Adjustable Back Positioner
- Tension Adjustable Back Straps
- 2-Point Pelvic Belt
- Seat Depth Growth of 2"
- 8" x 2" Front Ribbed, 12.5 x 2.5" Rear Pneumatic Knobby Tire
- Zippered Pockets in Seat & Back Upholstery
- Removable, Height Adjustable Swing-Away Footplates
- One Piece Height Adjustable Push Handle (Outside Mount)

**Coaster Scout Transit Additional Features:**

- 80° Front Rigging
- Attendant Handbrake
- Flip-Up Armrest
- Saddle Bags
- Clamshell Compact Fold
- Closure Strap
- Standard Convoid Upholstery
- Powder Coated Frame
- Made in the USA

- Transit Anchor Brackets
- Headrest Extension
- H-Harness with Padded Covers
- 3-Point Pelvic Belt
- Foot Positioners

**Standard Scout Models:**

- Coaster Scout CTR12 - C800123
- Coaster Scout CTR14 - C800143
- Coaster Scout CTR16 - C800163
- Coaster Scout CTR18 - C800183

**Price**  
**\$2,766**  
**\$2,766**  
**\$2,766**  
**\$2,766**

**Transit Included Models:**

- Coaster Scout Transit CTR12T - C800124
- Coaster Scout Transit CTR14T - C800144
- Coaster Scout Transit CTR16T - C800164
- Coaster Scout Transit CTR18T - C800184

**Price**  
**\$3,442**  
**\$3,442**  
**\$3,442**  
**\$3,442**

**Select Tilt and Recline Setting**

- Factory Default Setting if No Selection is Made: 10° Fixed Tilt and 105° Recline
- Fixed Tilt set at 10°  
Fixed Recline set at: 100 | 105 | 110 | 115
  - Fixed Tilt set at 20°  
Fixed Recline set at:  90 |  95 | 100 | 105
  - Fixed Tilt set at 30°  
Fixed Recline set at:  80 |  85 |  90 |  95

**Tilt and Recline Options**



**Upholstery Options**

**Colors:**

- Sassy Purple
- Cherry Red
- Princess Pink
- Royal Blue
- Forest Green
- Panther Black

C903851  
 C903855  
 C903759  
 C903850  
 C903857  
 C903854

**Upholstery Color Options**



- Foam Kit:**
- Positioning Wedges (Foam Kit for Zippered Upholstery Pockets)

**Price**  
**\$100**      **HCPCS**  
**K0108**

<b>CTR12</b>	<b>CTR14</b>	<b>CTR16</b>	<b>CTR18</b>
C905567	C905568	C905569	C905570

<u>Seat Depth Extension (Standard 2")</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> 3" Seat Depth Extension Tubes	\$100	E1297

<u>CTR12</u>	<u>CTR14</u>	<u>CTR16</u>	<u>CTR18</u>
C905540	C905540	C905540	C905540



<u>Seat Depth Ranges</u>		
<u>Seat Width</u>	<u>Standard 2"</u>	<u>Option +3"</u>
12"	10-12	12-15
14"	11-13	13-16
16"	12-14	14-17
18"	14.5-16.5	16.5-19.5

<u>Frame Options</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Height Adjustable Flip-Up Armrests	\$244	E0973
<input type="checkbox"/> Caster Lock for 2" Wide Front Tires <sup>1</sup>	\$100	K0073
<input type="checkbox"/> Rear Anti-Tip Tubes (Cannot use with Curb Tipper)	\$123	E0971
<input type="checkbox"/> Curb Tipper (Cannot use with Rear Anti-Tip Tubes)	\$105	

<sup>1</sup>Manufacturer Installation Required.

<u>CTR12</u>	<u>CTR14</u>	<u>CTR16</u>	<u>CTR18</u>
C905611	C905611	C905611	C905611
C905526	C905526	C905526	C905526
C905527	C905527	C905527	C905527
C510737	C510737	C510737	C510737

<u>Transit Required Items<sup>^*</sup></u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Transit Anchor Brackets <sup>1^</sup>	\$214	K0108
<input type="checkbox"/> Headrest Extension <sup>^</sup>	\$138	E0966
<input type="checkbox"/> H-Harness with Padded Covers (3-Point Pelvic Belt Required) <sup>^</sup>	\$105	E0960
<input type="checkbox"/> 3-Point Pelvic Belt <sup>^</sup>	\$52	E0978
<b>Foot Positioners - Select One Size Below:</b>		
<input type="checkbox"/> Small 7" - 10" <sup>^</sup>	\$167	K0108
<input type="checkbox"/> Medium 8" - 11" (Default if no size is selected) <sup>^</sup>	\$167	K0108
<input type="checkbox"/> Large 9" - 12" <sup>^</sup>	\$167	K0108

<sup>1</sup>Manufacturer Installation Required.  
<sup>^</sup>Transit required items are included in base price when selecting Coaster Transit Model  
 \*Crash Tested - WC19 compliant

<u>CTR12</u>	<u>CTR14</u>	<u>CTR16</u>	<u>CTR18</u>
C905531	C905531	C905531	C905531
C902863	C900484	C900138	C905571
C903988	C903988	C903988	C903988
C904017	C904017	C904017	C904017
C903415	C903415	C903415	C903415
C903416	C903416	C903416	C903416
C903417	C903417	C903417	C903417

<u>Transit Option</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Q'Straint Transit Lap Belt	\$198	E0978

<u>CTR12</u>	<u>CTR14</u>	<u>CTR16</u>	<u>CTR18</u>
C903094	C903094	C903094	C903094

<u>Tire Options</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> 8" x 2" Front Ribbed, 12.5" x 2" Rear Solid Knobby Tires	\$319	

<u>CTR12</u>	<u>CTR14</u>	<u>CTR16</u>	<u>CTR18</u>
C905545	C905545	C905545	C905545

<u>Head Supports</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Padded Headwing	\$260	E0955
<input type="checkbox"/> Occi Headwings	\$244	E0955
<input type="checkbox"/> Headrest Extension	\$138	E0966

<u>CTR12</u>	<u>CTR14</u>	<u>CTR16</u>	<u>CTR18</u>
C905572	C905573	C905574	C905575
C905270	C905271	C905271	C905272
C902863	C900484	C900138	C905571

<u>Trunk Positioning</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> H-Harness with Padded Covers (3-Point Pelvic Belt Required)	\$105	E0960
<input type="checkbox"/> Full Torso Support Vest	\$282	K0108

<u>CTR12</u>	<u>CTR14</u>	<u>CTR16</u>	<u>CTR18</u>
C903988	C903988	C903988	C903988
C905598	C905599	C905600	C905601

<u>Pelvic Positioning</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> 3-Point Pelvic Belt	\$52	E0978

<u>CTR12</u>	<u>CTR14</u>	<u>CTR16</u>	<u>CTR18</u>
C904017	C904017	C904017	C904017

<u>Foot &amp; Leg Positioning</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Medial Thigh Support, Abductor	\$173	E0957
<input type="checkbox"/> Lateral Thigh Support, Adductor	\$173	E0956
<input type="checkbox"/> Calf Panel*	\$66	E0995
<b>Front Rigging Options: (80° Fixed Front Riggings Standard)</b>		
<input type="checkbox"/> Elevating Leg Rest (90° - 180°)	\$267	K0990
<input type="checkbox"/> Padded Elevating Legrest Cover	\$100	K0108
<b>Footplate Options:</b>		
<input type="checkbox"/> Individual Angle Adjustable Footplates (Elevating Leg Rest Required)	\$170	K0040
<input type="checkbox"/> Elevating Leg Rest (90° - 180°)	\$267	K0990
<input type="checkbox"/> Angle Adjustable Footplates	\$170	K0040
<input type="checkbox"/> One Piece Angle Adjustable Flip-Up Footplate (-15° - +15°)	\$118	K0040
<input type="checkbox"/> Fully Padded Footbox (One Piece Footplate Required)	\$347	K0108
<input type="checkbox"/> Footplate Securement Strap	\$33	
<b>Foot Positioners - Select One Size Below:</b>		
<input type="checkbox"/> Small 7" - 10"	\$167	K0108
<input type="checkbox"/> Medium 8" - 11"	\$167	K0108
<input type="checkbox"/> Large 9" - 12"	\$167	K0108

\*Calf Panel cannot be used when footplate extension tube is rotated, and in the highest position.

<u>CTR12</u>	<u>CTR14</u>	<u>CTR16</u>	<u>CTR18</u>
C905576	C905577	C905578	C905579
C905580	C905581	C905582	C905583
C904580	C904580	C904580	C904580
C905541	C905542	C905543	C905544
C904619	C904619	C904619	C904619
C905624	C905625	C905626	C905627
C905541	C905542	C905543	C905544
C903445	C903268	C903268	C903269
C905557	C905558	C905559	
C904954	C904955	C904956	
C903584	C903584	C903584	C903584
C903415	C903415	C903415	C903415
C903416	C903416	C903416	C903416
C903417	C903417	C903417	C903417

<u>Medical Support Accessories</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Medical Necessity Storage Basket	\$315	
<input type="checkbox"/> IV Pole (Height Adjustable and Collapsible) <sup>^^</sup>	\$204	K0105
<input type="checkbox"/> Adj. Swing Away Hardware (Required for IV Pole)	\$147	E1028

<sup>^^</sup>Cannot be used with canopy.

<u>CTR12</u>	<u>CTR14</u>	<u>CTR16</u>	<u>CTR18</u>
C905587	C905588	C905588	C905589
C904566	C904566	C904566	C904566
C905612	C905612	C905612	C905612

<u>Accessories</u>	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Headrest Cover (Canopy)	\$244	K0108
<input type="checkbox"/> Extended Headrest (Canopy)	\$267	K0108
<input type="checkbox"/> Under Seat Storage Basket	\$122	
<input type="checkbox"/> Utility Bag	\$45	
<input type="checkbox"/> Rain Cover (Canopy Required)	\$161	
<input type="checkbox"/> Mosquito Net (Canopy Required)	\$161	

<u>CTR12</u>	<u>CTR14</u>	<u>CTR16</u>	<u>CTR18</u>
C905017	C905018	C905019	C905020
C905602	C905603	C905604	C905605
C905584	C905585	C905585	C905586
C904531	C904531	C904531	C904531
C905590	C905591	C905592	C905593
C905594	C905595	C905596	C905597

The Coaster transit option includes four manufacturer installed wheelchair transportation anchors. The Coaster transit model has been crash tested and conforms to RESNA WC4:2012 section 19 standards. Anthropomorphic Test Device weights are CTR12: 66 lbs./30kg; CTR14: 100 lbs./45.5kg; CTR16 and CTR18: 170 lbs./77kg.

The Convoid Coaster was tested equipped with a head support and wheelchair-anchored pelvic lap belt that comply with RESNA WC4:2012 section 18. The performance and safety of accessories and/or other seating options not tested by Convoid is unknown, and no warranty or claim about their performance is being, or can be, made.

*HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.*

**Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245 or email [convoidsales.us@etac.com](mailto:convoidsales.us@etac.com)**

Prices in USD & subject to change, ex works Torrance, California.

Effective: 4/24/2024

Convoid  
2830 California Street  
Torrance, CA 90503  
1 844-US-MOBILITY (844-876-6245)  
[www.etac.us.com](http://www.etac.us.com)

OF2610-ENG-REV.03

 **etac**®  
Creating Possibilities