



Etac Flex grab rail

Flexibility and freedom

 **etac**[®]

Etac Flex grab rail

Innovative design that creates possibilities

Flexible modular grab rails

Etac Flex grab rail is the perfect choice for those who need a little extra support in the bathroom. Thanks to its flexible modular function, it is easy to adjust the handle to suit both you and your bathroom environment. The pure Scandinavian design blends in with other products and you can choose between two colours - white and grey.

Support for good and bad days

The material of the grab rail and its well-designed shape and surface provide a stable grip even if your hands or the handles are wet. The triangular shape and the thumb support also make you feel safer.



Etac Flex fits perfectly in the bathroom environment.


Etac Flex screw mount		Item no.	Weight
- Extension Flex 30 cm white		81706110	0.35 kg per unit (0.77 lbs)
- Grab rail Flex 30 cm white		81707110	
- Grab rail Flex 60 cm white		81707120	Maximum user weight 100 kg
- Grab rail Flex 90 cm white		81707130	
- Extension Flex 30 cm grey		81706310	Material Polypropylene
- Grab rail Flex 30 cm grey		81707310	
- Grab rail Flex 60 cm grey		81707320	Colour White and Volcano grey
- Grab rail Flex 90 cm grey		81707330	
Etac Flex adhesive mount		Item no.	
- Extension Flex 30 cm white		81706210	
- Grab rail Flex 30 cm white		81707210	
- Grab rail Flex 60 cm white		81707220	
- Grab rail Flex 90 cm white		81707230	Clean the product with a solvent-free cleaning agent (pH level 5–9), or with a 70% disinfectant solution. Can be decontaminated at max 85 °C / 185 °F if so required.
- Extension Flex 30 cm grey		81706410	
- Grab rail Flex 30 cm grey		81707410	
- Grab rail Flex 60 cm grey		81707420	
- Grab rail Flex 90 cm grey		81707430	

Note: To add one or more grab rail to a 30 cm, 60 cm or 90 cm combination, it is necessary to use the Etac Flex Extensions.



Easy to adjust

The modular system can easily be assembled into a length and shape that suits you and your bathroom.




Secure grip

The well-designed shape and surface of the grab rail provide a stable grip even if your hands are wet.



Creating Possibilities

You can connect as many grab rails as you want and choose the angles that suit you. This creates possibilities.



Attach in two ways

Attach with screws or adhesive. The latter is perfect for those who cannot or do not want to drill in the bathroom.



Etac Flex offers a wide range of positioning options when connecting grab rails. Choose the angle between -90° to $+90^\circ$ that suits your need.





Etac is a world-leading developer and provider of ergonomic assistive devices and patient handling equipment. Our heart lies in the solutions that optimise quality of life for the individual, their family and caregivers.

For the latest news and continuously updated product information – visit: www.etac.com

Etac AB
Kista Science Tower
SE-164 51 Kista, Sweden
Tel +46 371 - 58 73 30
info@etac.se www.etac.com


Creating Possibilities