

etac[®]
Creating Possibilities



molift[®]

Molift Air 200

The compact, powerful hoist optimized for long term care

Molift Air 200

Molift Air 200 is an affordable, well-designed and highly functional ceiling hoist for users up to 205 kg. Combined with Molift Rail System and Molift slings, this hoist solves most day-to-day hoisting situations in health care and other patient handling environments.

Comfort and quality

Molift Air 200 combines the need for maximum comfort with efficient and safe hoisting situations. The design and user-friendliness bring an exceptional freedom, which enables a single caregiver to perform daily lifting situations, allowing more time for caring for the individual user. For the caregiver, a better working environment is possible without the need for manual lifting. The user experiences a dignified, safe and comfortable lifting.



The flexible hoist

Molift Air 200 is suitable for sitting and horizontal transfers, as well as standing and gait training situations.

Molift Air 200 is used with fixed ceiling rails but can also be used for portable free-standing systems. It is easy to handle during installation and servicing, due to a very low weight and the smart quick release function.

Molift Air 200 offers flexible solutions such as remote control and possibilities for room-to-room transfer. Add to this the safety features, the large hoisting interval and the multiple range of sling bars. This hoist is the complete solution in all health care segments.



Product model	Including	Item no.	Facts	
Molift Air 200 SLA Battery	Hoist, handcontrol, charger	26200	Expected lifetime	30 000 lifts SWL/10 years
Molift Air 200 Litium Battery	Hoist, handcontrol, charger	26201	Molift Service Software	Included
Facts			Emergency lowering	Mechanical and electrical
Maximum user weight (SWL)		205 kg	Emergency stop	Reachable from floor with string. Reset possible from floor
Total weight		8,6 kg	Protection class; hand control & hoist motor	IP24
Length		360 mm	 Clean / disinfect the product using a standard non-abrasive cleaning agent with a pH between 5 and 9, or with a 70% disinfectant solution.	
Width		190 mm		
Height		195 mm		
Battery SLA		12 V SLA 2,9 Ah		
Battery Litium		25,6 V LiFePO4 2,4 Ah		
Battery charger		100-240 V AC, 50-60 Hz		
Battery charge time		6 hrs		
Motor		24 V DC		
Hoisting interval		3000 mm		
Hoisting speed 75 kg		60 mm/sec		
				

Features



Small hoist - great lifting capacity

Molift Air 200 is very small and weighs only 8,6 kg. Nevertheless it has a lifting capacity (SWL) of 205 kg. The size and weight makes it very easy to operate and move for the carer during different transfer situations.



Silent and soft operation

Molift Air 200 works very quietly and has a soft start and stop operation. This results in a comfortable and safe feeling for the user during the transfer situation.



Quick release device for trolley

Molift Air 200 is easy to handle during installation and servicing, due to a very low weight and the smart quick release function. With a simple grip, the motor is disconnected from its trolley and removed.



Intuitive hand control

Molift Air 200 has a new intuitive handcontrol with attachment possibility on sling bar and optional docking station for wall.

Accessories



4-point Sling bar

Light weight aluminium, SWL 350 kg

- 4-pt Sling bar 350 (S)2530210
- 4-pt Sling bar 450 (M)2530211
- 4-pt Sling bar 550 (L)2530212
- 4-pt Sling bar 650 (XL)2530213



Molift Air Scale

This scale can calculate BMI (Body Mass Index) with a simple touch. Approved Class III.

- Scale SWL 300 kg.....1840000



2-point Sling bar

Light weight aluminium, SWL 350 kg

- 2-pt Sling bar 350 (S)2530115
- 2-pt Sling bar 450 (M)2530120
- 2-pt Sling bar 550 (L)2530125



Molift Stretcher MRI

Flexible scoop stretcher, ideal for horizontal hoisting situations at hospitals and institutions. SWL 300 kg

- Molift Stretcher MRI..... 2150108



2-point Pigtail Sling bar

Light weight aluminium, SWL 350 kg

- 2-pt Pigtail Sling bar 350 (S)1830033
- 2-pt Pigtail Sling bar 450 (M)1830032
- 2-pt Pigtail Sling bar 550 (L)1830031



Other accessories

- Room-to-room unit 2510295



8-point Sling bar

For transfer situations with a fabric stretcher. Steel, SWL 300 kg

- 8-pt Sling bar2140003



Etac is a world-leading developer of innovative assistive devices and patient handling equipment. Since 1973, we have been committed to improving quality of life for the individual, family members and caregivers.

For the intended purpose of the products and latest information, please visit www.etac.com

Etac

- +46 371 587 330
- info@etac.se
- www.etac.com

etac[®]
Creating Possibilities