


Avant

Valid from serial no: 17C 33000



Pos	Description	Item no
1	Handle strap black	25119
2	Brake handle compl (up to ser no 010543)	81636
	Brake handle compl (from ser no 010544)	82154
3	Knob and security lock kit	81639
4	Frame extension front	26203
5	Frame extension rear	81638
6	Wheel kit, 2 pcs	81947
7	Dust cover wheel, 4 pcs	81637
8	Seat compl wood (up to ser no 010543)	25673
	Seat compl plastic (from ser no 010544) ..	82030
9	Basket	25659-15
10	Crutch holder	24940

74077C

09-06-04

Etac AB, Box 203, SE-334 24 Anderstorp
Tel. +46/371-58 73 00, Fax. +46/371-58 73 90

www.etac.com

etac