


Tango

Valid from serial number: 27A 06614



Pos	Prescription	Art no	Pos	Prescription	Art no
1	Screw PT-KX K50x 18	11412	15	Cross pipe Tango right	25803-42-1
2	Handle bow Tango	25917		Cross pipe Tango left	25803-42-2
3	Brake tubes Tango compl	25842	16	Screw M 6x 55 MC6S 12.9	11351
4	Brake plate Tango kompl right	25843-1	17	Prifile washer Tango D22 d6.3	50437
	Brake plate Tango compl left	25843-2	18	Domed nut M 6 black	12005
5	Spring Tango brake lower	80928	19	Dust cover 22x1.5 black	50070
6	Knob Tango	80930	20	Screw B 8x 16 M6S RXS	11301
7	Slidebushing Tango	50457	21	Slidebar right Tango	50447
8	Tipper Tango	50443		Slidebar left Tango	50448
9	Axle wheel Tango back	80929	22	Side frame Tango right	25844-14-1
10	Wheel Tango 210 mm solid	50453		Side frame Tango left	25844-14-2
11	Starlock black chrome with cap	13073	23	Bushing Tango	50456
12	Washer BRB HB200 d13	13018	24	Bushing to tube Tango	50442
13	Front fork Tango compl	25845	25	Seat Tango	50451
14	Screw PT-KX K40x 8	11409	26	Brakeleash Tango	50459
			27	Plug Tango	50458
			28	Handle right Tango	50445
				Handle left Tango	50446

74170C

05-05-31

Etac AB, Box 203, SE-334 24 Anderstorp

Tel. +46/371-58 73 00, Fax. +46/371-58 73 90

www.etac.com