

80209240 / 75 mm, 84005033 / 125 mm

80209263 / 75 mm, 84001280 / 125 mm

80209264 / 75 mm, 84001281 / 125 mm


84005034 / 125 mm

Länkhjul Clean 75 / 125 mm

Castor Clean 75 / 125 mm

Lenkrad Clean 75 / 125 mm

78208F 2018-03-29



Etac Supply Center
Box 203
SE-334 24 Anderstorp
Tel +46 371 58 73 00
Fax +46 371 58 73 90
www.etac.com

pH
min 5
max 9

85 °C
185 °F
3 min

