


molift[®]
a part of Etac



Molift Service Tool

The software that interacts with the hoist

All Molift hoists have integrated service software that calculates when a service is necessary but it also records lifter use and lifting patterns. This provides a good platform both for technical service and maintenance but also an information data base for lifter usage. It is unique software that reassures user, caregiver and management that the hoist is safe in relation to applicable laws and regulations.



Ensuring correct service and maintenance

The automatic service notification software calculates when service is required and gives a warning signal. The automatic warning system will alert the end-user when service is due.

The technician can, via the Molift Service Tool 4, connect and interact with the lifter to know which service is necessary and for technical fault finding. The software allows also reading, checking and adjusting key parameters in the lifter, as well as receiving authorized service history and in printable format.

Measure, record and analyze hoist usage and patterns

Through the Molift Service Tool 4 the software can also be used to retrieve data about hoists usage, lifting patterns, etc.

This is very useful when analyzing staff "sick-leave" absence due to injury and whether the types of hoists are suitable for the patient group and their needs. It can also be a helpful notification system for the management that personnel are utilizing the equipment to reduce injury and maximize efficiency.

Features

- Automatic service notification
- Platform for technical service and management
- Measure, record and report usage and patterns
- Simplification of service planning
- Prevention of interrupted use

Application Areas

- Assessment of service and service planning
- Periodic inspection or recycling
- Documentation of service and inspection
- Management analysis of lifter usage, lifting patterns

For service and technicians

- Read, check hoist use and statistics
- Documentation of service and history
- Adjust key parameters

For Management

- Analysis use and lifting patterns
- Analysis solution efficiency
- Analysis equipment solution

Model	Item no.
Molift Service Tool 4, New user, incl. software pack, cables and new user registration	1299011
Molift Service Tool 4, Software upgrade, only available for existing registered user	1299013