Combination agreement



Wheelchairs

1A15211.8 21-11-17

The list of combination agreements is only an overview. For more detailed information regarding the agreements, please contact our customer service. The maximum user weight for the combination is always the same as for the product in the combination that has the lowest maximum user weight.

Manufacturer	Product	Etac Product
Alber	e-motion M15, M25	Cross 6, Elite, Cross, Cross XL, Cross 5, Cross 5 XL, M100, Twin, Prio, Crissy Swing-Away, Crissy Active
Alber	e-fix E25/E26	Cross 6, Cross, Cross XL, Cross 5, Cross 5 XL Transit, Twin, M100
Alber	Quix Q10	Cross, Cross XL, Transit, Twin
Alber	Scalamobil S30 equipped with Scalaport X7	Act, Elite, Next, Next Comfort, M100
Alber	Scalamobil S35/S38	Cross 6, Cross, Cross XL, Cross 5, Cross 5 XL, Transit, Twin, M100, Prio, Crissy Swing-Away, Crissy Active
Alber	Viamobil eco V14	Cross, Cross XL, Cross 5, Cross 5 XL, Twin, M100, Crissy Swing-Away, Crissy Active
Alber	Viamobil V15	Cross, Cross XL, Transit, Twin, M100
Alber	Viamobil V25	Cross 6, Cross, Cross XL, Cross 5, Cross 5 XL, M100, Twin, Transit, Prio, Crissy Swing-Away, Crissy Active
Alber	e-fix E35/E36	Cross 6, Cross, Cross XL, Cross 5, Cross 5 XL, M100, Twin, Transit, Prio, Crissy Swing-Away, Crissy Active
Alber	e-pilot P15	Cross 6, Crissy Swing Away, Crissy Active
Alber	SMOOV 010	Cross 6, Cross 5, Cross 5 XL, Crissy Swing Away, Crissy Active
Alu Rehab A/S	Netti Contour, UNO seat, Netti Sit, UNO Back, Netti Smart, Netti Stabil, Netti Super Stabil, Netti Kyfotic Back, Netti Nordic Back.	Cross, Cross XL, Cross 5, Cross 5 XL, Act, Elite, M100, M100 TR, Transit III, Next, Next Comfort, Etac Sting.
Anatomic Sitt	Proserie Headrest	Cross 6, Next, Cross, Cross XL, Twin, Transit III, Sting Act, Elite, Crissy Swing Away, Crissy Active, Cross 5, Cross 5 XL, Prio
Armbike/Jonny Borg	Armbike	Act, Elite
Bodypoint, Inc	Bodypoint 2-pt and 4-pt Hip Belts: NB200-series, HB100-series, HB200-series, HB400-series, EB200-series, LB200-series, LB400-series, LH600-series, SH100-series, SH200-series, SH300-series, SP110-series, SP111-series	Cross, Cross XL, Cross 5, Cross 5 XL, Cross 6, Twin, Transit, M100, M100 TR, Prio
Decon AB	Easy Wheel	Cross, Twin, Act, Elite, Prio
Decon AB	E-Move, E-Drive	Cross 5, Cross 5 XL, Cross 5 Active, Cross 6, Prio
Enkatsu Solutions AB	Enkatsu Climber	Cross, Twin, Transit, Ono Rollator, Clean och Mobil
G.Niklassons	Maxgrip Rear wheel	Next, Cross XL, Cross, Twin, Sting, Act, Elite
HEA Medical AB	HEA attachment for chest slings	Cross 5, Cross 5 XL, Cross, Cross XL, Cross 6, M100 with height adjustable push handles, Transit with height adjustable push handles
Janmek HB	Holder for oxygen tanks Holder for liquid oxygen tanks IV stand	Cross 6, Cross, Cross XL, Cross 5, Cross 5XL, Twin, Transit III, Elite, Act, Sting, Next, Next Comfort, HandiVipp Frame Only, HandiVipp, M100, M100TR, M100 VB

Combination agreement Wheelchairs

oetac[®]

1A15211.8 21-11-17

Manufacturer	Product	Etac Product
Max mobility	Smart Drive MX2+	Cross 5, Cross 5 XL, Cross, Cross XL
Medema	M-Drive power assist add on	Cross, Cross XL, Cross VB & Cross XL VB
Mercado Medic	Real 9000, Real 9000 Plus, Real 6100, Real 6100 Plus	Amputee leg support, Leg support- wide -narrow and -narrow short, Foot supports
Otto Bock Scandinavian AB	Head and Necks supports, Foot Buckets and Calf Pads, Firm Seating and Back systems and Cushions	Next. Cross XL, Cross, Twin, Transit III
Performance Health	Backrest V-Trak Products LTD	Cross, Cross XL, Next, Next Comfort, Twin, Transit
R82	Head and Necks support	Cross, Cross XL, Cross 5, Next, Next Comfort, Elite, Sting, M100, M100 TR
Sunrise Medical AB	Backrest Jay 2, standard, deep & Back Jay 3	Cross 6, Cross 5, Cross 5 XL, Prio, Cross, Cross XL, Next, Next Com- fort, Twin, Transit
Sunrise Medical AS	PowerTec F16	Cross
Sunrise Medical Europe	Quickie Wheeldrive	Cross 5
Särnmark Hjälpmedel AB	Hemi-armrest Positioning seat & Active badge	Etac's wheelchair
Tzora Active Systems	Combi Power, Combi Power Comfort	Cross, Cross XL, Transit III
UB-Let	EasyGo+	Cross XL, Twin XL