



Internal Product Information

July 24th, 2019 | Concerning: Product improvement | Molift 2 and 4-point sling bars

Product improvement for Molift 2 and 4-point sling bars

Etac will now implement the new zinc alloy hooks on all the 2 and 4-point sling bars, similar to the hooks implemented with 2-point sling bars late 2018.

The implementation will take place during the next couple of weeks.

New price lists will be distributed week 31 / 2019.

New material and design of the sling hook

The sling hook material has changed from mold cast composite, to a powder coated mold cast zinc-alloy. The new material, zinc alloy, provides improved features for the sling bars in terms of manufacturing, durability and mechanical features.

For sling bars to Molift Air (all versions), the SWL lifting capacity is increased from 300kg to 350 kg (770lbs) and the pivoting feature for the sling bars is replaced by a welded lifting eye.

This leaves us with a new set of item numbers for our sling bar portfolio. See complete table chart on the next page.

The new design is also contributing to a reduced risk of having the sling loops unattached unintendedly from the sling bar during lifting. There is no modification in terms of operating and using the new sling bars.

Outgoing design



New design



Do not hesitate to ask if any clarifications are necessary.

1/2



Internal Product Information

July 24th, 2019 | Concerning: Product improvement | Molift 2 and 4-point sling bars

Item	Description	Width (mm)	Swivel	SWL	Air (all versions)	Mover 180	Mover 205	Mover 300	Partner 255
1830043	2pt Sling Bar 350 (S) for Mover 180	350	Yes	300 kg		X			
1830042	2pt Sling bar 450 (M) for Mover 180	450	Yes	300 kg		X			
1830041	2pt Sling bar 550 (L) for Mover 180	550	Yes	300 kg		X			
2530115	2pt Sling bar 350 (S) for Air	350	No	350 kg	X				
2530120	2pt Sling bar 450 (M) for Air	450	No	350 kg	X				
2530125	2pt Sling bar 550 (L) for Air	550	No	350 kg	X				
1430115	2pt Sling bar 550 (S) for Mover 205	550	Yes	300 kg			X	(X)	
1430116	2pt Sling bar 450 (M) for Mover 205	450	Yes	300 kg			X	(X)	
1430117	2pt Sling bar 350 (L) for Mover 205	350	Yes	300 kg			X	(X)	
1830110	4pt Sling bar 350 (S) for Mover 180	350	Yes	300 kg		X			
1830111	4pt Sling bar 450(M) for Mover 180	450	Yes	300 kg		X			
1830112	4pt Sling bar 550 (L) for Mover 180	550	Yes	300 kg		X			
1830113	4pt Sling bar 650 (XL) for Mover 180	650	Yes	300 kg		X			
2530210	4pt Sling bar 350 (S) for Air	350	No	350 kg	X				
2530211	4pt Sling bar 450 (M) for Air	450	No	350 kg	X				
2530212	4pt Sling bar 550 (L) for Air	550	No	350 kg	X				
2530213	4pt Sling bar 650 (XL) Air	650	No	350 kg	X				
1530101	4pt Sling bar 450 (M) for Mover 205/300	450	Yes	300 kg			X	X	
1530102	4pt Sling bar 550 (L) for Mover 205/300	550	Yes	300 kg			X	X	
1530103	4pt Sling bar 650 (XL) for Mover 205/300	650	Yes	300 kg			X	X	
1330160	4pt Sling bar 450 (M) for Partner 255	450	No	255 kg					X
1330161	4pt Sling bar 550 (L) for Partner 255	550	No	255 kg					X
1330162	4pt Sling bar 650 (XL) for Partner 255	650	No	255 kg					X
2510435	Bariatric Sling bar (pig tale hook)	550x550	No	500 kg	X				

Information sent to:
 Distribution list for Internal Product Information.
 Please forward this information to relevant staff in your organisation

2019-07-18 Ver: A

2/2