

Product change information

Concerning: Molift Air | Molift Nomad Cart

Molift Air 350

molift[®]

We are now able to deliver also the Air 350 with IRC Power Supply.

Item number 2610119.

If you need wall brackets, the item number is 2610123.

Reduction of ceiling hoist models from the factory

From August 1, 2020, we will reduce the number of versions, delivered from the factory.

The following models will remain:

- Air 200
- Air 205
- Air 300
- Air 350
- Air 500
- Nomad 205
- Nomad 255

All other SWL can still be provided to the end-user, but you will have to adjust it yourselves by using the Service Tool. Your certified technician will know how to do this.

For more details please contact [Trond Kleiven](#) or [Tore Bergengen](#).

All Air models are now NRTL-certified, including Air 200, Air 350 and Air 500.

Click to view the [certificate](#).

Nomad cart

We have made a new version of the Molift cart, that will replace the previous one.

New item number: 1200105

- The color is changed from red to grey
- The castors are bigger - the same used on Smart150

