

Product change information

Concerning: Molift Duo | Molift Rail System

Leveling Kit for MRS Duo

molift[®]

We now have a kit to retrofit adjustable feet on MRS Duo. This gives you the possibility to adjust the feet from the top. Please note that you have to drill with M8 from the top.

The item number for the leveling kit is 1600026.

For instruction, click on the link [BM01699 new adjustable foot](#)

Rail bracket for plaster wall

New item number 1109890 for mounting brackets directly on 2-layer plaster wall. SWL 350 kg.

We recommend using Hilti HTB-2-S M6x60, on plaster walls.

For instruction click on the link [000430](#)

