

Product change information

Concerning: Molift Hi-Trac | Molift Quick Raiser 205 | Molift Air 200/350/500

Hi-Trac delivered prior to 2007

molift[®]

We are not able to provide any more hand controls for Molift Hi-Trac delivered prior to 2007. Item numbers 1160014 & 1160012, are discontinued.

Reset periodic inspection indicator on Molift Quick Raiser 205 and Molift Air 200/350/500

To turn off the yellow service indicator during periodic inspection, you have to do a reset, in accordance with the illustrations in links below. Note that the procedure has to be repeated every 12 months.

[Reset periodic inspection QR205](#)

[Reset periodic inspection Air](#)

Battery replacement on PCB on Quick Raiser 205 and Air 200/350/500

Please look at the enclosed photos for instructions.

