


# Product change information

Concerning: Molift Product | Periodic Inspection label

## Molift Periodic Inspection label new version



Please note that Periodic Inspection label on all floor hoists and Air motors have been updated.



NEW: You will need to mark date for **next periodic inspection** on the sticker and not the date for service conducted.

Note that labels might have changed location on the product. For instance, on all Air motors the label is now found on the outside of the cover and not as currently on the inside.

The implementation of the labels will take place during the coming months.

Item no. label safety check: 1100303

Item no. yellow sticker: 1100306

2022-09-01

1/1