

# Product Information


**Concerning:** Molift Product | Groin Straps for RgoSling Ambulating Vest

## Product Change



We have made an assessment and concluded that the optional groin straps for RgoSling Ambulating Vest need to be shortened, to better fit and provide support to users. Both the padding "B" and the strap "C" have been shortened:

Part Number	Size	Old length padding (B)	New length padding (B)	Old length strap (C)	New length strap (C)
1721700	XXS-XS	48cm	26cm	50cm	30cm
1721720	S-M	54cm	38cm	50cm	35cm
1721740	L-XL	60cm	48cm </td <td>48cm</td> <td>40cm</td>	48cm	40cm



All orders placed now will be according to the new measurements.