

Product launch

Concerning: Molift Rail System | New upright wall support, Slim

molift[®]

Upright wall support, Slim

As of week 34 our new “Upright wall support, Slim” for Molift Rail System installations is ready for orders.



The main feature of the new product is that it is slimmer than the existing upright support (150 mm compared to 270 mm). This makes it a nice option when there is a wish the posts blend in well with the environment and not stand out too much from the wall.

The existing upright supports will still be available as well.

The installation guide (BM1901) will soon be updated to clearly present the new version alongside the existing ones.

Name:	Upright wall supp. Slim, H-rail
Article number:	1109265
Safe Working Load:	205 kg
Length:	2600 mm
Depth:	150 mm
Width:	60 mm

ITEM	SCBP 2023 SEK	SCBP 2023 DKK	SCBP 2023 EUR	SCBP 2023 NOK	SCBP 2023 USD	SCBP 2023 GBP
1109265	326.2	227.1	30.5	300.0	28.9	26.2

Upright wall support Slim, H-rail kit
Ceiling Hoist
Description
Floor Support Profile Slim (L=2,6m)
Holeplug 17mm
Topbracket wall floorsupport Slim
Multibolt MRS
Washer Large Plain DIN9021
Nut, Lock
Foot floorsuport slim

