

# Combination agreements

## When combining R82 product with other product

The combination must not change the intended purpose of either products or modify the products in such a way that compliance with the applicable requirements may be affected. The combiner of the products must ensure, that applicable requirements are fulfilled. Applicable requirements examples: alignment of user weight/mass of the combined seat and base, approval for transportation in vehicles of both products (ISO16840-4 for seat and ISO7176-19 for base), material safety (biocompatibility, flammability), mechanical safety (stability, squeezing hazards, inter connection), electrical safety. The maximum user weight for the combination is always the same as for the product in the combination that has the lowest maximum user weight.

Manufacturer	Product	R82 product
Active design Ltd.	CAPS II or Mini CAPS	R82 Wheelchairs with standard or Passive Locking interface from Active Design.
Alber GmbH	e-fix E35/E36	Combi Frame:x (w/ frame modification) Kudu (w/ 0° camber)
	e-motion M25	Combi Frame:x (w/ frame modification) Kudu
	e-motion M15	Kudu
BlueSky Designs	Mounting Systems Wheelchair hardware	Multi Frame
Decon Wheel AB	E-Drive	Combi Frame:x 20" Combi Frame:x 22" Combi Frame:x 24"
	E-Drive	Cougar 20" Cougar 22" Cougar 24
	E-Drive Adventus	Kudu 16"
	E-Drive E-Move Easywheel Adventus	Kudu 20" Kudu 22" Kudu 24" with 7" castors
	E-Drive	Multi Frame 22"
Degonda Rehab SA	Turbo Twist	x:panda size, 1-3 Panda Futura sizes 1-4
DIETZ Power	Sango	x:panda x:panda shape Customized seating
Etac AB	Cross 6 strap back	Kudu
	26100 elevating footrest	Kudu
Etac Supply Center AB	Cross / Cross XL / Cross VB / Cross XL VB Cross 5 / Cross 5XL / Cross 5VB / Cross 5XL VB Next / Next Comfort Elite Sting M100 / M100 TR	Head and Necks support: 9920092-83 / 9920093-83 9923492-83 / 9923493-83 9945492-83 / 9945493-83 9967492-83 / 9967493-83 87580 9975692-85 / 9975693-85 9975792-85 / 9975793-85 9971692-85 / 9971693-85 9971792-85 / 9971793-85 9905145-00 / 9906145-00 9905167-00 / 9906167-00
	Swift Mobil Tripp	Head support 3 (9905145-0, 9905167-0) Head support 4 (9906145-0, 9906167-0)

Manufacturer	Product	R82 product
	Neck support: 80209041 80209447 80209070	Flamingo Swan
	Lille Viking 2000 Lille Viking 2000+	Cougar Panther Combi Frame
	Tidy Hygiene Cover – 63087-63	Kudu
	Hemiplegia armrest	Heron
	Cross 6	R82 cushion set
Handicare B.V.	Hippo Alex Ibis XP Junior	x:panda
Huka Bikes	Diaz wheelchair bike	x:panda Panda Futura/colour x:panda shape
Janmek AB	Holder for oxygen cylinder	Kudu Panther
Ki Mobility LLC product	Flip Tyke	x:panda x:panda shape
Lewis Seating systems	Lewis Seating Orthosis	Combi Frame X
Otto Bock Mobility Solutions GmbH	Discovery Mobility Base size 1-4	x:panda
Permobil Benelux BV	C300	x:panda
Rehadapt Engineering GmbH & Co	Rehadapt frame clamps and wheelchair mounts	R82 A/S product
ROLTEC el-kørestole A/S	Boa2 Viper CT4	Customized seating x:panda
Scout Mobility	X10-FWD X10-RWD X10-MWD	x:panda and customized seating
Siesta Systems S.L.	Headpod Mounting Bracket HP408-R82	Head support
Stricker R&E Stricker Reha Entwicklungen GmbH	Stricker Handbikes (manual, hybrid) Stricker Electric Power Assists	Kudu
Supporttec & Orthototaal	Supporttec & Orthototaal zit-orthoses (soft & ABS)	Kudu
Tyromotion GmbH	PONDUS	Caribou Gazelle PS Meerkat Rabbit UP Toucan
Vela	Blues: 210-300	Customized seating x:panda
	Blues, Swing	x:panda
Veldink4kids	Veldink Adapt wheelchair base	x:panda Customized seating

# Covenant – Responsibility compiled medical device

## Transport in Motor Vehicles

Manufacturer	Product	R82 product
Bert Stevens Orthesebouw B.V	Contour ortese, Subframe	Combi Frame:x, Combi Frame
Perteon	Itura - Zitorthese Itura - Orthesezitting Itura -Ortheserug Bora- Flex orthse	Combi Frame:x, Combi Frame
RZ- groep	Classic Classic flex Contour	Combi Frame:x
Supporttec & Orthototaal	Orthototaal Soft flex orthese	Combi Frame:x
	Orthototaal ABS <ul style="list-style-type: none"> <li>- Flex</li> <li>- Gedeeld</li> <li>- Heel</li> <li>- Zitting</li> <li>- Rug</li> </ul> Orthototaal Soft (rond of Recht hout uitvoering) <ul style="list-style-type: none"> <li>- Flex</li> <li>- Gedeeld</li> <li>- Heel</li> <li>- Zitting</li> <li>- Rug</li> </ul> Supporttec Hardshell Flex Supporttec Hardshell Partial Supporttec Hardshell Supporttec Hardshell Seat Supporttec Hardshell Backrest Supporttec Foam Flex Supporttec Foam Partial Supporttec Foam Supporttec Foam Seat Supporttec Foam Backrest	Combi Frame:x, Combi Frame
Summit	Summit Alpine Plus orthosis Summit Alpine Ac-tive orthosis Summit Alpine Curve orthosis	Combi Frame:x
Veldink4kids	Adapt	x:panda