

Rogue2

Solid og aktiv rullestol for optimal ytelse!

Beskrivelse	Lev. artnr	HMS art.nr.	Setebredde (cm)	Setedybde (cm)	Setehøyde uten bytte av deler (cm)	Setehøyde med bytte av deler (cm)	Rygghøyde (cm)	Totalbredde (cm)	Totalvekt (kg)	Transportvekt (kg)
Rogue2 sb31 sd36	KI23751214	317551	31	36	46-48	34-52	28-51 (standard 30-36)	55,5	10,4	6,6
Rogue2 sb31 sd38	KI23751215	317552	31	38	46-48	34-52	28-51 (standard 30-36)	55,5	10,5	6,7
Rogue2 sb31 sd41	KI23751216	317553	31	41	46-48	34-52	28-51 (standard 30-36)	55,5	10,6	6,8
Rogue2 sb31 sd43	KI23751217	317554	31	43	46-48	34-52	28-51 (standard 30-36)	55,5	10,6	6,8
Rogue2 sb31 sd46	KI13751218	317555	31	46	46-48	34-52	28-51 (standard 30-36)	55,5	10,7	6,9
Rogue2 sb33 sd36	KI23751314	317556	33	36	46-48	34-52	28-51 (standard 30-36)	57,5	10,6	6,8
Rogue2 sb33 sd38	KI23751315	317557	33	38	46-48	34-52	28-51 (standard 30-36)	57,5	10,7	6,9
Rogue2 sb33 sd41	KI23751316	317558	33	41	46-48	34-52	28-51 (standard 30-36)	57,5	10,7	6,9
Rogue2 sb33 sd43	KI23751317	317559	33	43	46-48	34-52	28-51 (standard 30-36)	57,5	10,8	7
Rogue2 sb33 sd46	KI23751318	317560	33	46	46-48	34-52	28-51 (standard 30-36)	57,5	10,8	7
Rogue2 sb33 sd48	KI23751319	317561	33	48	46-48	34-52	28-51 (standard 30-36)	57,5	10,9	7,1
Rogue2 sb36 sd36	KI23751414	317562	36	36	46-48	34-52	28-51 (standard 30-36)	60,5	10,9	7,1
Rogue2 sb36 sd38	KI23751415	317563	36	38	46-48	34-52	28-51 (standard 30-36)	60,5	10,9	7,1
Rogue2 sb36 sd41	KI23751416	301422	36	41	46-48	34-52	28-51 (standard 30-36)	60,5	11	7,2
Rogue2 sb36 sd43	KI23751417	317564	36	43	46-48	34-52	28-51 (standard 30-36)	60,5	11	7,2
Rogue2 sb36 sd46	KI23751418	317565	36	46	46-48	34-52	28-51 (standard 30-36)	60,5	11	7,2
Rogue2 sb36 sd48	KI23751419	317566	36	48	46-48	34-52	28-51 (standard 30-36)	60,5	11,1	7,3
Rogue2 sb36 sd51	KI23751420	317567	36	51	46-48	34-52	28-51 (standard 30-36)	60,5	11,1	7,3
Rogue2 sb38 sd36	KI23751514	317568	38	36	46-48	34-52	28-51 (standard 30-36)	62,5	10,9	7,1
Rogue2 sb38 sd38	KI23751515	317569	38	38	46-48	34-52	28-51 (standard 30-36)	62,5	11	7,2
Rogue2 sb38 sd41	KI23751516	317570	38	41	46-48	34-52	28-51 (standard 30-36)	62,5	11	7,2
Rogue2 sb38 sd43	KI23751517	317571	38	43	46-48	34-52	28-51 (standard 30-36)	62,5	11,1	7,3
Rogue2 sb38 sd46	KI23751518	317572	38	46	46-48	34-52	28-51 (standard 30-36)	62,5	11,2	7,4
Rogue2 sb38 sd48	KI23751519	317573	38	48	46-48	34-52	28-51 (standard 30-36)	62,5	11,2	7,4
Rogue2 sb38 sd51	KI23751520	317574	38	51	46-48	34-52	28-51 (standard 30-36)	62,5	11,2	7,4
Rogue2 sb41 sd36	KI23751614	317575	41	36	46-48	34-52	28-51 (standard 30-36)	65,5	11	7,2
Rogue2 sb41 sd38	KI23751615	317576	41	38	46-48	34-52	28-51 (standard 30-36)	65,5	11	7,2
Rogue2 sb41 sd41	KI23751616	317577	41	41	46-48	34-52	28-51 (standard 30-36)	65,5	11,1	7,3
Rogue2 sb41 sd43	KI23751617	317578	41	43	46-48	34-52	28-51 (standard 30-36)	65,5	11,1	7,3
Rogue2 sb41 sd46	KI23751618	317579	41	46	46-48	34-52	28-51 (standard 30-36)	65,5	11,2	7,4
Rogue2 sb41 sd48	KI23751619	317580	41	48	46-48	34-52	28-51 (standard 30-36)	65,5	11,2	7,4
Rogue2 sb41 sd51	KI23751620	317581	41	51	46-48	34-52	28-51 (standard 30-36)	65,5	11,3	7,5
Rogue2 sb43 sd36	KI23751714	317582	43	36	46-48	34-52	28-51 (standard 30-36)	67,5	11,1	7,3
Rogue2 sb43 sd38	KI23751715	317583	43	38	46-48	34-52	28-51 (standard 30-36)	67,5	11,2	7,4
Rogue2 sb43 sd41	KI23751716	317584	43	41	46-48	34-52	28-51 (standard 30-36)	67,5	11,2	7,4
Rogue2 sb43 sd43	KI23751717	301295	43	43	46-48	34-52	28-51 (standard 30-36)	67,5	11,2	7,4
Rogue2 sb43 sd46	KI23751718	317585	43	46	46-48	34-52	28-51 (standard 30-36)	67,5	11,3	7,5
Rogue2 sb43 sd48	KI23751719	317586	43	48	46-48	34-52	28-51 (standard 30-36)	67,5	11,3	7,5

Rogue2

Solid og aktiv rullestol for optimal ytelse!

Beskrivelse	Lev. artnr	HMS art.nr.	Setebredde (cm)	Setedybde (cm)	Setehøyde uten bytte av deler (cm)	Setehøyde med bytte av deler (cm)	Rygghøyde (cm)	Totalbredde (cm)	Totalvekt (kg)	Transportvekt (kg)
Rogue2 sb43 sd51	KI23751720	317587	43	51	46-48	34-52	28-51 (standard 30-36)	67,5	11,4	7,6
Rogue2 sb46 sd38	KI23751815	317588	46	38	46-48	34-52	28-51 (standard 30-36)	70,5	11,2	7,4
Rogue2 sb46 sd41	KI23751816	317589	46	41	46-48	34-52	28-51 (standard 30-36)	70,5	11,2	7,4
Rogue2 sb46 sd43	KI23751817	317590	46	43	46-48	34-52	28-51 (standard 30-36)	70,5	11,3	7,5
Rogue2 sb46 sd46	KI23751818	317591	46	46	46-48	34-52	28-51 (standard 30-36)	70,5	11,3	7,5
Rogue2 sb46 sd48	KI23751819	317592	46	48	46-48	34-52	28-51 (standard 30-36)	70,5	11,4	7,6
Rogue2 sb46 sd51	KI23751820	317593	46	51	46-48	34-52	28-51 (standard 30-36)	70,5	11,4	7,6
Rogue2 sb48 sd41	KI23751916	317594	48	41	46-48	34-52	28-51 (standard 30-36)	72,5	11,2	7,4
Rogue2 sb48 sd43	KI23751917	317595	48	43	46-48	34-52	28-51 (standard 30-36)	72,5	11,3	7,5
Rogue2 sb48 sd46	KI23751918	317596	48	46	46-48	34-52	28-51 (standard 30-36)	72,5	11,3	7,5
Rogue2 sb48 sd48	KI23751919	317597	48	48	46-48	34-52	28-51 (standard 30-36)	72,5	11,3	7,5
Rogue2 sb48 sd51	KI23751920	317598	48	51	46-48	34-52	28-51 (standard 30-36)	72,5	11,4	7,6
Rogue2 sb51 sd41	KI23752016	317599	51	41	46-48	34-52	28-51 (standard 30-36)	75,5	11,3	7,5
Rogue2 sb51 sd43	KI23752017	317600	51	43	46-48	34-52	28-51 (standard 30-36)	75,5	11,3	7,5
Rogue2 sb51 sd46	KI23752018	317601	51	46	46-48	34-52	28-51 (standard 30-36)	75,5	11,4	7,6
Rogue2 sb51 sd48	KI23752019	317602	51	48	46-48	34-52	28-51 (standard 30-36)	75,5	11,4	7,6
Rogue2 sb51 sd51	KI23752020	317603	51	51	46-48	34-52	28-51 (standard 30-36)	75,5	11,5	7,7

