

# Convaid Carrot 3 Child Restraint System

## Canadian Retail Price List / Order Form

Account #:

Purchase Order #:

Date:

Contact:

Mark For:

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email:	Phone:

**Carrot Standard Features:**

- Recline Adjustment (90° - 180°)
- Head Support: 8.5" x 10.5" (H x L)
- Lower Anchors Connector System
- 5-Point Positioning Harness
- Trunk Positioning Inserts (Body Balance Pads)
- Seat Positioning Inserts
- Hip Positioning Inserts
- Chest Support: 9" Width

**Standard Models:**

<input type="checkbox"/> Carrot 3 CT3BLACK - CT3000-10	<b>Price</b> <b>\$2,943</b>
<input type="checkbox"/> Carrot 3 CT3BLUE - CT3000-20	<b>\$2,943</b>
<input type="checkbox"/> Carrot 3 CT3RED - CT3000-30	<b>\$2,943</b>

All prices in CAD

**Carrot 3 Recommended Item**

<input type="checkbox"/> Tether strap (Black)	<b>Price</b> <b>\$122</b>
<input type="checkbox"/> Top Securement Strap *	<b>\$73</b>
<input type="checkbox"/> Lower Securement Strap *	<b>\$89</b>

*\*For vehicles not equipped with anchor points.*

Black	Blue	Red
CT3TS0-10	CT3TS0-10	CT3TS0-10
CT3HSF-10	CT3HSF-10	CT3HSF-10
CT3LFS-10	CT3LFS-10	CT3LFS-10

**Pelvic Positioning**

<input type="checkbox"/> Rodded Seat Extension [Large/Curved (extends seat depth 7")] ^	<b>Price</b> <b>\$286</b>
<input type="checkbox"/> Rodded Seat Extension [Small/90° (extends seat depth 2")] ^	<b>\$233</b>
<input type="checkbox"/> Seat Extension [2" (Extends seat depth to 12")] ^ *	<b>\$128</b>
<input type="checkbox"/> Seat Extension [4" (Extends seat depth to 14")] ^ ^v	<b>\$155</b>
<input type="checkbox"/> Seat Angle Adjustment Wedge (Height: 4.25" / Angle: 20°)	<b>\$155</b>
<input type="checkbox"/> Pommel - Black *	<b>\$161</b>
<input type="checkbox"/> Anti-Adduction/Knee Separator Pad	<b>\$205</b>

*^ Growth specifications with listed extensions:**Seat depth: 9.8" - 22.3" / Backrest height: 11.8" - 22.8"**\* Rodded Seat Extension or Rodded Footrest required. 2" & 4" Seat Extensions cannot be used simultaneously when using Rodded Seat Extension.**v 7" Rodded Seat Extension required.*

Black	Blue	Red
CT3RSL-10	CT3RSL-20	CT3RSL-30
CT3RSS-10	CT3RSS-20	CT3RSS-30
CT3SE2-10	CT3SE2-20	CT3SE2-30
CT3SE4-10	CT3SE4-20	CT3SE4-30
CT3SA0-10	CT3SA0-20	CT3SA0-30
CT3PM0-10	CT3PM0-10	CT3PM0-10
CT3AAK-00	CT3AAK-00	CT3AAK-00

**Trunk Positioning**

<input type="checkbox"/> Backrest Extension - 2" (Extends back height to 13.8") *	<b>Price</b> <b>\$92</b>
<input type="checkbox"/> Backrest Extension - 4" (Extends back height to 15.8") *	<b>\$122</b>
<input type="checkbox"/> Shoulder Protector Wings (Extends back height 6")	<b>\$709</b>
<input type="checkbox"/> Firm Lateral Trunk Supports - Black^**	<b>\$251</b>

*\* The 2" & 4" Backrest Extensions cannot be used simultaneously.**^ Width adjustment from 10" - 12", and 2" of height adjustment.**\*\*Manufacturer installation required.*

Black	Blue	Red
CT3BE2-10	CT3BE2-20	CT3BE2-30
CT3BE4-10	CT3BE4-20	CT3BE4-30
CT3SP0-10	CT3SP0-20	CT3SP0-30
CT3LS0-10	CT3LS0-10	CT3LS0-10

**Foot & Leg Positioning**

<input type="checkbox"/> Rodded Footrest - Small [Gray (Adjusts from 10"-12")] *	<b>Price</b> <b>\$431</b>
<input type="checkbox"/> Rodded Footrest - Medium [Gray (Adjusts from 12"-16")] *	<b>\$431</b>
<input type="checkbox"/> Rodded Footrest - Large [Gray (Adjusts from 14"-18")] *	<b>\$431</b>

*\* Not available with Rodded Seat Extension*

Black	Blue	Red
CT3FRS-00	CT3FRS-00	CT3FRS-00
CT3FRM-00	CT3FRM-00	CT3FRM-00
CT3FRL-00	CT3FRL-00	CT3FRL-00

**Small Support Kit**

<input type="checkbox"/> Small Support Kit	<b>Price</b> <b>\$628</b>
• Small Head Support: 7.25" x 8.5" (H x L)	<b>\$0</b>
• Small Chest Support: Width 6.9"	<b>\$0</b>

Black	Blue	Red
CT3SSK-10	CT3SSK-20	CT3SSK-30

**Accessories**

<input type="checkbox"/> Support Tray (crash tested)	<b>Price</b> <b>\$239</b>
<input type="checkbox"/> 3D Body Balance Pads L/R Set	<b>\$206</b>
<input type="checkbox"/> Air-Mesh Cushions/Pads Set	<b>\$375</b>
<input type="checkbox"/> Pelvis Lock Harness (Small)	<b>\$342</b>
<input type="checkbox"/> Pelvis Lock Harness (Medium)	<b>\$342</b>
<input type="checkbox"/> Pelvis Lock Harness (Large)	<b>\$342</b>
<input type="checkbox"/> Anti-Escape Vest (Small)	<b>\$218</b>
<input type="checkbox"/> Anti-Escape Vest (Medium)	<b>\$218</b>
<input type="checkbox"/> Anti-Escape Vest (Large)	<b>\$218</b>
<input type="checkbox"/> Removeable Incontinence Cover	<b>\$274</b>
<input type="checkbox"/> Hinge Lock Belt	<b>\$57</b>
<input type="checkbox"/> FAA Fixing Strap to Aircraft Seat	<b>\$274</b>
<input type="checkbox"/> Booster Kit	<b>\$1,571</b>

Black	Blue	Red
CT3ST0-10	CT3ST0-20	CT3ST0-30
CT3BBP-00	CT3BBP-00	CT3BBP-00
CT3AMC-10	CT3AMC-10	CT3AMC-10
CT3PLS-10	CT3PLS-10	CT3PLS-10
CT3PLM-10	CT3PLM-10	CT3PLM-10
CT3PLL-10	CT3PLL-10	CT3PLL-10
CT3AVS-00	CT3AVS-00	CT3AVS-00
CT3AVM-00	CT3AVM-00	CT3AVM-00
CT3AVL-00	CT3AVL-00	CT3AVL-00
CT3RIR-00	CT3RIR-00	CT3RIR-00
CT3HLB-00	CT3HLB-00	CT3HLB-00
CT3FSA-00	CT3FSA-00	CT3FSA-00
CRT3BK-19	CRT3BK-29	CRT3BK-39

**Standard Shipping Rates: Ground \$69 • For expedited shipping rates, please call 1-844-876-6245**

Prices in CAD &amp; subject to change, ex works Torrance, California.

Effective: 1/1/2022

Convaid  
2830 California Street  
Torrance, CA 90503  
1 844-US-MOBILITY (844-876-6245)  
www.etac.us.com


OF0202-ENG-REV.05



Creating Possibilities

11/16/2021

CLIENT MEASUREMENTS:



Carrot 3 Base Model Specifications:

Seat Width:	11.5"
Seat Depth Adjustment:	9.8" - 11.3"
Backrest Height Adjustment:	11.8" - 12.8"
Head Support Height:	8.5"
Recline Adjustment:	90 - 180°
Weight:	19.1 lbs.

User Measurement Guide:

Minimum User Weight:	30 lbs.
Maximum User Weight:	108 lbs.
Minimum User Height:	37"
Maximum User Height:	60"
Recommended User Age:	Ages 3-15

Carrot 3 Specifications with Extensions:

Seat Depth Adjustment:	9.8" - 22.3"
Backrest Height Adjustment:	11.8" - 22.8"
Shoulder Wing Height:	6"
Shoulder Wing Width Adjustment:	15.25" - 19.25"

Base Model Measurements: \* ^

Seat Depth:	9.8" - 11.3"
Backrest Height:	11.8" - 12.8"
Head Support Height:	8.5"

\* Measurements do not include head support & growth extensions  
^ Base model includes 1.5" of seat depth growth (A) & 1" of backrest growth (B)

Growth Extensions:

	Inches
1a. Rodded seat extension - large	7"
1b. Rodded seat extension - small	2"
2. Seat extension - 4"	4"
3. Seat extension - 2"	2"
4. Backrest extension - 2"	2"
5. Backrest extension - 4"	4"
6. Shoulder protector wings	6"

Note: The same type of 2" & 4" extensions cannot be used simultaneously.

