

# Convoid Metro

## Canadian Retail Price List / Order Form

Account #: \_\_\_\_\_ Purchase Order #: \_\_\_\_\_ Date: \_\_\_\_\_  
 Contact: \_\_\_\_\_ Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email:	Phone:

**Metro Standard Features:**

- 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tire (ME12/14/16)
- Quick release Rolko rear wheels (ME12/14/16)
- 8" Front, 11" Rear Solid Tire (ME18)
- Wheel Locks
- Adjustable Swing-Away Footplates
- 30° Fixed Tilt
- Silver Powder Coated Frame
- 2-Pt Positioning Belt
- Made in the USA

<u>Standard Models:</u>	<u>Price</u>
<input type="checkbox"/> Metro ME12 - 901402	<b>\$1,443</b>
<input type="checkbox"/> Metro ME14 - 901400	<b>\$1,443</b>
<input type="checkbox"/> Metro ME16 - 901335	<b>\$1,443</b>
<input type="checkbox"/> Metro ME18 - 903127	<b>\$1,443</b>

**All prices in CAD**

<u>Upholstery Colors</u>	<u>Price</u>	<u>ME12</u>	<u>ME14</u>	<u>ME16</u>	<u>ME18</u>
<input type="checkbox"/> Cherry Red	<b>\$0</b>	903855	903855	903855	903855
<input type="checkbox"/> Royal Blue	<b>\$0</b>	903850	903850	903850	903850

<u>Transit Required Items</u>	<u>Price</u>	<u>ME12</u>	<u>ME14</u>	<u>ME16</u>	<u>ME18</u>
<input type="checkbox"/> Transit Option (Manufacturer Installation Required)	<b>\$275</b>	901403	901401	901336	903128
<input type="checkbox"/> Headrest Extension	<b>\$179</b>	903263	903264	903265	903266
<input type="checkbox"/> 3-Pt Positioning Belt	<b>\$68</b>	904017	904017	904017	904017
Foot Positioners - Select Size Below:		903415	903415	903415	903415
<input type="checkbox"/> Foot Positioners - Small 7" - 10" (pair)	<b>\$215</b>	903416	903416	903416	903416
<input type="checkbox"/> Foot Positioners - Medium 8" - 11" (pair)	<b>\$215</b>	903417	903417	903417	903417
<input type="checkbox"/> Foot Positioners - Large 9" - 12" (pair)	<b>\$215</b>				

<u>Transit Options</u>	<u>Price</u>	<u>ME12</u>	<u>ME14</u>	<u>ME16</u>	<u>ME18</u>
<input type="checkbox"/> Q'Straint Transit Lap Belt	<b>\$123</b>	903094	903094	903094	903094
<input type="checkbox"/> Sure-Lok Transit Lap Belt	<b>\$123</b>	903232	903232	903232	903232

<u>Tire Options</u>	<u>Price</u>	<u>ME12</u>	<u>ME14</u>	<u>ME16</u>	<u>ME18</u>
<input type="checkbox"/> 12.5" Rear Solid Knobby Tire	<b>\$272</b>	904117	904117	904117	
<input type="checkbox"/> 7.5" x 2" Front, 12.5" Rear Pneumatic Knobby Tire	<b>\$409</b>	904118	904118	904118	
<input type="checkbox"/> 8"x2" Front Ribbed, 12.5" Rear Solid Knobby Tire	<b>\$409</b>				903084
<input type="checkbox"/> 8"x2" Front Ribbed, 12.5" Rear Pneumatic Knobby Tire	<b>\$409</b>				903083

<u>Accessories</u>	<u>Price</u>	<u>ME12</u>	<u>ME14</u>	<u>ME16</u>	<u>ME18</u>
<input type="checkbox"/> 3-Pt Positioning Belt	<b>\$68</b>	904017	904017	904017	904017
<input type="checkbox"/> Wheelchair Travel Bag	<b>\$282</b>	902022	902022	902022	902022
<input type="checkbox"/> Under Seat Storage Basket - Fixed Tilt	<b>\$93</b>	904046	904047	904047	904048

The Metro Transit option includes four manufacturer installed wheelchair transportation anchors. The Metro transit model has been crash tested and conforms to RESNA WC4:2012 section 19 standards. Anthropomorphic Test Device weights are ME12: 66 lbs./30kg; ME14: 100 lbs./45.5kg/ME16 and ME18: 170 lbs./77 kg.

The Convoid Metro was tested equipped with a head support and wheelchair-anchored pelvic lap belt that comply with RESNA WC4:2012 section 18. The performance and safety of accessories and/or other seating options not tested by Convoid is unknown, and no warranty or claim about their performance is being, or can be, made.

**Standard Shipping Rates: Ground \$69 • For expedited shipping rates, please call 1-844-876-6245**

Prices in CAD & subject to change, ex works Torrance, California.

Effective: 1/1/2022

Convoid  
 2830 California Street  
 Torrance, CA 90503  
 1 844-US-MOBILITY (844-876-6245)  
 www.etac.us.com

