

Convoid Flyer - HCPCS E1233/E1161

Retail Price List / Order Form

Account #: _____

Purchase Order #: _____ Date: _____

Contact: _____

Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email:	Phone:

Flyer Base Standard Features:

- CG Tilt-in-Space (-5° to +40°)
- 30° back angle adjustment (80° to 110°)
- Quick Release Rear Wheels
- Single Foot Activated Wheel Lock
- HCPCS code E1161 is only valid on chairs 15" or greater
- Removable Front Riggings with Fixed Angled Footplate (Choose Option Below)
- Flip-up armrests - Long (Standard on FL15 and FL16)
- Choice of 3 Mobility Base colors
- One Piece Push Handle
- Choice of three back cane heights
- Made in the USA

Standard Models:

<input type="checkbox"/> Flyer FL12 - 905088	Price
<input type="checkbox"/> Flyer FL13 - 905089	\$2,919
<input type="checkbox"/> Flyer FL14 - 905090	\$2,919
<input type="checkbox"/> Flyer FL15 - 905091	\$2,919
<input type="checkbox"/> Flyer FL16 - 905092	\$2,919

Frame Colors

- Cherry Red
- Royal Blue
- Cosmic Silver
- Plum Purple

Price	HCPCS
N/C	
N/C	
N/C	
N/C	

FL12	FL13	FL14	FL15	FL16
905133	905133	905133	905133	905133
905134	905134	905134	905134	905134
905135	905135	905135	905135	905135
905461	905461	905461	905461	905461

Tire Selection

- 11.5" Rear Tire, 6" x 2" Front Caster (15" Front Seat to Floor Height)
Select One:
 - 11.5" Rear Fixed Wheels
 - 11.5" Quick Release Rear Wheels
- 16" Rear Tire, 7.5" x 2" Front Caster (17" Front Seat to Floor Height)
Select One:
 - 16" Rear Fixed Wheels
 - 16" Quick Release Rear Wheels

Price	HCPCS
N/C	
N/C	
N/C	
N/C	
N/C	
N/C	

FL12	FL13	FL14	FL15	FL16
905303	905304	905305	905306	905307
905384	905384	905384	905384	905384
905387	905387	905387	905387	905387
905059	905060	905061	905062	905063
905385	905385	905385	905385	905385
905388	905388	905388	905388	905388

Push Handle Option

- One-Piece Height Adjustable Push Handle Option (+ 9" height adjustable from back canes height) (Select One):**
- Outside Mount
 - Inside Mount

Price	HCPCS
N/C	
N/C	

905108	905109	905110	905111	905112
905098	905099	905100	905101	905102

Seat Depth Range Options (see chart on last page of order form)

- Position 1: 11"-15"
 - Position 2: 13"-17"
 - Position 3: 15"-19"
 - Position 4: 17"- 21"
- *Seat Depth is reduced by 1 1/2" when used with Solid Back System with Cushions and requires tools.

Price	HCPCS
N/C	
N/C	
N/C	
N/C	

FL12	FL13	FL14	FL15	FL16
905197	905197	905197	905197	905197
905198	905198	905198	905198	905198
905199	905199	905199	905199	905199
905200	905200	905200	905200	905200

Back Options

- Back Cane Fixed Heights* (7/8" diameter)- with Push Handles (Select One):**
- Short: 19" - (Supports back height range: 13" - 19")**** (Cannot be used with IV Pole)
 - Medium: 22" - (Supports back height range 13" - 22")***
 - Tall: 25" - (Supports back height range 13" - 25")****
- Height Adjustable Back Pans (4 Pieces)
 - Back Cushions & Cover (Includes Height Adjustable Back Pans) (13" - 25")****
- Cushion Accent Color (cushion is black, color is for piping along sides) Select One:**
- Royal Blue
 - Forest Green
 - Cherry Red
 - Sassy Purple
 - Princess Pink
 - Panther Black
 - Chocolate Brown
 - Aqua Blue
 - Apple Green

Price	HCPCS
N/C	
N/C	
N/C	
\$157	
\$547	E2291
N/C	
N/C	
N/C	
N/C	
N/C	
N/C	
N/C	
N/C	
N/C	
N/C	
N/C	
N/C	
N/C	

FL12	FL13	FL14	FL15	FL16
905054	905054	905054	905054	905054
905055	905055	905055	905055	905055
905056	905056	905056	905056	905056
905058	905058	905058	905058	905058
905093	905094	905095	905096	905097
903850	903850	903850	903850	903850
903857	903857	903857	903857	903857
903855	903855	903855	903855	903855
903851	903851	903851	903851	903851
903759	903759	903759	903759	903759
903854	903854	903854	903854	903854
904189	904189	904189	904189	904189
904753	904753	904753	904753	904753
904754	904754	904754	904754	904754

***Measured from pan to push handle bend.
****Height Adjustable Solid Back System w/Back Cushion HCPCS code only applies to pediatric sizes 12, 13, and 14.

Matrix Back Support Options

- Matrix E2 (FLwidth x 10"h) Color Black Meshtex¹
 - Matrix E2 (FLwidth x 12h) Color Black Meshtex¹
 - Matrix E2 (FLwidth x 14"h) Color Black Meshtex¹
 - Matrix E2 (FLwidth x 16"h) Color Black Meshtex¹
 - Matrix E2 (FLwidth x 18"h) Color Black Meshtex¹
 - Matrix E2 (FLwidth x 20"h) Color Black Meshtex¹
- ¹Must be ordered with Back Canes Only option.

Price	HCPCS
\$625	E2615
\$625	E2615
\$625	E2615
\$625	E2615
\$645	E2615
\$645	E2615

FL12	FL13	FL14	FL15	FL16
		E2S1410	E2S1510	E2S1610
		E2S1412	E2S1512	E2S1612
		E2S1414	E2S1514	E2S1614
		E2S1416	E2S1516	E2S1616
		E2S1418	E2S1518	E2S1618
		E2S1420	E2S1520	E2S1620

Matrx Back Support Options - Deep

- Matrx E2 Deep (FLwidth x 10"h) Color Black Meshtex¹
Matrx E2 Deep (FLwidth x 12"h) Color Black Meshtex¹
Matrx E2 Deep (FLwidth x 14"h) Color Black Meshtex¹
Matrx E2 Deep (FLwidth x 16"h) Color Black Meshtex¹
Matrx E2 Deep (FLwidth x 18"h) Color Black Meshtex¹
Matrx E2 Deep (FLwidth x 20"h) Color Black Meshtex¹

Price
\$755
\$755
\$755
\$755
\$775
\$775

HCPCS
E2620
E2620
E2620
E2620
E2620
E2620

Table with 4 columns: E2D1410, E2D1510, E2D1610, E2D1412, E2D1512, E2D1612, E2D1414, E2D1514, E2D1614, E2D1416, E2D1516, E2D1616, E2D1418, E2D1518, E2D1618, E2D1420, E2D1520, E2D1620

Matrx Back Support Options - Mini

- Matrx Mini ((FLwidth-1)"w x 7"h) Color Kidabra¹
Matrx Mini ((FLwidth-1)"w x 10"h) Color Kidabra¹
Matrx Mini ((FLwidth-1)"w x 12"h) Color Kidabra¹
Matrx Mini ((FLwidth-1)"w x 14"h) Color Kidabra¹
Matrx Mini ((FLwidth-1)"w x 7"h) Color Pink¹
Matrx Mini ((FLwidth-1)"w x 10h) Color Pink¹
Matrx Mini ((FLwidth-1)"w x 12"h) Color Pink¹
Matrx Mini ((FLwidth-1)"w x 14"h) Color Pink¹
Matrx Mini ((FLwidth-1)"w x 7"h) Color Blue¹
Matrx Mini ((FLwidth-1)"w x 10"h) Color Blue¹
Matrx Mini ((FLwidth-1)"w x 12"h) Color Blue¹
Matrx Mini ((FLwidth-1)"w x 14"h) Color Blue¹
Matrx Mini ((FLwidth-1)"w x 7"h) Color Black Startex¹
Matrx Mini ((FLwidth-1)"w x 10"h) Color Black Startex¹
Matrx Mini ((FLwidth-1)"w x 12"h) Color Black Startex¹
Matrx Mini ((FLwidth-1)"w x 14"h) Color Black Startex¹

Price
\$625
\$625
\$625
\$625
\$625
\$625
\$625
\$625
\$625
\$625
\$625
\$625
\$625
\$625
\$625
\$625

HCPCS
E2615
E2615
E2615
E2615
E2615
E2615
E2615
E2615
E2615
E2615
E2615
E2615
E2615
E2615
E2615
E2615

Table with 3 columns: MINI1107-K, MINI1207-K, MINI1307-K, MINI1110-K, MINI1210-K, MINI1310-K, MINI1112-K, MINI1212-K, MINI1312-K, MINI1114-K, MINI1214-K, MINI1314-K, MINI1107-P, MINI1207-P, MINI1307-P, MINI1110-P, MINI1210-P, MINI1310-P, MINI1112-P, MINI1212-P, MINI1312-P, MINI1114-P, MINI1214-P, MINI1314-P, MINI1107-B, MINI1207-B, MINI1307-B, MINI1110-B, MINI1210-B, MINI1310-B, MINI1112-B, MINI1212-B, MINI1312-B, MINI1114-B, MINI1214-B, MINI1314-B, MINI1107-S, MINI1207-S, MINI1307-S, MINI1110-S, MINI1210-S, MINI1310-S, MINI1112-S, MINI1212-S, MINI1312-S, MINI1114-S, MINI1214-S, MINI1314-S

Seat Options

- Planar Seat Cushion (Includes Solid Seat Pan)** (Same Color as Selected Back Cushion)
Seat Pan: 12"-16" Width, 4" Depth Growth ^^

Price
\$339
\$157

HCPCS
E2292
E2231

Table with 5 columns: FL12, FL13, FL14, FL15, FL16. Rows: 905083, 905084, 905085, 905086, 905087; 905057, 905057, 905057, 905057, 905057

Matrx Back Support Options - Libra

- Matrx Libra ((FLwidth)"w x 12"d) Color Black Startex²
Matrx Libra ((FLwidth)"w x 14"d) Color Black Startex²
Matrx Libra ((FLwidth)"w x 16"d) Color Black Startex²
Matrx Libra ((FLwidth)"w x 18"d) Color Black Startex²
Matrx Libra ((FLwidth)"w x 20"d) Color Black Startex²

Price
\$485
\$485
\$485
\$485
\$485

HCPCS
E2624
E2624
E2624
E2624
E2624

Table with 5 columns: FL12, FL13, FL14, FL15, FL16. Rows: LC1412, LC1414, LC1416, LC1516, LC1616, LC1518, LC1618, LC1520, LC1620

Matrx Back Support Options - Vi

- Matrx Kidabra-Vi ((FL width)x12d Color Kidabra²
Matrx Kidabra-Vi ((FL width)x14d Color Kidabra²
Matrx-Vi ((FL width)x15d Color Black Meshtex²
Matrx-Vi ((FL width)x16d Color Black Meshtex²
Matrx-Vi ((FL width)x17d Color Black Meshtex²
Matrx-Vi ((FLwidth)x18d Color Black Meshtex²
Matrx-Vi ((FLwidth)x20d Color Black Meshtex²

Price
\$425
\$425
\$425
\$425
\$425
\$425
\$425

HCPCS
E2607
E2607
E2607
E2607
E2607
E2607
E2607

Table with 5 columns: FL12, FL13, FL14, FL15, FL16. Rows: KA1212-Vi, KA1412-Vi, KA1214-Vi, KA1414-Vi, MA1415-Vi, MA1417-Vi, MA1616-Vi, MA1618-Vi, MA1620-Vi

Matrx Back Support Options - PS

- Matrx PS ((FLwidth)"w x 12"d) Color Kidabra²
Matrx PS ((FLwidth)"w x 14"d) Color Kidabra²
Matrx PS ((FLwidth)"w x 16"d) Color Kidabra²
Matrx PS ((FLwidth)"w x 12"d) Color Pink²
Matrx PS ((FLwidth)"w x 14"d) Color Pink²
Matrx PS ((FLwidth)"w x 16"d) Color Pink²
Matrx PS ((FLwidth)"w x 12"d) Color Blue²
Matrx PS ((FLwidth)"w x 14"d) Color Blue²
Matrx PS ((FLwidth)"w x 16"d) Color Blue²
Matrx PS ((FLwidth)"w x 12"d) Black Startex Reversible Cover²
Matrx PS ((FLwidth)"w x 14"d) Black Startex Reversible Cover²
Matrx PS ((FLwidth)"w x 16"d) Black Startex Reversible Cover²
Matrx PS ((FLwidth)"w x 18"d) Black Startex Reversible Cover²
Matrx PS ((FLwidth)"w x 20"d) Black Startex Reversible Cover²

Price
\$345
\$345
\$345
\$345
\$345
\$345
\$345
\$345
\$345
\$345
\$345
\$345
\$345
\$345

HCPCS
E2605
E2605
E2605
E2605
E2605
E2605
E2605
E2605
E2605
E2605
E2605
E2605
E2605
E2605

Table with 5 columns: FL12, FL13, FL14, FL15, FL16. Rows: PS1212-K, PS1312-K, PS1412-K, PS1214-K, PS1314-K, PS1414-K, PS1216-K, PS1316-K, PS1416-K, PS1212-P, PS1312-P, PS1412-P, PS1214-P, PS1314-P, PS1414-P, PS1216-P, PS1316-P, PS1416-P, PS1212-B, PS1312-B, PS1412-B, PS1214-B, PS1314-B, PS1414-B, PS1216-B, PS1316-B, PS1416-B, PS1212, PS1312, PS1412, PS1214, PS1314, PS1414, PS1216, PS1316, PS1416, PS1616, PS1618, PS1620

Matrx Back Support Options - PSP

- Matrx PSP ((FLwidth)"w x 12"d) Color Kidabra²
Matrx PSP ((FLwidth)"w x 14"d) Color Kidabra²
Matrx PSP ((FLwidth)"w x 16"d) Color Kidabra²
Matrx PSP ((FLwidth)"w x 12"d) Color Pink²
Matrx PSP ((FLwidth)"w x 14"d) Color Pink²
Matrx PSP ((FLwidth)"w x 16"d) Color Pink²
Matrx PSP ((FLwidth)"w x 12"d) Color Blue²
Matrx PSP ((FLwidth)"w x 14"d) Color Blue²
Matrx PSP ((FLwidth)"w x 16"d) Color Blue²
Matrx PSP ((FLwidth)"w x 12"d) Black Startex Reversible Cover²
Matrx PSP ((FLwidth)"w x 14"d) Black Startex Reversible Cover²
Matrx PSP ((FLwidth)"w x 16"d) Black Startex Reversible Cover²
Matrx PSP ((FLwidth)"w x 18"d) Black Startex Reversible Cover²
Matrx PSP ((FLwidth)"w x 20"d) Black Startex Reversible Cover²

Price
\$395
\$395
\$395
\$395
\$395
\$395
\$395
\$395
\$395
\$395
\$395
\$395
\$395
\$395

HCPCS
E2607
E2607
E2607
E2607
E2607
E2607
E2607
E2607
E2607
E2607
E2607
E2607
E2607
E2607

Table with 5 columns: FL12, FL13, FL14, FL15, FL16. Rows: PSP1212-K, PSP1312-K, PSP1412-K, PSP1214-K, PSP1314-K, PSP1414-K, PSP1216-K, PSP1316-K, PSP1416-K, PSP1212-P, PSP1312-P, PSP1412-P, PSP1214-P, PSP1314-P, PSP1414-P, PSP1216-P, PSP1316-P, PSP1416-P, PSP1212-B, PSP1312-B, PSP1412-B, PSP1214-B, PSP1314-B, PSP1414-B, PSP1216-B, PSP1316-B, PSP1416-B, PSP1212, PSP1312, PSP1412, PSP1214, PSP1314, PSP1414, PSP1216, PSP1316, PSP1416, PSP1616, PSP1618, PSP1620

Transit

- Transit Option (Manufacturer Installation Required)± ^
±See bottom of order form regarding the specific configuration in which the chair was tested for WC-19.
^Pelvic positioning belt required.

Price
\$209

HCPCS
K0108

Table with 5 columns: FL12, FL13, FL14, FL15, FL16. Row: 905082, 905082, 905082, 905082, 905082

Transit Options

- Q'Straint Transit Lap Belt
Sure-Lok Transit Lap Belt

Price
\$94
\$94

HCPCS
E0978
E0978

Table with 5 columns: FL12, FL13, FL14, FL15, FL16. Rows: 903094, 903094, 903094, 903094, 903094; 903232, 903232, 903232, 903232, 903232

Pelvic Positioning		Price	HCPCS	FL12	FL13	FL14	FL15	FL16
<input type="checkbox"/> 2-Point Pelvic Positioning Belt		\$47	E0978	903931	903931	903931	904559	904559
4-Point Pelvic Positioning Belt by R82 (Camlock Kit and Seat Pan Required)								
<input type="checkbox"/> Hip Width: 6.75" - 12.5" (Size XS)		\$140	K0108	905148	905148	905148	905148	905148
<input type="checkbox"/> Hip Width: 8.25" - 13.75" (Size S)		\$140	K0108	905149	905149	905149	905149	905149
<input type="checkbox"/> Hip Width: 9.75" - 16.25" (Size M)		\$140	K0108	905150	905150	905150	905150	905150
<input type="checkbox"/> Hip Width: 12.5" - 17.25" (Size L)		\$140	K0108	905151	905151	905151	905151	905151
<input type="checkbox"/> Hip Width: 16.25" - 20.5" (Size XL)		\$140	K0108	905475	905475	905475	905475	905475
<input type="checkbox"/> Camlock Kit for 4-Point Belt		N/C		905393	905393	905393	905393	905393
Left Width, Height Adj. Lateral Pelvic Support (Hip Guide)								
<input type="checkbox"/> Small - H: 4.5" - 5.5"		\$125	E0956	905320	905320	905320	905320	905320
<input type="checkbox"/> Large - H: 5.25" - 6.25"		\$125	E0956	905322	905322	905322	905322	905322
Right Width, Height Adj. Lateral Pelvic Support (Hip Guide)								
<input type="checkbox"/> Small - H: 4.5" - 5.5"		\$125	E0956	905321	905321	905321	905321	905321
<input type="checkbox"/> Large - H: 5.25" - 6.25"		\$125	E0956	905323	905323	905323	905323	905323
Medial Thigh Support								
<input type="checkbox"/> Medium (2.6" front width, 4.1" depth, 2" back width, 2.7" height)		\$169	E0957	905389	905389	905389	905389	905389
<input type="checkbox"/> Large (3.3" front width, 5" depth, 2.3" back width, 2.9" height)		\$169	E0957	905390	905390	905390	905390	905390
<input type="checkbox"/> X-Large (4" front width, 6" depth, 3" back width, 3" height)		\$169	E0957	905391	905391	905391	905391	905391

*Height is reduced by 2" when used with Cushion.

Foot & Leg Positioning		Price	HCPCS	FL12	FL13	FL14	FL15	FL16
Front Riggings (Footplates attach to the Front Riggings)								
<input type="checkbox"/> 65° Removable Front Riggings - Short (5" - 11" Seat Pan to Footrest Height)		N/C	E0990	905466	905466	905466	905466	905466
<input type="checkbox"/> 65° Removable Front Riggings - Long (5" - 13" Seat Pan to Footrest Height)		N/C	E0990	905465	905465	905465	905465	905465
<input type="checkbox"/> 90° Removable Front Riggings - (5.7" - 11" Seat Pan to Footrest Height)		N/C	E0990	905467	905467	905467	905467	905467
<input type="checkbox"/> Angle Adjustable Elevating Leg Rest (ELR) - Pair ***		\$260	K0195	905077	905077	905078	905078	905136
<input type="checkbox"/> Legrest Extension +2" - Pair (ELR Required)		\$57		904095	904095	904096	904096	904097
Footplates								
<input type="checkbox"/> One-Piece Angle Adjustable Footplates (Flip Up and Reversible) (ELR Required)		\$115	K0040	905160	905161	905162	905163	905164
<input type="checkbox"/> One-Piece Fixed Angle Footplates (Flip-up and Reversible Footplate) (ELR Required)		N/C		905072	905073	905074	905075	905076
<input type="checkbox"/> Individual Angle Adjustable Footplates - Pair (Compatible with std fixed front riggings or ELR)		\$186	K0040	905165	905165	905166	905166	905166
Foot Positioners (Pair)								
<input type="checkbox"/> Small 7" - 10"		\$163	K0108	903415	903415	903415	903415	903415
<input type="checkbox"/> Medium 8" - 11"		\$163	K0108	903416	903416	903416	903416	903416
<input type="checkbox"/> Large 9" - 12"		\$163	K0108	903417	903417	903417	903417	903417
Therafin Shoe Holders w/Padded Straps (Pair)								
<input type="checkbox"/> Small		\$196	K0108	904679	904679	904679	904679	904679
<input type="checkbox"/> Medium		\$196	K0108	904680	904680	904680	904680	904680
<input type="checkbox"/> Fully Padded Footbox (one-piece footplate required)		\$339	K0108	904954	904954	904955	904955	904956
<input type="checkbox"/> Calf Panel		\$64		904580	904580	904580	904580	904580
<input type="checkbox"/> Padded Elevating Legrest Cover (Pair)		\$15	K0108	904619	904619	904619	904619	904619

*2" Legrest Extension interferes with caster swivel when seat is in upright position. Seat must be tilted back at least 5° OR Legrests must be elevated to allow caster swivel.

**Seat to footplate range varies from 5.25" - 18" depending on chair size.

***Recommended for full recline.

Trunk Positioning		Price	HCPCS	FL12	FL13	FL14	FL15	FL16
Neoprene Vest by R82 (Camlock Kit and Back Pan Required)								
<input type="checkbox"/> LT Shoulder to RT Hip: 7.75" - 10.75" (Size XS)		\$209	E0960	905476	905476	905476	905476	905476
<input type="checkbox"/> LT Shoulder to RT Hip: 10.75" - 12.5" (Size S)		\$209	E0960	905152	905152	905152	905152	905152
<input type="checkbox"/> LT Shoulder to RT Hip: 12.5" - 16" (Size M)		\$209	E0960	905153	905153	905153	905153	905153
<input type="checkbox"/> LT Shoulder to RT Hip: 16" - 18.5" (Size L)		\$209	E0960	905154	905154	905154	905154	905154
<input type="checkbox"/> LT Shoulder to RT Hip: 18.5" - 20.5" (Size XL)		\$209	E0960	905480	905480	905480	905480	905480
<input type="checkbox"/> LT Shoulder to RT Hip: 20.5" - 23.5" (Size XXL)		\$209	E0960	905481	905481	905481	905481	905481
<input type="checkbox"/> Camlock Kit for Vest		N/C		905392	905392	905392	905392	905392
Chest Belt, Upholstered by R82 (Camlock Kit and Back Pan Required)								
<input type="checkbox"/> W: 11.25" x D:7.5" x H: 9.25-12" (Size 1)		\$282	E0960	905155	905155	905155	905155	905155
<input type="checkbox"/> W: 14" X D: 10.25" X H: 10.25-12.75" (Size 2)		\$282	E0960	905156	905156	905156	905156	905156
<input type="checkbox"/> W: 16.25" X D: 12.5" X H: 10.25-12.75" (Size 3)		\$282	E0960	905157	905157	905157	905157	905157
<input type="checkbox"/> W: 20.25" X D: 16" X H: 11.25-14.5" (Size 4)		\$282	E0960	905158	905158	905158	905158	905158
<input type="checkbox"/> Camlock Kit for Vest		N/C		905392	905392	905392	905392	905392
Cross Vest by R82 (Camlock Kit and Back Pan Required)								
<input type="checkbox"/> LT Shoulder to RT Hip: 7.75" - 12.5" (Size XS)		\$198	E0960	905176	905176	905176	905176	905176
<input type="checkbox"/> LT Shoulder to RT Hip: 10.25" - 13.5" (Size S)		\$198	E0960	905177	905177	905177	905177	905177
<input type="checkbox"/> LT Shoulder to RT Hip: 13.5" - 17.75" (Size M)		\$198	E0960	905178	905178	905178	905178	905178
<input type="checkbox"/> LT Shoulder to RT Hip: 17.75" - 20.75" (Size L)		\$198	E0960	905179	905179	905179	905179	905179
<input type="checkbox"/> LT Shoulder to RT Hip: 20.75" - 24.5" (Size XL)		\$198	E0960	905180	905180	905180	905180	905180
<input type="checkbox"/> Camlock Kit for Vest		N/C		905392	905392	905392	905392	905392
Swing-Away Lateral Support Bracket by R82 (Hrdw Req)								
<input type="checkbox"/> Swing-Away Lateral Support Bracket, Right (Ea.)		\$235	E01028	905338	905338	905338	905338	905338
<input type="checkbox"/> Swing-Away Lateral Support Bracket, Left (Ea.)		\$235	E01028	905337	905337	905337	905337	905337
Swing-Away Hardware								
<input type="checkbox"/> Bar, Straight ¾", Right (Ea.)		\$0		905327	905327	905327	905327	905327
<input type="checkbox"/> Bar, Straight ¾", Left (Ea.)		\$0		905327	905327	905327	905327	905327
<input type="checkbox"/> Bar, Offset, 1", Right (Ea.)		\$0		905328	905328	905328	905328	905328
<input type="checkbox"/> Bar, Offset, 1", Left (Ea.)		\$0		905328	905328	905328	905328	905328
<input type="checkbox"/> Bar, Offset 1¼", Right (Ea.)		\$0		905329	905329	905329	905329	905329
<input type="checkbox"/> Bar, Offset 1¼", Left (Ea.)		\$0		905329	905329	905329	905329	905329
<input type="checkbox"/> Bar, Offset 2½", Right (Ea.)		\$0		905330	905330	905330	905330	905330
<input type="checkbox"/> Bar, Offset 2½", Left (Ea.)		\$0		905330	905330	905330	905330	905330
Swing-Away Lateral Support Curved Cushion (Ea.)								
<input type="checkbox"/> Curved Cushion 3½" x 7½", Right (Ea.)		\$110	E0956	905331	905331	905331	905331	905331
<input type="checkbox"/> Curved Cushion 3½" x 7½", Left (Ea.)		\$110	E0956	905331	905331	905331	905331	905331
<input type="checkbox"/> Curved Cushion 2½" x 6¾", Right (Ea.)		\$110	E0956	905332	905332	905332	905332	905332
<input type="checkbox"/> Curved Cushion 2½" x 6¾", Left (Ea.)		\$110	E0956	905332	905332	905332	905332	905332
<input type="checkbox"/> Curved Cushion 2½" x 6", Right (Ea.)		\$110	E0956	905333	905333	905333	905333	905333
<input type="checkbox"/> Curved Cushion 2½" x 6", Left (Ea.)		\$110	E0956	905333	905333	905333	905333	905333
Swing-Away Lateral Support Standard Cushion (Ea.)								
<input type="checkbox"/> Standard Cushion 4" x 4", Right (Ea.)		\$110	E0956	905334	905334	905334	905334	905334
<input type="checkbox"/> Standard Cushion 4" x 4", Left (Ea.)		\$110	E0956	905334	905334	905334	905334	905334
<input type="checkbox"/> Standard Cushion 2½" x 6¾", Right (Ea.)		\$110	E0956	905335	905335	905335	905335	905335
<input type="checkbox"/> Standard Cushion 2½" x 6¾", Left (Ea.)		\$110	E0956	905335	905335	905335	905335	905335
<input type="checkbox"/> Protective Cushion for Laterals, Black (Pr.) (Laterals by R82 only)		\$96		905336	905336	905336	905336	905336
<input type="checkbox"/> 1" Thick Padded Lateral Support Cover (X Small only) (Contoured Laterals only)		\$52	E0960	905411	905411	905411	905411	905411

*Pad Size.

^Depth of Laterals is reduced by 1.5" with back upholstery.

**Requires height adjustable solid back system.

	<u>Customs*</u>	<u>Price</u>	<u>HCPCS</u>	<u>FL12</u>	<u>FL13</u>	<u>FL14</u>	<u>FL15</u>	<u>FL16</u>
<input type="checkbox"/> Two Piece Push Handle *30 Day Standard Delivery.		N/C		905311	905311	905311	905311	905311

***The Flyer Transit option includes four manufacturer installed black wheelchair transportation anchors. The Flyer Transit model has been crash tested and conforms to RESNA WC19 standards.

The Convaid Flyer Mobility Base was tested with the Flyer Seating Module equipped with a head support and wheelchair anchored pelvic lap belt. The performance and safety of accessories and/or other seating options not tested by Convaid is unknown, and no warranty or claim about their performance is or can being made.

Transit and Non Transit weight capacity: Back Cane Position 1: 85 lbs/39 kg and Back Cane Position 2-4: 170 lbs/77 kg.

Furthermore, when being transported in a vehicle it is Convaid's position that the wheelchair user transfer from the wheelchair into appropriate vehicle seating during transportation. Please note that at this time, the Department of Transportation has not officially approved any tie-down system for a user while in a wheelchair, therefore Convaid cannot recommend a transit system.

For more information on RESNA WC-19 and WC-20 standards please visit <http://www.recrwts.pitt.edu/WC19.html>.

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245


Prices in USD & subject to change, ex works Torrance, California.

Effective: 1/1/2022

Convaid
2830 California Street
Torrance, CA 90503
1 844-US-MOBILITY (844-876-6245)
www.etac.us.com



Seat Depth Adjustable Positions:



Position Option	1	2	3	4
Seating System Seat Depth Adjustment Range				
Back Canes	11" – 15"	13" – 17"	15" – 19"	17" – 21"
Solid Back*	9.5" – 13.5"	11.5" 15.5"	13.5" – 17.5"	15.5" – 19.5"
*Seat Depth is reduced by 1 1/2" when used with solid back seating system				

Footplate Dimensions:

Seat to Footplate Dimensions	Mounting	12" Flyer	13" Flyer	14" Flyer	15" Flyer	16" Flyer
One-Piece Angle Adjustable Footplates	Flip-Up	9in-14in	9in-14in	11in-16in	11in-16in	13in-18in
	Reversible	5.25in-11.25in	5.25in-11.25in	7.25in-13.25in	7.25in-13.25in	9.25in-15.25in
One-Piece Fixed Angle Footplates	Flip-Up	9in-14in	9in-14in	11in-16in	11in-16in	13in-18in
	Reversible	5.25in-11.25in	5.25in-11.25in	7.25in-13.25in	7.25in-13.25in	9.25in-15.25in
Individual Angle Adjustable Footplates - Pair	Flip-Up	7in-13in	7in-13in	9.5in-15.5in	9.5in-15.5in	12in-18in
	Reversible	4.5in-9.25in	4.25in-11in	4.25in-11in	4.25-11in	7.25in-13.25in

Height Adjustable Removable Armrest Length:

Recommended Height Adjustable Removable Armrest Length per Seat Depth Position				
	Seat Depth Position 1	Seat Depth Position 2	Seat Depth Position 3	Seat Depth Position 4
Height Adjustable Removable Armrest - Short 9.00" Length	●	●		
Height Adjustable Removable Armrest - Short 9.00" Length			●	●