

# R82 Caribou

## CAD Retail Price List/ Order Form


Account #: \_\_\_\_\_ Purchase Order No.: \_\_\_\_\_ Date: \_\_\_\_\_  
 Contact: \_\_\_\_\_ Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	
Phone:	City/State/Zip:
Email	Phone:

### US Configuration


Effective: 01.01.22 - Rev 7

- Caribou
- Top Pad/Shoulder Support
- Cross Straps/Back Support
- Chest Strap
- Head Support
- Extra Cushion
- Hip/Side Support
- Knee Support/Pommel with Strap
- Tray



- A (Total Height) \_\_\_\_\_
- B (Shoulder Height) \_\_\_\_\_
- C (Chest Height) \_\_\_\_\_
- D (Top of Pelvis) \_\_\_\_\_
- D1 (Buttock Height) \_\_\_\_\_
- D2 (Below Knee) \_\_\_\_\_
- E (Chest Width) \_\_\_\_\_
- F (Shoulder Width) \_\_\_\_\_
- G (Pelvis Width) \_\_\_\_\_
- Age \_\_\_\_\_
- Weight \_\_\_\_\_

### Measurements



	Size 1 Inches	Size 2 Inches	Size 3 Inches	Size 4 Inches
Frame width (G)	20"	22½"	25½"	29¼"
Frame length (H)	29¼"	36"	43¼"	43¼"
Height, vertical (I)	18-22½"	20¾-25¼"	26¾-31"	34¾-39"
Height, horizontal (I1)	27½"	27½"	31½"	31½"
Height from floor to frame (I2)	5"	5"	5"	5"
Basis plate width, top-bottom (J-J1)	9¾-12½"	11¾-15"	13¾-18"	15¾-21¾"
Basis plate height (K)	17"	19¾"	25½"	33½"
Angling	0-87°	0-87°	0-87°	0-87°
Foot plate angling	-10-10°	-10-10°	-10-10°	-10-10°
Weight, wo/power	45.2 lbs.	53 lbs.	67.2 lbs.	77.2 lbs.
Weight, w/power	49.6 lbs.	56.2 lbs.	71.7 lbs.	81.6 lbs.
Max load	88.2 lbs.	132.3 lbs.	176.4 lbs.	220½"
Approximate total user height	22-42"	32-52"	42-62"	52-72"

**Caribou Standard Features**



- Sizes 1, 2, 3 and 4
- 4" Locking Castors
- Contoured Top Plate
- Tilt: From Vertical to Horizontal
- Height and Angle Adjustable 2-Piece Footrest

Description	Size	Part No.	CAD Price
-------------	------	----------	-----------

**Important: For sizes 3 and 4 it is recommended to use the power version to ease the tilting; see below**

<input type="checkbox"/> Caribou Base w/ 2-Piece Footrest	1	860101-FP	\$3,551
<input type="checkbox"/> Caribou Base w/ 2-Piece Footrest	2	860102-FP	\$3,680
<input type="checkbox"/> Caribou Base w/ 2-Piece Footrest	3	860113-FP	\$4,077
<input type="checkbox"/> Caribou Base w/ 2-Piece Footrest	4	860114-FP	\$4,741

Description	Size	Part No.	CAD Price
-------------	------	----------	-----------

<input type="checkbox"/> Caribou Base with Power Actuator ●●	1	860101-ELFP	\$5,540
<input type="checkbox"/> Caribou Base with Power Actuator ●●	2	860102-ELFP	\$5,671
<input type="checkbox"/> Caribou Base with Power Actuator	3	860103-ELFP	\$6,067
<input type="checkbox"/> Caribou Base with Power Actuator	4	860104-ELFP	\$6,730

●●Custom order; Approx. lead time: 6-8 weeks. Non-returnable item (Sz 1 and 2)

Description	Measurements	Size	Part No.	CAD Price
-------------	--------------	------	----------	-----------

**Required, must select one:**

<input type="checkbox"/> Top Pad ♦	W1:6, W2:7, H:7¾"	1	861101-1	\$0
<input type="checkbox"/> Top Pad ♦	W1:7½, W2:9, H:9¾"	2	861101-2	\$0
<input type="checkbox"/> Top Pad ♦	W1:8¾, W2:10½, H:12½"	3	861101-3	\$0
<input type="checkbox"/> Top Pad ♦	W1:9¾, W2:12¼, H:15¾"	4	861101-4	\$0

♦Top Pad size must match Shoulder Support size

<input type="checkbox"/> Head Pad	W: 6", H: 7½"	1, 2	861116-1	\$657
<input type="checkbox"/> Head Pad	W: 8", H: 9½"	3, 4	861116-2	\$657

Description	Measurements	Size	Part No.	CAD Price
-------------	--------------	------	----------	-----------

<input type="checkbox"/> Extra Cushion (Ea.)	W:9¾" x H:3"	One Size	861106-1	\$299
<input type="checkbox"/> Extra Cushion (Ea.)	W:11¾ x H:4"	One Size	861106-2	\$299
<input type="checkbox"/> Extra Cushion (Ea.)	W:13¾ x H:4"	One Size	861106-3	\$326
<input type="checkbox"/> Extra Cushion (Ea.)	W:15¾ x H:4"	One Size	861106-4	\$326


Accessories	Description	Measurements	Size	Part No.	CAD Price
-------------	-------------	--------------	------	----------	-----------

■ XX = 60 (Black)














<input type="checkbox"/> Fix-Locks, Black, Bag of 6 Pcs ♦			One Size	862109-35	\$55
---	--	--	----------	-----------	------












♦Required for use of Vests













<input type="checkbox"/> Cross Vest (Ea.)	S-P: 7¾-12½"	XS	85720-XS	\$247
<input type="checkbox"/> Cross Vest (Ea.)	S-P 10¼-13½"	S	85720-S	\$247
<input type="checkbox"/> Cross Vest (Ea.)	S-P: 13½ - 17¾"	M	85720-M	\$247
<input type="checkbox"/> Cross Vest (Ea.)	S-P: 17¾-20¾"	L	85720-L	\$247
<input type="checkbox"/> Cross Vest (Ea.)	S-P: 20¾ - 24 ½"	XL	85720-XL	\$247



Accessories	Description	Measurements	Size	Part No.	CAD Price
	<input type="checkbox"/> Torso vest (Ea.)	S-P: 7 <sup>3</sup> / <sub>4</sub> -10 <sup>3</sup> / <sub>4</sub> "	XS	85705-XS	\$260
	<input type="checkbox"/> Torso vest (Ea.)	S-P: 10 <sup>3</sup> / <sub>4</sub> -12 <sup>1</sup> / <sub>2</sub> "	S	85705-S	\$260
	<input type="checkbox"/> Torso vest (Ea.)	S-P: 12 <sup>1</sup> / <sub>2</sub> -16"	M	85705-M	\$260
	<input type="checkbox"/> Torso vest (Ea.)	S-P: 16-18 <sup>1</sup> / <sub>2</sub> "	L	85705-L	\$267
	<input type="checkbox"/> Torso vest (Ea.)	S-P: 18 <sup>1</sup> / <sub>2</sub> -20 <sup>1</sup> / <sub>2</sub> "	XL	85705-XL	\$260
	<input type="checkbox"/> Torso vest (Ea.)	S-P: 20 <sup>1</sup> / <sub>2</sub> -23 <sup>1</sup> / <sub>2</sub> "	XXL	85705-XXL	\$260
	<input type="checkbox"/> Chest Belt, Upholstered (Ea.)		1	85428-1	\$351
	<input type="checkbox"/> Chest Belt, Upholstered (Ea.)		2	85428-2	\$351
	<input type="checkbox"/> Chest Belt, Upholstered (Ea.)		3	85428-3	\$351
<b>Other types of vests and belts available, please see Vest and Belts Order Form</b>					
	<input type="checkbox"/> Chest Strap (Ea.)	W: 16 <sup>1</sup> / <sub>2</sub> -24 <sup>3</sup> / <sub>4</sub> "	1, 2	862101-2	\$156
	<input type="checkbox"/> Chest Strap (Ea.)	W: 19-30 <sup>1</sup> / <sub>4</sub> "	3	862101-3	\$156
	<input type="checkbox"/> Chest Strap (Ea.)	W: 22 <sup>3</sup> / <sub>4</sub> -37 <sup>3</sup> / <sub>4</sub> "	4	862101-4	\$156
<b>Note: A chest strap is not included as standard, quantity of two straps is recommended.</b>					
	<input type="checkbox"/> Cross Straps (Pr.)	15-19 <sup>1</sup> / <sub>4</sub> "	1	862102-1	\$364
	<input type="checkbox"/> Cross Straps (Pr.)	17-21 <sup>1</sup> / <sub>2</sub> "	2	862102-2	\$364
	<input type="checkbox"/> Cross Straps (Pr.)	20 <sup>3</sup> / <sub>4</sub> -27 <sup>1</sup> / <sub>2</sub> "	3	862102-3	\$384
	<input type="checkbox"/> Cross Straps (Pr.) ●●	24 <sup>3</sup> / <sub>4</sub> -33 <sup>3</sup> / <sub>4</sub> "	4	862102-4	\$384
<b>●●Custom order; Approx. lead time: 6-8 weeks (Size 4)</b>					
	<input type="checkbox"/> Side Supports, Fixed, Type 7 (Pr.) Pad: 44.33 x 4.72"		1, 2	861102-10	\$326
	<input type="checkbox"/> Side Supports, Fixed, Type 5 (Pr.) Pad: 6.3 x 5.9"		3, 4	861102-20	\$390
<input type="checkbox"/> Adjustable Plate for Fixed Side Supports (Pr.)			One Size	861102-PL	\$182
	<input type="checkbox"/> Side Supports, size 1 (Pr.)	L: 5 <sup>1</sup> / <sub>2</sub> X H:5 <sup>1</sup> / <sub>2</sub> x W:3 <sup>3</sup> / <sub>4</sub> "	One Size	861102-1	\$370
	<input type="checkbox"/> Side Supports, size 2 (Pr.)	L: 6 <sup>1</sup> / <sub>2</sub> X H:7 x W:5"	One Size	861102-2	\$370
	<input type="checkbox"/> Side Supports, size 3 (Pr.)	L: 7 <sup>3</sup> / <sub>4</sub> H:7 x W:5"	One Size	861102-3	\$370
<b>Width between laterals on torso pad (min-max)</b>					
<b>Size 1: 4" to 9"</b>		<b>Size 2: 6<sup>1</sup>/<sub>2</sub>" to 14"</b>		<b>Size 3: 12" to 18"</b>	
	<input type="checkbox"/> Lateral Support Strap (Ea.)	10-12 <sup>1</sup> / <sub>2</sub> "	1	862103-1	\$66
	<input type="checkbox"/> Lateral Support Strap (Ea.)	13-16"	2	862103-2	\$66
	<input type="checkbox"/> Lateral Support Strap (Ea.)	13 <sup>1</sup> / <sub>2</sub> -19 <sup>3</sup> / <sub>4</sub> "	3	862103-3	\$66
	<input type="checkbox"/> Lateral Support Strap (Ea.)	17 <sup>3</sup> / <sub>4</sub> -23 <sup>1</sup> / <sub>2</sub> "	3	862103-4	\$66
	<input type="checkbox"/> Shoulder Support (Pr.) ♦	W:3, H:7 <sup>1</sup> / <sub>2</sub> , D:2"	1	861126-1	\$332
	<input type="checkbox"/> Shoulder Support (Pr.) ♦	W:3 <sup>1</sup> / <sub>4</sub> H:9 <sup>3</sup> / <sub>4</sub> , D:2"	2	861126-2	\$332
	<input type="checkbox"/> Shoulder Support (Pr.) ♦	W:4, H:12 <sup>1</sup> / <sub>2</sub> , D:2"	3	861126-3	\$332
	<input type="checkbox"/> Shoulder Support (Pr.) ♦	W:4 <sup>1</sup> / <sub>4</sub> H:15 <sup>1</sup> / <sub>4</sub> D:2"	4	861126-4	\$332
<b>♦ Shoulder Support size must match Top Pad size.</b>					
	<input type="checkbox"/> Upper Arm Support, Small		1, 2	861109-10	\$487
	<input type="checkbox"/> Upper Arm Support, Large		3, 4	861109-20	\$520
	<input type="checkbox"/> Wedges for Side Slots (Pr.) ●●		1	861107-1	\$21
	<input type="checkbox"/> Wedges for Side Slots (Pr.) ●●		2	861107-2	\$21
	<input type="checkbox"/> Wedges for Side Slots (Pr.) ●●		3	861107-3	\$21

Accessories	Description	Measurements	Size	Part No.	CAD Price
	<input type="checkbox"/> Back Support (Ea.) ●●		1, 2	861115-2	\$384
	<input type="checkbox"/> Back Support (Ea.)		3, 4	861115-3	\$390
	<b>●●Custom order; Approx. lead time: 6-8 weeks (Sizes 1 &amp; 2)</b>				
	<input type="checkbox"/> Turtle Bar System ◆	A: 10", B: 4-8½"	Long	98619-26	\$299
	<b>◆Required with all Head Supports</b>				
	<input type="checkbox"/> Anatomic Head Support ◆	A:5½ x B:5½ x C:9½"	Small	5099001	\$332
	<input type="checkbox"/> Head Support Cover (Required)		Small	9923-83	N/C
	<input type="checkbox"/> Anatomic Head Support ◆	A:6 x B:5½ x C:10"	Medium	5099002	\$332
	<input type="checkbox"/> Head Support Cover (Required)		Medium	9945-83	N/C
	<input type="checkbox"/> Anatomic Head Support ◆	A:6½ x B:5½ x C:12½"	Large	5099003	\$332
	<input type="checkbox"/> Head Support Cover (Required)		Large	9967-83	N/C
	<b>◆Turtle Bar System Short or Long Required; 98619-25 or 98619-26</b>				
<b>Other types of headrest available, please see Head Supports Order Form</b>					
Description	Measurements	Size	Part No.	CAD Price	
	<input type="checkbox"/> Pommel, Size 1 (Ea.)	W1:2,W2:2¾,H:6,D:3½"	One Size	861103-1	\$176
	<input type="checkbox"/> Pommel, Size 2 (Ea.)	W1:6,W2:2¾,H:6,D:6¾"	One Size	861103-2	\$176
	<input type="checkbox"/> Strap for Pommel, Sz 1 (Ea.) ●●	L:9.84 - 11.81"	One Size	862104-1	\$188
	<input type="checkbox"/> Strap for Pommel, Sz 2 (Ea.) ●●	L:10.63 - 18½"	One Size	862104-2	\$188
	<b>●●Custom order; Approx. lead time: 6-8 weeks</b>				
	<input type="checkbox"/> Knee Support, Fixed, sz 1		1, 2	861144-2	\$501
	<input type="checkbox"/> Knee Support, Hinged, Fixed, sz 1 ●●		1, 2	861144-20	\$696
	<input type="checkbox"/> Knee Support, Fixed, sz 2 **		3, 4	861144-3	\$501
	<input type="checkbox"/> Knee Support, Hinged, Fixed, sz 2 ●●		3, 4	861144-30	\$696
	<b>●●Custom order; Approx. lead time: 6-8 weeks</b>				
	<input type="checkbox"/> Knee Support, Swing-Away, sz 1		1, 2	861154-2	\$865
	<input type="checkbox"/> Knee Support, Hinged, Swing-Away, sz 1		1, 2	861154-20	\$1,060
	<input type="checkbox"/> Knee Support, Swing-Away, sz 2		3, 4	861154-3	\$865
	<input type="checkbox"/> Knee Support, Hinged, Swing-Away, sz 2		3, 4	861154-30	\$1,060
	<input type="checkbox"/> Formed Footplates (Pr.)	L: 8½"	1, 2	86147	\$215
	<input type="checkbox"/> Formed Footplates (Pr.)	L: 9¼"	3, 4	86148	\$215
	<input type="checkbox"/> Foot Straps, Plastic Lock (Pr.)	Max: 22", A: 13½-20½"	One Size	85096	\$73
	<input type="checkbox"/> Foot Straps, Velcro (Pr.)		One Size	85102	\$66
	<input type="checkbox"/> Heel Support for Parted Foot Support (Pr.)		1, 2	861251	\$73
	<input type="checkbox"/> Heel Support for Parted Foot Support (Pr.)		3, 4	861252	\$73
	<input type="checkbox"/> Upholstered Ankle Strap (Pr.)	A: 5½-7 x B:8½-10"	1	85426-1	\$240
	<input type="checkbox"/> Upholstered Ankle Strap (Pr.)	A: 6¼-7¾ x B:10-11¼"	2	85426-2	\$240
	<input type="checkbox"/> Upholstered Ankle Strap (Pr.)	A: 7-8½ x B:11-12½"	3	85426-3	\$240

Accessories	Description	Measurements	Size	Part No.	CAD Price
	<input type="checkbox"/> Neoprene Sandals (Pr.)		One Size	85696-1	\$169
<hr/>					
	<b>Description</b>	<b>Measurements</b>	<b>Size</b>	<b>Part No.</b>	<b>CAD Price</b>
<hr/>					
<b>Both Trays can be used in Prone and Supine Position</b>					
	<input type="checkbox"/> Tray, Plastic, with Bowl & Cover (Ea.)		1, 2	862130-2	\$729
	<input type="checkbox"/> Tray, Plastic, with Bowl & Cover (Ea.)		3, 4	862130-3	\$729
<hr/>					
<b>For Supine Use Only</b>					
	<input type="checkbox"/> Low Mount Kit for Plastic Tray ♦		One Size	862131-KIT	\$267
♦ <b>Plastic Tray required (862130-2 or 862130-3)</b>					
	<input type="checkbox"/> Armrest for Plastic Tray		One Size	861131	\$209
	<input type="checkbox"/> Toybar for Plastic Tray		One Size	861135	\$312
	<input type="checkbox"/> Grip Handle for Tray ••		One Size	8910750	\$527
•• <b>Custom order; Approx. lead time: 6-8 weeks</b>					
	<input type="checkbox"/> Helper's Handle ♦		One Size	861108-U	\$182
♦ <b>For moving and tilting the Caribou</b>					

Spare Parts	Description	Measurements	Size	Part No.	CAD Price
	<input type="checkbox"/> Power Kit, Complete (Ea.)		1, 2	861140	\$3,017
	<input type="checkbox"/> Power Kit, Complete (Ea.)		3	861150	\$2,066
	<input type="checkbox"/> Power Kit, Complete (Ea.)		4	861160	\$2,066
	<input type="checkbox"/> Plug for USA/Canada (Ea.)		One Size	2222166	\$34
	<input type="checkbox"/> Batteries for Caribou with Power (Pr.) ♦		One Size	2222012	\$397
♦ <b>Both batteries must be charged at the same time</b>					

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that all appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Prices in CAD are subject to change without notice

For Quotes and Orders please send completed Order Form to: [sales.us@etac.com](mailto:sales.us@etac.com) or fax: 1-800-760-0537

Etac  
13137 Bleinheim Lane  
Matthews, North Carolina 28105  
+1 844-US-MOBILITY (844-876-6245)  
[www.etac.us.com](http://www.etac.us.com)

OF5002-ENG-REV.07

  
Creating Possibilities