

Convoid Carrot 3 Child Restraint System

Canadian Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____

Contact: _____ ATP/RTS: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email:	

<ul style="list-style-type: none">● Recline Adjustment (90° - 180°)● Head Support: 8.5" x 10.5" (H x L)● Lower Anchors Connector System● 5-Point Positioning Harness● Tether Strap (Black) <p>Standard Models:</p> <table><tr><td><input type="checkbox"/> Carrot 3 CT3BLACK - CCT3000-100</td><td>Price</td></tr><tr><td><input type="checkbox"/> Carrot 3 CT3BLUE - CCT3000-200</td><td>\$3,368</td></tr><tr><td><input type="checkbox"/> Carrot 3 CT3RED - CCT3000-300</td><td>\$3,368</td></tr></table>	<input type="checkbox"/> Carrot 3 CT3BLACK - CCT3000-100	Price	<input type="checkbox"/> Carrot 3 CT3BLUE - CCT3000-200	\$3,368	<input type="checkbox"/> Carrot 3 CT3RED - CCT3000-300	\$3,368	<p>Carrot 3 Standard Features:</p> <ul style="list-style-type: none">● Trunk Positioning Inserts (Body Balance Pads)● Seat Positioning Inserts● Hip Positioning Inserts● Chest Support: 9" Width <p>All prices in CAD</p>
<input type="checkbox"/> Carrot 3 CT3BLACK - CCT3000-100	Price						
<input type="checkbox"/> Carrot 3 CT3BLUE - CCT3000-200	\$3,368						
<input type="checkbox"/> Carrot 3 CT3RED - CCT3000-300	\$3,368						

- Carrot 3 Recommended Item**
- | | |
|--|--------------|
| <input type="checkbox"/> Top Securement Strap* | Price |
| <input type="checkbox"/> Lower Securement Strap* | \$79 |
| | \$99 |
- *For vehicles not equipped with anchor points.*

- Trunk Positioning**
- | | |
|--|--------------|
| <input type="checkbox"/> Backrest Extension - 2" (Extends back height to 13.8")* | Price |
| <input type="checkbox"/> Backrest Extension - 4" (Extends back height to 15.8")* | \$103 |
| <input type="checkbox"/> Shoulder Protector Wings (Extends back height 6") | \$137 |
| <input type="checkbox"/> Firm Lateral Trunk Supports - Black** | \$795 |
| | \$282 |
- *The 2" & 4" Backrest Extensions cannot be used simultaneously.*
**Width adjustment from 10" - 12", and 2" of height adjustment.*
***Manufacturer installation required.*

- Pelvic Positioning**
- | | |
|---|--------------|
| <input type="checkbox"/> Rodded Seat Extension [Large/Curved (extends seat depth 7")]*^ | Price |
| <input type="checkbox"/> Rodded Seat Extension [Small/90° (extends seat depth 2")]*^ | \$313 |
| <input type="checkbox"/> Seat Extension [2" (Extends seat depth to 12")]*^ | \$255 |
| <input type="checkbox"/> Seat Extension [4" (Extends seat depth to 14")]*^v | \$141 |
| <input type="checkbox"/> Seat Angle Adjustment Wedge (Height: 4.25" / Angle: 20°) | \$171 |
| <input type="checkbox"/> Pommel - Black* (Medial Thigh Separator) | \$171 |
| <input type="checkbox"/> Anti-Adduction/Knee Separator Pad | \$179 |
| <input type="checkbox"/> Pelvis Lock Harness (Small) (20.5" - 33.75") | \$227 |
| <input type="checkbox"/> Pelvis Lock Harness (Medium) (23" - 39.5") | \$377 |
| <input type="checkbox"/> Pelvis Lock Harness (Large) (25" - 46") | \$377 |

**Growth specifications with listed extensions:*

Seat depth: 9.8" - 22.3" / Backrest height: 11.8" - 22.8".

**Rodded Seat Extension or Rodded Footrest required. 2" & 4" Seat Extensions cannot be used simultaneously when using Rodded Seat Extension.*
v7" Rodded Seat Extension required.

- Foot & Leg Positioning**
- | | |
|--|--------------|
| <input type="checkbox"/> Rodded Footrest - Small [Gray (Adjusts from 10"-12")]* | Price |
| <input type="checkbox"/> Rodded Footrest - Medium [Gray (Adjusts from 12"-16")]* | \$484 |
| <input type="checkbox"/> Rodded Footrest - Large [Gray (Adjusts from 14"-18")]* | \$484 |
- *Not available with Rodded Seat Extension.*

- Small Support Kit**
- | | |
|--|--------------|
| <input type="checkbox"/> Small Support Kit | Price |
| <ul style="list-style-type: none">● Small Head Support: 7.25" x 8.5" (H x L)● Small Chest Support: Width 6.9" | \$694 |

- Accessories**
- | | |
|---|----------------|
| <input type="checkbox"/> Upper Extremity Support (Tray) (Crash Tested) | Price |
| <input type="checkbox"/> 3D Body Balance Pads L/R Set | \$261 |
| <input type="checkbox"/> Air-Mesh Cushions/Pads Set | \$227 |
| <input type="checkbox"/> Anti-Escape Vest (Small) (Circum: 24.75" Height: 10.75") | \$415 |
| <input type="checkbox"/> Anti-Escape Vest (Medium) (Circum: 31.5" Height: 12") | \$241 |
| <input type="checkbox"/> Anti-Escape Vest (Large) (Circum: 37.5" Height: 13") | \$241 |
| <input type="checkbox"/> Removeable Incontinence Cover | \$302 |
| <input type="checkbox"/> Hinge Lock Belt | \$64 |
| <input type="checkbox"/> FAA Fixing Strap to Aircraft Seat | \$302 |
| <input type="checkbox"/> Booster Kit | \$1,736 |

Black	Blue	Red
CCT3HSF-10	CCT3HSF-10	CCT3HSF-10
CCT3LFS-10	CCT3LFS-10	CCT3LFS-10

Black	Blue	Red
CCT3BE2-10	CCT3BE2-20	CCT3BE2-30
CCT3BE4-10	CCT3BE4-20	CCT3BE4-30
CCT3SP0-10	CCT3SP0-20	CCT3SP0-30
CCT3LS0-10	CCT3LS0-10	CCT3LS0-10

Black	Blue	Red
CCT3RSL-10	CCT3RSL-20	CCT3RSL-30
CCT3RSS-10	CCT3RSS-20	CCT3RSS-30
CCT3SE2-10	CCT3SE2-20	CCT3SE2-30
CCT3SE4-10	CCT3SE4-20	CCT3SE4-30
CCT3SA0-10	CCT3SA0-20	CCT3SA0-30
CCT3PM0-10	CCT3PM0-10	CCT3PM0-10
CCT3AAK-00	CCT3AAK-00	CCT3AAK-00
CCT3PLS-10	CCT3PLS-10	CCT3PLS-10
CCT3PLM-10	CCT3PLM-10	CCT3PLM-10
CCT3PLL-10	CCT3PLL-10	CCT3PLL-10

Black	Blue	Red
CCT3FRS-00	CCT3FRS-00	CCT3FRS-00
CCT3FRM-00	CCT3FRM-00	CCT3FRM-00
CCT3FRL-00	CCT3FRL-00	CCT3FRL-00

Black	Blue	Red
CCT3SSK-10	CCT3SSK-20	CCT3SSK-30

Black	Blue	Red
CCT3ST0-10	CCT3ST0-20	CCT3ST0-30
CCT3BBP-00	CCT3BBP-00	CCT3BBP-00
CCT3AMC-10	CCT3AMC-10	CCT3AMC-10
CCT3AVS-00	CCT3AVS-00	CCT3AVS-00
CCT3AVM-00	CCT3AVM-00	CCT3AVM-00
CCT3AVL-00	CCT3AVL-00	CCT3AVL-00
CCT3RIR-00	CCT3RIR-00	CCT3RIR-00
CCT3HLB-00	CCT3HLB-00	CCT3HLB-00
CCT3FSA-00	CCT3FSA-00	CCT3FSA-00
CCRT3BK-199	CCRT3BK-299	CCRT3BK-399

Standard Shipping Rate: Ground \$71.95 • For expedited shipping rates, please call 1-844-876-6245 or email convoidsales.us@etac.com

Prices in CAD & subject to change, ex works Torrance, California.

Effective: 2/1/2025

Convoid
2830 California Street
Torrance, CA 90503
1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

OF0101-ENG-REV.01



Carrot 3 Base Model Specifications:

Seat Width:	11.5"
Seat Depth Adjustment:	9.8" - 11.3"
Backrest Height Adjustment:	11.8" - 12.8"
Head Support Height:	8.5"
Recline Adjustment:	90 - 180°
Weight:	19.1 lbs.

User Measurement Guide:

Minimum User Weight:	30 lbs.
Maximum User Weight:	108 lbs.
Minimum User Height:	37"
Maximum User Height:	60"
Recommended User Age:	Ages 3-15

Carrot 3 Specifications with Extensions:

Seat Depth Adjustment:	9.8" - 22.3"
Backrest Height Adjustment:	11.8" - 22.8"
Shoulder Wing Height:	6"
Shoulder Wing Width Adjustment:	15.25" - 19.25"

Base Model Measurements:**

Seat Depth:	9.8" - 11.3"
Backrest Height:	11.8" - 12.8"
Head Support Height:	8.5"

*Measurements do not include head support & growth extensions
*Base model includes 1.5" of seat depth growth (A) & 1" of backrest growth (B)

Growth Extensions:


- 1a. Rodded seat extension - large
- 1b. Rodded seat extension - small
- 2.a Seat extension - 4"
- 3.b Seat extension - 2"

- 4.a Backrest extension - 2"
- 5.a Backrest extension - 4"
- 6.a Shoulder protector wings


Note: The same type of 2" & 4" extensions cannot be used simultaneously.

Inches

- 7"
- 2"
- 4"
- 2"
- 2"
- 4"
- 6"



Convaid Carrot 3 Seat Depth Growth:



Position 1 (Inch)	9.8"	11.8"	13.8"	16.8"	18.8	20.8"
Position 2 (Inch)	11.3"	13.3"	15.3"	18.3"	20.3	22.3"

NOTE: 2" & 4" Seat Depth Extensions cannot be used simultaneously when using Rodded Seat Depth Extension

2830 California St., Torrance, CA 90503 • 1-844-876-6245, 310-618-0111 • Fax: 310-618-2166