

# Convoid Safari Tilt - WMK5CO040

## Canadian Retail Price List / Order Form

Account #: \_\_\_\_\_ Purchase Order #: \_\_\_\_\_ Date: \_\_\_\_\_

Contact: \_\_\_\_\_ ATP/RTS#: \_\_\_\_\_ Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email:	

### Safari Tilt Standard Features:

- 5° to 45° of Adjustable Tilt
- Silver Powder Coated Frame
- 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tire (SFT16)
- 6" Front, 11" Rear Solid Tires (SFT18)
- Quick Release Rear Wheels (SFT16)
- Foot Operated Wheel Locks
- 2-Point Pelvic Belt
- Heavy-Duty Back Upholstery
- One Piece, Non-Folding/Non-Elevating Footplate
- Support Cushion with Solid Seat Insert

### Standard Models:

- |   |                                |
|---|--------------------------------|
| <input type="checkbox"/> Safari SFT16 - C900526 | <b>Price</b><br><b>\$4,305</b> |
| <input type="checkbox"/> Safari SFT18 - C902432 | <b>\$4,305</b>                 |

All prices in CAD

### Upholstery Colors

- Panther Black  
 Royal Blue  
 Sassy Purple

**Price**  
N/C  
N/C  
N/C

SFT16	SFT18
C903854	C903854
C903850	C903850
C903851	C903851

### Frame Options

- Anatomic Back<sup>1</sup>  
 Recline Option<sup>1</sup>  
 Height Adjustable Flip-Up Armrests  
 Attendant Hand Brakes\* (Knobby Tires Required)  
 Curb Tipper

**Price**  
\$622  
\$858  
\$425  
\$439  
\$148

SFT16	SFT18
C901319	C902476
C901321	
C901322	C902477
C904064	C904064
C904063	C904063

<sup>1</sup>Manufacturer Installation Required.

### Transit Required Items

- Transit Anchor Brackets<sup>1\*</sup>  
 Headrest Extension\*  
 H-Harness with Padded Covers\* (Requires 3-Point Pelvic Belt)  
 3-Point Pelvic Belt\*  
 Elevating Legrests with Swing-Away Footplates\*  
 Foot Positioners\* - Select Size Below:  
 Small 7" - 10"  
 Medium 8" - 11" (Default if size is not selected)  
 Large 9" - 12"

**Price**  
\$373  
\$195  
\$148  
\$73  
\$502  
\$236  
\$236  
\$236

SFT16	SFT18
C900527	C902435
C901310	C902471
C903988	C903988
C904017	C904017
C901299	C902460
C903415	C903415
C903416	C903416
C903417	C903417

<sup>1</sup>Manufacturer Installation Required.

\*Crash Tested - WC19 and ISO 7176-19 Compliant

### Transit Options

- Q'Straint Transit Lap Belt

**Price**  
\$279

SFT16	SFT18
C903095	C903095

### Head Supports

- Padded Headwings  
 Occi Headwings  
 Headrest Extension

**Price**  
\$367  
\$344  
\$195

SFT16	SFT18
C905637	C905637
C905272	C905272
C901310	C902471

### Trunk Positioning

- H-Harness with Padded Covers (Requires 3-Point Pelvic Belt)  
 Full Torso Support Vest  
 Reducer Seat Insert (Reduces Seat Width 2" & Seat Depth 1")

**Price**  
\$148  
\$491  
\$577

SFT16	SFT18
C903988	C903988
C901308	C902469
C904110	C904110

### Pelvic Positioning

- 3-Point Pelvic Belt  
 Position Cushion, Medial Thigh Support/Anti Thrust<sup>1\*</sup>  
 Align Cushion, Lateral Pelvis/Lateral Thigh Support<sup>1\*</sup>  
 Incontinence Cushion Liner<sup>1</sup> (Cushion Required)

**Price**  
\$73  
\$703  
\$703  
\$109

SFT16	SFT18
C904017	C904017
C903590	C903591
C903594	C903595
C902324	C902324

<sup>1</sup>Manufacturer Installation Required.

\*Solid Seat Insert included with Cushion.

**Foot & Leg Positioning**

	<u>Price</u>
<input type="checkbox"/> Medial Thigh Support, Abductor	<b>\$301</b>
<input type="checkbox"/> Lateral Thigh Support, Adductor	<b>\$301</b>
<input type="checkbox"/> Elevating Legrests with Swing-Away Footplates	<b>\$502</b>
<input type="checkbox"/> 3" Legrest Extension*^	<b>\$103</b>
<input type="checkbox"/> Calf Panel	<b>\$93</b>
<input type="checkbox"/> Angle Adjustable Footplates^	<b>\$332</b>
Foot Positioners - Select Size Below:	
<input type="checkbox"/> Small 7" - 10"	<b>\$236</b>
<input type="checkbox"/> Medium 8" - 11"	<b>\$236</b>
<input type="checkbox"/> Large 9" - 12"	<b>\$236</b>
<input type="checkbox"/> Footplate Securement Strap	<b>\$47</b>

SFT16	SFT18
C901297	C902458
C901298	C902459
C901299	C902460
C901324	C902478
C904580	C904580
C903271	C903271
C903415	C903415
C903416	C903416
C903417	C903417
C903584	C903584

**\*3" Legrest Extension interferes with caster swivel when seat is in upright position.  
 Seat must be tilted back at least 5 degrees OR Legrests must be elevated to allow caster swivel.  
 ^Elevating Legrests with Swing-Away Footplates required.**

**Tire Options**

	<u>Price</u>
<input type="checkbox"/> 12.5" x 2" Rear Solid Knobby Tires	<b>\$301</b>
<input type="checkbox"/> 7.5" x 2" Front, 12.5" x 2" Rear Pneumatic Knobby Tires	<b>\$450</b>
<input type="checkbox"/> 12.5" x 2" Rear Solid Knobby Tires	<b>\$325</b>

SFT16	SFT18
C904117	
C904118	
	C903420

**Accessories**

	<u>Price</u>
<input type="checkbox"/> Headrest Cover (Canopy)	<b>\$370</b>
<input type="checkbox"/> Extended Headrest Cover (Canopy)	<b>\$404</b>
<input type="checkbox"/> Upper Extremity Support Surface Tray (Clear Tray)^	<b>\$181</b>
<input type="checkbox"/> Adj. Removeable Hardware (Required for Clear Tray)	<b>\$368</b>
<input type="checkbox"/> Solid Side Panels - (N/A w/ Recline Option)	<b>\$392</b>
<input type="checkbox"/> Wheelchair Travel Bag	<b>\$421</b>
<input type="checkbox"/> Utility Bag	<b>\$64</b>
<input type="checkbox"/> Mosquito Net (Canopy Required)	<b>\$227</b>

SFT16	SFT18
C901313	C902475
C904734	C904735
C903403	C903404
C904519	C904519
C904392	C904393
C902022	C902022
C904531	C904531
C904656	C904656

**^Flip-Up Armrests required.**

The Safari transit option includes four manufacturer installed wheelchair transportation anchors. The Safari transit model has been crash tested and conforms to RESNA WC4:2017 section 19 and ISO 7176-19:2022 standards.

*HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.*

**Standard Shipping Rates: Ground \$71.95 • For expedited shipping rates, please call 1-844-876-6245 or email [convaidsales.us@etac.com](mailto:convaidsales.us@etac.com)**

Prices in CAD & subject to change, ex works Torrance, California.

Convaid  
 2830 California Street  
 Torrance, CA 90503  
 1 844-US-MOBILITY (844-876-6245)  
[www.etac.us.com](http://www.etac.us.com)

Effective: 2/1/2025

OF0601-ENG-REV.01

