

Convoid Carrot 3 Child Restraint System

Retail Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____

Contact: _____ ATP/RTS#: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email:	

Carrot 3 Standard Features:

- Recline Adjustment (90° - 180°)
- Head Support: 8.5" x 10.5" (H x L)
- Lower Anchors Connector System
- 5-Point Positioning Harness
- Tether Strap (Black)
- Trunk Positioning Inserts (Body Balance Pads)
- Seat Positioning Inserts
- Hip Positioning Inserts
- Chest Support: 9" Width

Standard Models:

- | | |
|---|----------------|
| <input type="checkbox"/> Carrot 3 CT3 BLACK - CCT3000-100 | Price |
| <input type="checkbox"/> Carrot 3 CT3 BLUE - CCT3000-200 | \$2,434 |
| <input type="checkbox"/> Carrot 3 CT3 RED - CCT3000-300 | \$2,434 |

Carrot 3 Recommended Item

- | | |
|--|--------------|
| <input type="checkbox"/> Top Securement Strap* | Price |
| <input type="checkbox"/> Lower Securement Strap* | \$57 |
| \$71 | |

*For vehicles not equipped with anchor points.

Trunk Positioning

- | | |
|--|--------------|
| <input type="checkbox"/> Backrest Extension - 2" (Extends back height to 13.8")* | Price |
| <input type="checkbox"/> Backrest Extension - 4" (Extends back height to 15.8")* | \$74 |
| <input type="checkbox"/> Shoulder Protector Wings (Extends back height 6") | \$99 |
| <input type="checkbox"/> Firm Lateral Trunk Supports - Black*** | \$574 |
| | \$204 |

*The 2" & 4" Backrest Extensions cannot be used simultaneously.

**Width adjustment from 10" - 12", and 2" of height adjustment.

***Manufacturer installation required.

Pelvic Positioning

- | | |
|---|--------------|
| <input type="checkbox"/> Rodded Seat Extension [Large/Curved (extends seat depth 7")] | Price |
| <input type="checkbox"/> Rodded Seat Extension [Small/90° (extends seat depth 2")] | \$226 |
| <input type="checkbox"/> Seat Extension [2" (Extends seat depth to 12")] | \$185 |
| <input type="checkbox"/> Seat Extension [4" (Extends seat depth to 14")] | \$102 |
| <input type="checkbox"/> Seat Angle Adjustment Wedge (Height: 4.25" / Angle: 20°) | \$123 |
| <input type="checkbox"/> Pommel - Black* (Medial Thigh Separator) | \$123 |
| <input type="checkbox"/> Anti-Adduction/Knee Separator Pad | \$130 |
| <input type="checkbox"/> Pelvis Lock Harness (Small) (20.5" - 33.75") | \$164 |
| <input type="checkbox"/> Pelvis Lock Harness (Medium) (23" - 39.5") | \$272 |
| <input type="checkbox"/> Pelvis Lock Harness (Large) (25" - 46") | \$272 |

^Growth specifications with listed extensions:

Seat depth: 9.8" - 22.3" / Backrest height: 11.8" - 22.8".

*Rodded Seat Extension or Rodded Footrest required. 2" & 4" Seat Extensions cannot be used simultaneously when using Rodded Seat Extension.

v7" Rodded Seat Extension required.

Foot & Leg Positioning

- | | |
|---|--------------|
| <input type="checkbox"/> Rodded Footrest - Small [Gray (Adjusts from 10"-12")] | Price |
| <input type="checkbox"/> Rodded Footrest - Medium [Gray (Adjusts from 12"-16")] | \$350 |
| <input type="checkbox"/> Rodded Footrest - Large [Gray (Adjusts from 14"-18")] | \$350 |

*Not available with Rodded Seat Extension.

Small Support Kit

- | | |
|--|--------------|
| <input type="checkbox"/> Small Support Kit | Price |
| • Small Head Support: 7.25" x 8.5" (H x L) | \$502 |
| • Small Chest Support: Width 6.9" | N/C |
| | N/C |

Accessories

- | | |
|---|----------------|
| <input type="checkbox"/> Upper Extremity Support (Tray) (Crash Tested) | Price |
| <input type="checkbox"/> 3D Body Balance Pads L/R Set | \$189 |
| <input type="checkbox"/> Air-Mesh Cushions/Pads Set | \$164 |
| <input type="checkbox"/> Anti-Escape Vest (Small) (Circum: 24.75" Height: 10.75") | \$300 |
| <input type="checkbox"/> Anti-Escape Vest (Medium) (Circum: 31.5" Height: 12") | \$174 |
| <input type="checkbox"/> Anti-Escape Vest (Large) (Circum: 37.5" Height: 13") | \$174 |
| <input type="checkbox"/> Removeable Incontinence Cover | \$174 |
| <input type="checkbox"/> Hinge Lock Belt | \$218 |
| <input type="checkbox"/> FAA Fixing Strap to Aircraft Seat | \$46 |
| <input type="checkbox"/> Booster Kit | \$218 |
| | \$1,255 |

Black	Blue	Red
CCT3HSF-10	CCT3HSF-10	CCT3HSF-10
CCT3LFS-10	CCT3LFS-10	CCT3LFS-10

Black	Blue	Red
CCT3BE2-10	CCT3BE2-20	CCT3BE2-30
CCT3BE4-10	CCT3BE4-20	CCT3BE4-30
CCT3SP0-10	CCT3SP0-20	CCT3SP0-30
CCT3LS0-10	CCT3LS0-10	CCT3LS0-10

Black	Blue	Red
CCT3RSL-10	CCT3RSL-20	CCT3RSL-30
CCT3RSS-10	CCT3RSS-20	CCT3RSS-30
CCT3SE2-10	CCT3SE2-20	CCT3SE2-30
CCT3SE4-10	CCT3SE4-20	CCT3SE4-30
CCT3SA0-10	CCT3SA0-20	CCT3SA0-30
CCT3PM0-10	CCT3PM0-10	CCT3PM0-10
CCT3AAK-00	CCT3AAK-00	CCT3AAK-00
CCT3PLS-10	CCT3PLS-10	CCT3PLS-10
CCT3PLM-10	CCT3PLM-10	CCT3PLM-10
CCT3PLL-10	CCT3PLL-10	CCT3PLL-10

Black	Blue	Red
CCT3FRS-00	CCT3FRS-00	CCT3FRS-00
CCT3FRM-00	CCT3FRM-00	CCT3FRM-00
CCT3FRL-00	CCT3FRL-00	CCT3FRL-00

Black	Blue	Red
CCT3SSK-10	CCT3SSK-20	CCT3SSK-30

Black	Blue	Red
CCT3ST0-10	CCT3ST0-20	CCT3ST0-30
CCT3BBP-00	CCT3BBP-00	CCT3BBP-00
CCT3AMC-10	CCT3AMC-10	CCT3AMC-10
CCT3AVS-00	CCT3AVS-00	CCT3AVS-00
CCT3AVM-00	CCT3AVM-00	CCT3AVM-00
CCT3AVL-00	CCT3AVL-00	CCT3AVL-00
CCT3RIR-00	CCT3RIR-00	CCT3RIR-00
CCT3HLB-00	CCT3HLB-00	CCT3HLB-00
CCT3FSA-00	CCT3FSA-00	CCT3FSA-00
CCRT3BK-199	CCRT3BK-299	CCRT3BK-399

Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245 or email convoidsales.us@etac.com

Prices in USD & subject to change, ex works Torrance, California.

Effective: 2/1/2025

Convoid
2830 California Street
Torrance, CA 90503
1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

OF0100-ENG-REV.01


Creating Possibilities

Carrot 3 Base Model Specifications:

Seat Width:	11.5"
Seat Depth Adjustment:	9.8" - 11.3"
Backrest Height Adjustment:	11.8" - 12.8"
Head Support Height:	8.5"
Recline Adjustment:	90 - 180°
Weight:	19.1 lbs.

User Measurement Guide:

Minimum User Weight:	30 lbs.
Maximum User Weight:	108 lbs.
Minimum User Height:	37"
Maximum User Height:	60"
Recommended User Age:	Ages 3-15

Carrot 3 Specifications with Extensions:

Seat Depth Adjustment:	9.8" - 22.3"
Backrest Height Adjustment:	11.8" - 22.8"
Shoulder Wing Height:	6"
Shoulder Wing Width Adjustment:	15.25" - 19.25"

Base Model Measurements:**

Seat Depth:	9.8" - 11.3"
Backrest Height:	11.8" - 12.8"
Head Support Height:	8.5"


***Measurements do not include head support & growth extensions
 Base model includes 1.5" of seat depth growth (A) & 1" of backrest growth (B)

Growth Extensions:

- 1a. Rodded seat extension - large
- 1b. Rodded seat extension - small
- 2.a Seat extension - 4"
- 3.b Seat extension - 2"
- 4.a Backrest extension - 2"
- 5.a Backrest extension - 4"
- 6.a Shoulder protector wings


Inches

- 7"
- 2"
- 4"
- 2"
- 2"
- 4"
- 6"



Note: The same type of 2" & 4" extensions cannot be used simultaneously.

Convaid Carrot 3 Seat Depth Growth:



Position 1 (inch)	9.8"	11.8"	13.8"	16.8"	18.8"	20.8"
Position 2 (inch)	11.3"	13.3"	15.3"	18.3"	20.3"	22.3"

NOTE: 2" & 4" Seat Depth Extensions cannot be used simultaneously when using Rodded Seat Depth Extension

2830 California St., Torrance, CA 90503 • 1-844-876-6245, 310-618-0111 • Fax: 310-618-2166