

Convoid Coaster Scout- HCPCS E1236

Retail Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____

Contact: _____ ATP/RTS#: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email:	

Coaster Scout Standard Features:

- Three fixed tilt setting options 10°, 20°, or 30° (Adjustable with Tools)
- Four Fixed Recline Setting Options (Adjustable with Tools)
- Height and Tension Adjustable Back Positioner
- 2-Point Pelvic Belt
- Seat Depth Growth of 2"
- Standard Convoid Upholstery
- Zippered Pockets in Seat & Back Upholstery
- Removable, Height Adjustable Swing-Away Footplates
- One Piece Height Adjustable Push Handle

Standard Models:

- Coaster Scout CTR12 - C800123
- Coaster Scout CTR14 - C800143
- Coaster Scout CTR16 - C800163
- Coaster Scout CTR18 - C800183

Price

- \$2,821**
- \$2,821**
- \$2,821**
- \$2,821**

Transit Models Additional Features:

- Transit Anchor Brackets

Transit Included Models:

- Coaster Scout Transit CTR12T - C800124
- Coaster Scout Transit CTR14T - C800144
- Coaster Scout Transit CTR16T - C800164
- Coaster Scout Transit CTR18T - C800184

Price

- \$3,039**
- \$3,039**
- \$3,039**
- \$3,039**

Select Tilt and Recline Setting

Factory Default Setting if No Selection is Made: 10° Fixed Tilt and 105° Recline

- Fixed Tilt set at 10°**
Fixed Recline set at: 100 | 105 | 110 | 115
- Fixed Tilt set at 20°**
Fixed Recline set at: 90 | 95 | 100 | 105
- Fixed Tilt set at 30°**
Fixed Recline set at: 80 | 85 | 90 | 95

Upholstery Options

Colors:

- Panther Black C903854
- Royal Blue C903850
- Sassy Purple C903851
- Princess Pink C903759
- Forest Green C903857
- Cherry Red C903855

Foam Kit:

- Positioning Wedges (Foam Kit for Zippered Upholstery Pockets)

Price
\$102

HCPCS
K0108

CTR12	CTR14	CTR16	CTR18
C905567	C905568	C905569	C905570

Tilt and Recline Options



Seat Depth Extensions[^]

	Price	HCPCS
<input type="checkbox"/> 3" Seat Depth Extension	\$102	K0108

CTR12	CTR14	CTR16	CTR18
C905540	C905540	C905638	C905638

[^]The Coaster will ship with the +3 Seat depth extension that will only accomplish the +3 seat depth ranges (shown below).

Future Seat Depth Growth Kit^{*}

	Price	HCPCS
<input type="checkbox"/> 3" Seat Depth Extension with Fixed Front Riggings (Growth Kit) [^]	\$102	K0108
<input type="checkbox"/> 3" Seat Depth Extension with Elevating Legrest (Growth Kit) [^]	\$102	K0108


CTR12	CTR14	CTR16	CTR18
C905639	C905639	C905640	C905640
C905641	C905642	C905643	C905644

^{*}The Coaster will ship with standard seat depth and include a seat depth growth kit for future growth that will accomplish all seat depth ranges (shown below).

[^]Seat depth growth kit selection is dependent on the front rigging's selection.

Seat Depth Ranges

Seat Width	Standard	Option +3"
12"	10-12	13-15
14"	11-13	14-16
16"	12-15	15-18
18"	14.5-17.5	17.5-20.5



Fixed Front Riggings
Standard



Elevating Legrest
Option

Frame Options

	Price	HCPCS
<input type="checkbox"/> Caster Lock for 2" Wide Front Tires ¹	\$102	K0073
<input type="checkbox"/> Rear Anti-Tip Tubes (Cannot use with Curb Tipper)	\$125	E0971
<input type="checkbox"/> Curb Tipper (Cannot use with Rear Anti-Tip Tubes)	\$107	

CTR12	CTR14	CTR16	CTR18
C905526	C905526	C905526	C905526
C905527	C905527	C905527	C905527
C510737	C510737	C510737	C510737

¹Manufacturer Installation Required.

Transit Required Items^{}**

	Price	HCPCS
<input type="checkbox"/> Transit Anchor Brackets ¹	\$218	K0108
<input type="checkbox"/> Headrest Extension	\$141	E0966
<input type="checkbox"/> H-Harness with Padded Covers (Requires 3-Point Pelvic Belt)	\$107	E0960
<input type="checkbox"/> 3-Point Pelvic Belt	\$53	E0978
Foot Positioners - Select One Size Below:		
<input type="checkbox"/> Small 7" - 10"	\$170	K0108
<input type="checkbox"/> Medium 8" - 11" (Default if no size is selected)	\$170	K0108
<input type="checkbox"/> Large 9" - 12"	\$170	K0108

CTR12	CTR14	CTR16	CTR18
C905531	C905531	C905531	C905531
C902863	C900484	C900138	C905571
C903988	C903988	C903988	C903988
C904017	C904017	C904017	C904017
C903415	C903415	C903415	C903415
C903416	C903416	C903416	C903416
C903417	C903417	C903417	C903417

¹Manufacturer Installation Required.

^{*}Crash Tested - WC19 and ISO 7176-19 compliant

[^]The Transit Included Model excludes the price of the transit required items except the Transit Anchor Brackets.

Transit Options

	Price	HCPCS
<input type="checkbox"/> Q'Straint Transit Lap Belt	\$202	E0978

CTR12	CTR14	CTR16	CTR18
C903094	C903094	C903094	C903094

Tire Options

	Price	HCPCS
<input type="checkbox"/> 8" x 2" Front Ribbed, 12.5" x 2" Rear Solid Knobby Tires	\$325	

CTR12	CTR14	CTR16	CTR18
C905545	C905545	C905545	C905545

Head Supports

	Price	HCPCS
<input type="checkbox"/> Padded Headwing	\$265	E0955
<input type="checkbox"/> Occi Headwings	\$249	E0955
<input type="checkbox"/> Headrest Extension	\$141	E0966

CTR12	CTR14	CTR16	CTR18
C905635	C905636	C905637	C905637
C905271	C905272	C905272	C905272
C902863	C900484	C900138	C905571

Trunk Positioning

	Price	HCPCS
<input type="checkbox"/> H-Harness with Padded Covers (Requires 3-Point Pelvic Belt)	\$107	E0960
<input type="checkbox"/> Full Torso Support Vest	\$288	K0108
<input type="checkbox"/> Reducer Seat Insert (Reduces Seat Width 2" & Seat Depth 1") [*]	\$417	K0108
<input type="checkbox"/> Side Panels	\$231	K0108

CTR12	CTR14	CTR16	CTR18
C903988	C903988	C903988	C903988
C905598	C905599	C905600	C905601
C904211	C904110	C904110	C904110
C905552	C905553	C905554	C905555

^{*}Requires Side Panels

Pelvic Positioning

	Price	HCPCS
<input type="checkbox"/> 3-Point Pelvic Belt	\$53	E0978

CTR12	CTR14	CTR16	CTR18
C904017	C904017	C904017	C904017

Foot & Leg Positioning

	Price	HCPCS
<input type="checkbox"/> Medial Thigh Support, Abductor	\$176	E0957
<input type="checkbox"/> Lateral Thigh Support, Adductor	\$176	E0956
<input type="checkbox"/> Calf Panel [*]	\$67	E0995

CTR12	CTR14	CTR16	CTR18
C905576	C905577	C905578	C905579
C905580	C905581	C905582	C905583
C904580	C904580	C904580	C904580

Front Rigging Options: (80° Fixed Front Riggings Standard)

<input type="checkbox"/> Elevating Leg Rest (90° - 180°)	\$272	E0990
<input type="checkbox"/> Padded Elevating Legrest Cover	\$102	K0108

C905541	C905542	C905543	C905544
C904619	C904619	C904619	C904619

Footplate Options:

<input type="checkbox"/> Individual Angle Adjustable Footplates (Elevating Leg Rest Required)	\$173	K0040
<input type="checkbox"/> Elevating Leg Rest (90° - 180°)	\$272	E0990
<input type="checkbox"/> Angle Adjustable Footplates	\$173	K0040
<input type="checkbox"/> One Piece Angle Adjustable Flip-Up Footplate (-15° - +15°)	\$120	K0040
<input type="checkbox"/> Fully Padded Footbox (One Piece Footplate Required)	\$354	K0108
<input type="checkbox"/> Footplate Securement Strap	\$34	

C905624	C905625	C905626	C905627
C905541	C905542	C905543	C905544
C903445	C903268	C903268	C903269
C905557	C905558	C905559	
C904954	C904955	C904956	
C903584	C903584	C903584	C903584

Foot Positioners - Select One Size Below:

<input type="checkbox"/> Small 7" - 10"	\$170	K0108
<input type="checkbox"/> Medium 8" - 11"	\$170	K0108
<input type="checkbox"/> Large 9" - 12"	\$170	K0108

C903415	C903415	C903415	C903415
C903416	C903416	C903416	C903416
C903417	C903417	C903417	C903417

^{*}Calf Panel cannot be used when footplate extension tube is rotated, and in the highest position.

Medical Support Accessories

- Medical Necessity Storage Basket
- IV Pole (Height Adjustable and Collapsible)** *
- Adj. Swing Away Hardware (Required for IV Pole)

Price	HCPCS
\$321	
\$208	K0105
\$150	E1028

CTR12	CTR14	CTR16	CTR18
C905587	C905588	C905588	C905589
C904566	C904566	C904566	C904566
C905612	C905612	C905612	C905612

^^Cannot be used with canopy.
***Rear Anti-tip Tubes Required.**

Accessories

- Headrest Cover (Canopy)
- Extended Headrest Cover (Canopy)
- Under Seat Storage Basket
- Utility Bag
- Rain Cover (Canopy Required)
- Mosquito Net (Canopy Required)

Price	HCPCS
\$249	K0108
\$272	K0108
\$124	
\$46	
\$164	
\$164	

CTR12	CTR14	CTR16	CTR18
C905017	C905018	C905019	C905020
C905602	C905603	C905604	C905605
C905584	C905585	C905585	C905586
C904531	C904531	C904531	C904531
C905590	C905591	C905592	C905593
C905594	C905595	C905596	C905597

The Coaster transit option includes four manufacturer installed wheelchair transportation anchors. The Coaster transit model has been crash tested and conforms to RESNA WC4:2017 section 19 and ISO 7176-19:2022 standards. Anthropomorphic Test Device weights are CTR12: 66

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245 or email convoidsales.us@etac.com

Prices in USD & subject to change, ex works Torrance, California.

Effective: 2/1/2025

Convoid
 2830 California Street
 Torrance, CA 90503
 1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

