

Convaid Coaster - HCPCS E1236

Retail Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____

Contact: _____ ATP/RTS#: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email:	

Coaster Standard Features:

- Three fixed tilt setting options 10°, 20°, or 30° (Adjustable with Tools)
- Four Fixed Recline Setting Options (Adjustable with Tools)
- Tension Adjustable Back Straps
- 2-Point Pelvic Belt
- Seat Depth Growth of 2"
- Standard Convaid Upholstery
- Zippered Pockets in Seat & Back Upholstery
- Removable, Height Adjustable Swing-Away Footplates
- One Piece Height Adjustable Push Handle

Standard Models:

- Coaster CTR12 - C800121
- Coaster CTR14 - C800141
- Coaster CTR16 - C800161
- Coaster CTR18 - C800181

- 80° Fixed Front Riggings
- 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tires
- Quick Release Rolko Rear Wheels
- Single Foot Operated Wheel Locks
- Clamshell Compact Fold
- Closure Strap
- Saddle Bags
- Silver Powder Coated Frame
- Height Adjustable Flip Up Armrest

Transit Models Additional Features:

- Transit Anchor Brackets

Transit Included Models:

- Coaster Transit CTR12T - C800122
- Coaster Transit CTR14T - C800142
- Coaster Transit CTR16T - C800162
- Coaster Transit CTR18T - C800182

Price
\$2,005
\$2,005
\$2,005
\$2,005

Price
\$2,223
\$2,223
\$2,223
\$2,223

Select Tilt and Recline Setting

Factory Default Setting if No Selection is Made: 10° Fixed Tilt and 105° Recline

Fixed Tilt set at 10°

Fixed Recline set at: 100 | 105 | 110 | 115

Fixed Tilt set at 20°

Fixed Recline set at: 90 | 95 | 100 | 105

Fixed Tilt set at 30°

Fixed Recline set at: 80 | 85 | 90 | 95

Upholstery Options

Colors:

- Panther Black C903854
- Royal Blue C903850
- Sassy Purple C903851
- Princess Pink C903759
- Forest Green C903857
- Cherry Red C903855

Foam Kit:

- Positioning Wedges (Foam Kit for Zippered Upholstery Pockets)

Lumbar Support:

- Height and Tension Adjustable Back Positioner

Price
\$102
\$388

HCPCS
K0108
K0108

CTR12	CTR14	CTR16	CTR18
C905567	C905568	C905569	C905570
C905538	C905618	C905539	C905619

Tilt and Recline Options



- Seat Depth Extension**
- 3" Seat Depth Extension (Option)[^] **Price** **HCPCS**
\$102 **K0108**
- [^]The Coaster will ship with the +3 Seat depth extension that will only accomplish the +3 seat depth ranges (shown below).

CTR12	CTR14	CTR16	CTR18
C905540	C905540	C905638	C905638

- Future Seat Depth Growth Kit***
- 3" Seat Depth Extension with Fixed Front Riggings (Growth Kit)[^] **Price** **HCPCS**
\$102 **K0108**
- 3" Seat Depth Extension with Elevating Legrest (Growth Kit)[^] **Price** **HCPCS**
\$102 **K0108**

CTR12	CTR14	CTR16	CTR18
C905639	C905639	C905640	C905640
C905641	C905642	C905643	C905644

[^]The Coaster will ship with standard seat depth and include a seat depth growth kit for future growth that will accomplish all seat depth ranges (shown below).
[^]Seat depth growth kit selection is dependent on the front rigging's selection.

Seat Depth Ranges		
Seat Width	Standard	Option +3"
12"	10-12	13-15
14"	11-13	14-16
16"	12-15	15-18
18"	14.5-17.5	17.5-20.5



Fixed Front Riggings
Standard



Elevating Legrest
Option

- Frame Options**
- Caster Lock for 2" Wide Front Tires¹ **Price** **HCPCS**
\$102 **K0073**
- Rear Anti-Tip Tubes (Cannot use with Curb Tipper) **Price** **HCPCS**
\$125 **E0971**
- Curb Tipper (Cannot use with Rear Anti-Tip Tubes) **Price**
\$107
- ¹Manufacturer Installation Required.

CTR12	CTR14	CTR16	CTR18
C905526	C905526	C905526	C905526
C905527	C905527	C905527	C905527
C510737	C510737	C510737	C510737

- Transit Required Items****
- Transit Anchor Brackets¹ **Price** **HCPCS**
\$218 **K0108**
- Headrest Extension **Price** **HCPCS**
\$141 **E0966**
- H-Harness with Padded Covers (3-Point Pelvic Belt Required) **Price** **HCPCS**
\$107 **E0960**
- 3-Point Pelvic Belt **Price** **HCPCS**
\$53 **E0978**
- Foot Positioners - Select One Size Below:**
- Small 7" - 10" **Price** **HCPCS**
\$170 **K0108**
- Medium 8" - 11" (Default if no size is selected) **Price** **HCPCS**
\$170 **K0108**
- Large 9" - 12" **Price** **HCPCS**
\$170 **K0108**

CTR12	CTR14	CTR16	CTR18
C905531	C905531	C905531	C905531
C902863	C900484	C900138	C905571
C903988	C903988	C903988	C903988
C904017	C904017	C904017	C904017
C903415	C903415	C903415	C903415
C903416	C903416	C903416	C903416
C903417	C903417	C903417	C903417

¹Manufacturer Installation Required.

[^]Crash Tested - WC19 compliant

[^]The Transit Included Model excludes the price of the transit required items except the Transit Anchor Brackets.

- Transit Options**
- Q'Straint Transit Lap Belt **Price** **HCPCS**
\$202 **E0978**

CTR12	CTR14	CTR16	CTR18
C903094	C903094	C903094	C903094

- Tire Options**
- 7.5" x 2" Front Ribbed, 12.5" x 2" Rear Solid Knobby Tires **Price**
\$325
- 7.5" x 2" Front, 12.5" x 2" Rear Pneumatic Tires **Price**
\$325

CTR12	CTR14	CTR16	CTR18
C905535	C905535	C905535	C905535
C905561	C905561	C905561	C905561

- Head Supports**
- Padded Headwing **Price** **HCPCS**
\$265 **E0955**
- Occi Headwings **Price** **HCPCS**
\$249 **E0955**
- Headrest Extension **Price** **HCPCS**
\$141 **E0966**

CTR12	CTR14	CTR16	CTR18
C905635	C905636	C905637	C905637
C905271	C905272	C905272	C905272
C902863	C900484	C900138	C905571

- Trunk Positioning**
- H-Harness with Padded Covers (Requires 3-Point Pelvic Belt) **Price** **HCPCS**
\$107 **E0960**
- Full Torso Support Vest **Price** **HCPCS**
\$288 **K0108**
- Reducer Seat Insert (Reduces Seat Width 2" & Seat Depth 1")^{*} **Price** **HCPCS**
\$417 **K0108**
- Side Panels **Price** **HCPCS**
\$231 **K0108**

CTR12	CTR14	CTR16	CTR18
C903988	C903988	C903988	C903988
C905598	C905599	C905600	C905601
C904211	C904110	C904110	C904110
C905552	C905553	C905554	C905555

^{*}Requires Side Panels

- Pelvic Positioning**
- 3-Point Pelvic Belt **Price** **HCPCS**
\$53 **E0978**

CTR12	CTR14	CTR16	CTR18
C904017	C904017	C904017	C904017

- Foot & Leg Positioning**
- Medial Thigh Support, Abductor **Price** **HCPCS**
\$176 **E0957**
- Lateral Thigh Support, Adductor **Price** **HCPCS**
\$176 **E0956**
- Calf Panel^{*} **Price** **HCPCS**
\$67 **E0995**

CTR12	CTR14	CTR16	CTR18
C905576	C905577	C905578	C905579
C905580	C905581	C905582	C905583
C904580	C904580	C904580	C904580

- Front Rigging Options: (80° Fixed Front Riggings Standard)**
- Elevating Leg Rest (90° - 180°) **Price** **HCPCS**
\$272 **E0990**
- Padded Elevating Legrest Cover **Price** **HCPCS**
\$102 **K0108**

CTR12	CTR14	CTR16	CTR18
C905541	C905542	C905543	C905544
C904619	C904619	C904619	C904619

- Footplate Options:**
- Individual Angle Adjustable Footplates (Elevating Leg Rest Required) **Price** **HCPCS**
\$173 **K0040**
- Elevating Leg Rest (90° - 180°) **Price** **HCPCS**
\$272 **E0990**
- Angle Adjustable Footplates **Price** **HCPCS**
\$173 **K0040**
- One Piece Angle Adjustable Flip-Up Footplate (-15° - +15°) **Price** **HCPCS**
\$120 **K0040**
- Fully Padded Footbox (One Piece Footplate Required) **Price** **HCPCS**
\$354 **K0108**
- Footplate Securement Strap **Price**
\$34

CTR12	CTR14	CTR16	CTR18
C905624	C905625	C905626	C905627
C905541	C905542	C905543	C905544
C903445	C903268	C903268	C903269
C905557	C905558	C905559	
C904954	C904955	C904956	
C903584	C903584	C903584	C903584

- Foot Positioners - Select One Size Below:**
- Small 7" - 10" **Price** **HCPCS**
\$170 **K0108**
- Medium 8" - 11" **Price** **HCPCS**
\$170 **K0108**
- Large 9" - 12" **Price** **HCPCS**
\$170 **K0108**

CTR12	CTR14	CTR16	CTR18
C903415	C903415	C903415	C903415
C903416	C903416	C903416	C903416
C903417	C903417	C903417	C903417

^{*}Calf Panel cannot be used when footplate extension tube is rotated, and in the highest position.

Medical Support Accessories

- Medical Necessity Storage Basket
- IV Pole (Height Adjustable and Collapsible)**
- Adj. Swing Away Hardware (Required for IV Pole)

Price	HCPCS
\$321	
\$208	K0105
\$150	E1028

CTR12	CTR14	CTR16	CTR18
C905587	C905588	C905588	C905589
C904566	C904566	C904566	C904566
C905612	C905612	C905612	C905612

****Cannot be used with canopy.
*Rear Anti-tip Tubes Required.**

Accessories

- Headrest Cover (Canopy)
- Extended Headrest Cover (Canopy)
- Under Seat Storage Basket
- Utility Bag
- Rain Cover (Canopy Required)
- Mosquito Net (Canopy Required)

Price	HCPCS
\$249	K0108
\$272	K0108
\$124	
\$46	
\$164	
\$164	

CTR12	CTR14	CTR16	CTR18
C905017	C905018	C905019	C905020
C905602	C905603	C905604	C905605
C905584	C905585	C905585	C905586
C904531	C904531	C904531	C904531
C905590	C905591	C905592	C905593
C905594	C905595	C905596	C905597

The Coaster transit option includes four manufacturer installed wheelchair transportation anchors. The Coaster transit model has been crash tested and conforms to RESNA WC4:2017 section 19 and ISO 7176-19:2022 standards.

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245 or email convoidsales.us@etac.com

Prices in USD & subject to change, ex works Torrance, California.

Effective: 2/1/2025

Convoid
2830 California Street
Torrance, CA 90503
1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

OF0200-ENG-REV.01

