

Convoid Rodeo Tilt - HCPCS E1232/E1161

Retail Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____

Contact: _____ ATP/RTS#: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email:	

Rodeo Tilt Standard Features:

- 5° to 45° of Adjustable Tilt
- Silver Powder Coated Frame
- One Piece Height Adjustable Push Handle
- One Piece Angle Adjustable Footplate
- 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tire
- Quick Release Rolko Rear Tires
- Single Foot Activated Wheel Locks

- Zippered Pockets in Seat Bottom Upholstery
- 2-Point Pelvic Belt
- Seat Depth Growth of 5"
- Elevating Legrest
- Tool-less Upholstery Removal

Transit Models Additional Features:

- Transit Anchor Brackets

Standard Models:

- Rodeo RD10 - C904340
- Rodeo RD12 - C902689
- Rodeo RD14 - C902691
- Rodeo RD16 - C903134

Price

- \$2,691**
- \$2,691**
- \$2,691**
- \$2,691**

Transit Included Models:

- Rodeo Transit RD10T - C904860 **\$2,909**
- Rodeo Transit RD12T - C904861 **\$2,909**
- Rodeo Transit RD14T - C904862 **\$2,909**
- Rodeo Transit RD16T - C904863 **\$2,909**

Upholstery Colors

- Panther Black
- Royal Blue
- Sassy Purple
- Princess Pink
- Forest Green
- Cherry Red

Price

- N/C
- N/C
- N/C
- N/C
- N/C
- N/C

RD10	RD12	RD14	RD16
C903854	C903854	C903854	C903854
C903850	C903850	C903850	C903850
C903851	C903851	C903851	C903851
C903759	C903759	C903759	C903759
C903857	C903857	C903857	C903857
C903855	C903855	C903855	C903855

Frame Option

- Recline Option¹
- Height Adjustable Flip-Up Armrests
- Attendant Hand Brakes¹
- Caster Lock for 2" Wide Front Tires¹
- Curb Tipper (Cannot use with Rear Anti-Tip Tubes)
- Rear Anti-Tip Tubes (Cannot use with Curb Tipper)

Price

- \$503**
- \$249**
- \$257**
- \$102**
- \$107**
- \$125**

HCPCS

- E1225**
- E0973**
-
- K0073**
-
- E0971**

RD10	RD12	RD14	RD16
C904342	C904343	C904344	C904345
C904346	C904346	C904346	C904346
		C904370	C904371
C904650	C904650	C904650	C904650
C905310	C905310	C905310	C905310
C904649	C904649	C904649	C904649

Transit Required Items^{***}

- Transit Anchor Brackets^{1A}
- Headrest Extension^A
- H-Harness with Padded Covers (Requires 3-Point Pelvic Belt)^A
- 3-Point Pelvic Belt^A

Price

- \$218**
- \$141**
- \$107**
- \$53**

HCPCS

- K0108**
- E0966**
- E0960**
- E0978**

RD10	RD12	RD14	RD16
C904341	C903167	C903168	C903169
C903332	C902863	C900484	C900138
C903988	C903988	C903988	C903988
C904017	C904017	C904017	C904017
C903415	C903415	C903415	C903415
C903416	C903416	C903416	C903416
C903417	C903417	C903417	C903417

Foot Positioners - Select Size Below:

- Small 7" - 10"^A
- Medium 8" - 11"^A
- Large 9" - 12"^A

- \$170**
- \$170**
- \$170**

- K0108**
- K0108**
- K0108**

¹Manufacturer Installation Required.

^AThe Transit Included Model excludes the price of the transit required items except the Transit Anchor Brackets.

^{**}Crash Tested - WC19 and ISO 7176-19 compliant.

Transit Options

- Q'Straint Transit Lap Belt

Price

- \$202**

HCPCS

- E0978**

RD10	RD12	RD14	RD16
C903094	C903094	C903094	C903094

Head Supports

- Padded Headwings
- Occi Headwings
- Headrest Extension

Price

- \$265**
- \$249**
- \$141**

HCPCS

- E0955**
- E0955**
- E0966**

RD10	RD12	RD14	RD16
C905634	C905635	C905635	C905637
C905270	C905271	C905271	C905271
C903332	C902863	C900484	C900138

Trunk Positioning

- H-Harness with Padded Covers (Requires 3-Point Pelvic Belt)
- Right Adj Firm Lateral Trunk Sup - Contoured (Hrdw Required)**
 - X Small= 2.75"x 3.75"^{**}
 - Small = 3"x 5"
 - Medium = 3.75"x 5.5"
 - Large = 4.75"x 6"
- Left Adj Firm Lateral Trunk Sup - Contoured (Hrdw Required)**
 - X Small= 2.75"x 3.75"^{**}
 - Small = 3"x 5"
 - Medium = 3.75"x 5.5"
 - Large = 4.75"x 6"
- Right S/A Adj. Hardware for Lateral Trunk Support
- Left S/A Adj. Hardware for Lateral Trunk Support
- Full Torso Support Vest
- Reducer Seat Insert (Reduces Seat Width 2" & Seat Depth 1")

Price

- \$107**
- \$129**
- \$129**
- \$129**
- \$129**
- \$129**
- \$129**
- \$129**
- \$129**
- \$129**
- \$259**
- \$259**
- \$288**
- \$417**

HCPCS

- E0960**
- E0956**
- E0956**
- E0956**
- E0956**
- E0956**
- E0956**
- E0956**
- E0956**
- E0956**
- E1028**
- E1028**
- K0108**
- K0108**

RD10	RD12	RD14	RD16
C903988	C903988	C903988	C903988
C904385	C904385	C904385	C904385
	C904071	C904071	C904071
	C904072	C904072	C904072
	C904073	C904073	C904073
C904465	C904465	C904465	C904465
	C904466	C904466	C904466
	C904467	C904467	C904467
	C904468	C904468	C904468
C904362	C904362	C904362	C904362
C904363	C904363	C904363	C904363
C904364	C903261	C903260	C903259
C904211	C904211	C904110	C904110

^{**}RD10 requires X Small laterals.

- Pelvic Positioning**
- 3-Point Pelvic Belt
 - Support Cushion, General Use*1
 - Position Cushion, Medial Thigh Support/Anti Thrust**1
 - Align Cushion, Lateral Pelvis/Lateral Thigh Support**1
 - Incontinence Cushion Liner (Cushion Required)*

Price	HCPCS
\$53	E0978
\$262	E2601
\$508	E2607
\$508	E2607
\$79	E2619

RD10	RD12	RD14	RD16
C904017	C904017	C904017	C904017
C904347	C904350	C904353	C904356
C904348	C904351	C904354	C904357
C904349	C904352	C904355	C904358
C902324	C902324	C902324	C902324

***Manufacturer Installation Required.**
***Solid Seat Insert included with all Cushions.**

- Foot & Leg Positioning**
- Medial Thigh Support, Abductor
 - Lateral Thigh Support, Adductor
 - Padded Elevating Legrest Cover
 - 2" Legrest Extension**1
 - Calf Panel
 - Individual Angle Adjustable Footplates*1
 - Full Padded Footbox (One-Piece Footplate Required)

Price	HCPCS
\$176	E0957
\$176	E0956
\$102	K0108
\$114	
\$67	E0995
\$375	K0040
\$354	K0108

RD10	RD12	RD14	RD16
C904359	C902667	C902647	C903140
C904360	C902668	C902648	C903141
C904619	C904619	C904619	C904619
C904095	C904095	C904096	C904097
C904361	C904361	C904361	C904361
	C904957	C904958	C904959
C904953	C904954	C904955	C904956
C903415	C903415	C903415	C903415
C903416	C903416	C903416	C903416
C903417	C903417	C903417	C903417

Foot Positioners - Select Size Below:

- Small 7" - 10"^^
- Medium 8" - 11"^^
- Large 9" - 12"^^

***Manufacturer installation required if model was made prior to June 2019.**
****2" Legrest Extension interferes with caster swivel when seat is in upright position.**
Seat must be tilted back at least 5 degrees OR Legrests must be elevated to allow caster swivel.

- Tire Options**
- 12.5" x 2" Rear Solid Knobby Tires
 - 7.5" x 2" Front, 12.5" x 2" Rear Pneumatic Knobby Tires

Price
\$217
\$325

RD10	RD12	RD14	RD16
C904117	C904117	C904117	C904117
C904118	C904118	C904118	C904118

- Medical Support Accessories**
- Medical Necessity Storage Basket
 - IV Pole (height adjustable and collapsible)*
 - Adj. Swing Away Hardware (Required for IV Pole)
 - Oxygen Tank Bag

Price	HCPCS
\$321	
\$208	K0105
\$150	E1028
\$191	E2208

RD10	RD12	RD14	RD16
C904476	C904476	C904477	C904477
C904566	C904566	C904566	C904566
C905612	C905612	C905612	C905612
C902870	C902870	C902870	C902870

***Rear Anti-tip Tubes Required**

- Accessories**
- Headrest Cover (Canopy)
 - Extended Headrest Cover (Canopy)
 - Upper Extremity Support Surface Tray (Clear Tray)**^
 - Adj. Removable Hardware (Required for Clear Tray)
 - Under Seat Storage Basket
 - Utility Bag
 - Solid Side Panels
 - Mesh Side Panels
 - Rain Cover (Canopy Required)
 - Mosquito Net (Canopy Required)

Price	HCPCS
\$249	K0108
\$272	K0108
\$121	E0950
\$248	E1028
\$124	
\$46	
\$231	K0108
\$162	K0108
\$164	
\$164	

RD10	RD12	RD14	RD16
C904366	C902679	C902659	C903152
C904728	C904729	C904730	C904731
C904367	C903405	C903401	C903402
C903369	C903369	C903369	C903369
C904529	C904529	C904530	C904530
C904531	C904531	C904531	C904531
C904388	C902672	C902652	C903145
C904389	C902680	C902681	C903146
C904939	C904939	C904939	C904940
C904655	C904655	C904655	C904656

***Flip-Up Armrests and Hardware required.**
^Cannot be used during Transit.

The Rodeo Transit option includes four manufacturer installed black wheelchair transportation anchors. The Rodeo Transit model has been crash tested and conforms to RESNA WC19 and ISO 7176-19:2022 standards.

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245 or email convaidsales.us@etac.com
 2830 California St., Torrance, CA 90503 • 1-844-876-6245, 310-618-0111 • Fax: 310-618-2166

Prices in USD & subject to change, ex works Torrance, California.

Effective: 2/1/2025

Convaid
 2830 California Street
 Torrance, CA 90503
 1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

OF0500-ENG-REV.01

 Creating Possibilities