

Convoid Carrot 3 Booster Seat - Child Restraint System

Retail Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____
 Contact: _____ ATP/RTS: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email:	

Carrot 3 Booster Seat Standard Features:

- Carrot 3 Backrest
- Recline Adjustment (90° - 180°)
- Head Support: 8.5" x 10.5" (H x L)
- 4 Point Positioning Harness
- H-Harness Attaching Bracket
- Booster Seat Base
- Lower Anchors Connector System
- 2" Backrest Extension
- Shoulder Protector Wings
- Tether Strap (Black)

Standard Models:

	Price
<input type="checkbox"/> Carrot 3 Booster Seat CRT3BLACK - CCRT3BS-10	\$3,096
<input type="checkbox"/> Carrot 3 Booster Seat CRT3BLUE - CCRT3BS-20	\$3,096
<input type="checkbox"/> Carrot 3 Booster Seat CRT3RED - CCRT3BS-30	\$3,096

Carrot Booster Seat Recommended Item

- | | Price |
|--|-------------|
| <input type="checkbox"/> Top Securement Strap* | \$57 |
| <input type="checkbox"/> Lower Securement Strap* | \$71 |

**For vehicles not equipped with anchor points.*

Black	Blue	Red
CCT3HSF-10	CCT3HSF-10	CCT3HSF-10
CCT3LFS-10	CCT3LFS-10	CCT3LFS-10

Accessories

- | | Price |
|---|--------------|
| <input type="checkbox"/> Booster Support Tray (Crash Tested) | \$316 |
| <input type="checkbox"/> 3D Body Balance Pads L/R Set | \$164 |
| <input type="checkbox"/> Seat Angle Adjustment Wedge (Height: 4.25" / Angle: 20°) | \$123 |
| <input type="checkbox"/> Backrest Extension - 4" (Extends back height to 32.8") | \$99 |
| <input type="checkbox"/> Large Crotch Belt | \$259 |

Black	Blue	Red
CCRT3BT-19	CCRT3BT-29	CCRT3BT-39
CCT3BBP-00	CCT3BBP-00	CCT3BBP-00
CCT3SA0-10	CCT3SA0-20	CCT3SA0-30
CCT3BE4-10	CCT3BE4-20	CCT3BE4-30
CCT3CBX-00	CCT3CBX-00	CCT3CBX-00

Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245 or Email convoidsales.us@etac.com

Prices in USD & subject to change, ex works Torrance, California.


Effective: 1/20/2025

Convoid
 2830 California Street
 Torrance, CA 90503
 1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

OF1100-ENG-REV.00



CLIENT MEASUREMENTS:



- A** = _____ in.
- B** = _____ in.
- C** = _____ in.
- D** = _____ in.
- E** = _____ in.
- F** = _____ in.

A = Seat back height D = Shoulder width
 B = Seat depth E = Chest width
 C = Seat to floor F = Hip width

Carrot 3 Base Model Specifications:

Seat Width:	16.5"
Seat Depth:	15" - 16.5" *
Backrest Height Adjustment:	28" - 33"
Head Support Height:	8.5"
Backrest Height:	29" - 34"
Backrest Extension Height (Standard):	2"
Backrest Extension Height (Option):	4"
Shoulder Wing Height (Standard):	6"
Shoulder Wing Width Adjustment:	15.25" - 19.25"
Product Weight:	21 lbs.
Recline Adjustment:	90° - 180°
Product Overall Width:	20" - 22"

User Measurement Guide:

Minimum User Weight:	79 lbs.
Maximum User Weight:	165 lbs.
Minimum User Height:	55"
Maximum User Height:	70"

*Factory Setting

2830 California St., Torrance, CA 90503 • 1-844-876-6245, 310-618-0111 • Fax: 310-618-2166