

# Convoid Coaster - WM0002143

## Canadian Price List / Order Form

Account #: \_\_\_\_\_ Purchase Order #: \_\_\_\_\_ Date: \_\_\_\_\_  
 Contact: \_\_\_\_\_ ATP/RTS: \_\_\_\_\_ Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email:	

### Coaster Standard Features:

- Three fixed tilt setting options 10°, 20°, or 30° (Adjustable with Tools)
- Four Fixed Recline Setting Options (Adjustable with Tools)
- Tension Adjustable Back Straps
- 2-Point Pelvic Belt
- Seat Depth Growth
- Standard Convoid Upholstery
- Zippered Pockets in Seat & Back Upholstery
- Removable, Height Adjustable Swing-Away Footplates
- One Piece Height Adjustable Push Handle

### **Standard Models:**

- ☐ Coaster CTR12 - C800121
- ☐ Coaster CTR14 - C800141
- ☐ Coaster CTR16 - C800161
- ☐ Coaster CTR18 - C800181

### **Price**

**\$2,831**  
**\$2,831**  
**\$2,831**  
**\$2,831**

### Coaster Scout All Terrain Additional Features:

- Height and Tension Adjustable Back Positioner
  - Height Adjustable, Flip-Up Armrest
  - Attendant Hand Brake
  - 8" x 2" Front Ribbed, 12.5" Rear Pneumatic Tires
- Total for Scout Standard Features \$1,152**

**All Prices in CAD**

### **Coaster Scout Models:**

- ☐ Coaster Scout CTR12 - C800123
- ☐ Coaster Scout CTR14 - C800143
- ☐ Coaster Scout CTR16 - C800163
- ☐ Coaster Scout CTR18 - C800183

### **Price**

**\$3,983**  
**\$3,983**  
**\$3,983**  
**\$3,983**

### Select Tilt and Recline Setting

Factory Default Setting if No Selection is Made: 10° Fixed Tilt and 105° Recline

#### ☐ **Fixed Tilt set at 10°**

Fixed Recline set at: ☐ 100 | ☐ 105 | ☐ 110 | ☐ 115

#### ☐ **Fixed Tilt set at 20°**

Fixed Recline set at: ☐ 90 | ☐ 95 | ☐ 100 | ☐ 105

#### ☐ **Fixed Tilt set at 30°**

Fixed Recline set at: ☐ 80 | ☐ 85 | ☐ 90 | ☐ 95

### Upholstery Options

#### **Colors:**

- ☐ Panther Black C903854
- ☐ Royal Blue C903850
- ☐ Sassy Purple C903851
- ☐ Princess Pink C903759
- ☐ Forest Green C903857
- ☐ Cherry Red C903855

#### **Foam Kit:**

- ☐ Positioning Wedges (Foam Kit for Zippered Upholstery Pockets)

#### **Lumbar Support:**

- ☐ Height and Tension Adjustable Back Positioner

### **Price**

**\$144**

**\$548**

### Tilt and Recline Options



CTR12	CTR14	CTR16	CTR18
C905567	C905568	C905569	C905570
C905538	C905618	C905539	C905619

**Seat Depth Extension**

- ☐
- 3" Seat Depth Extension (Option)
- <sup>^</sup>

**Price****\$144**<sup>^</sup>The Coaster will ship with the +3 Seat depth extension that will only accomplish the +3 seat depth ranges (shown below).

CTR12	CTR14	CTR16	CTR18
C905540	C905540	C905638	C905638

**Future Seat Depth Growth Kit\***

- ☐
- 3" Seat Depth Extension with Fixed Front Riggings (Growth Kit)
- <sup>^</sup>
- 
- ☐
- 3" Seat Depth Extension with Elevating Legrest (Growth Kit)
- <sup>^</sup>

**Price****\$144****\$144**<sup>^</sup>The Coaster will ship with standard seat depth and include a seat depth growth kit for future growth that will accomplish all seat depth ranges (shown below).<sup>^</sup>Seat depth growth kit selection is dependent on the front rigging's selection.

CTR12	CTR14	CTR16	CTR18
C905639	C905639	C905640	C905640
C905641	C905642	C905643	C905644

**Seat Depth Ranges**

Seat Width	Standard	Option +3"
12"	10-12	13-15
14"	11-13	14-16
16"	12-15	15-18
18"	14.5-17.5	17.5-20.5

**Fixed Front Riggings**  
Standard**Elevating Legrest**  
Option**Frame Options**

- ☐
- Height Adjustable Flip-Up Armrests
- 
- ☐
- Caster Lock for 2" Wide Front Tires
- <sup>1</sup>
- 
- ☐
- Rear Anti-Tip Tubes (Cannot use with Curb Tipper)
- 
- ☐
- Curb Tipper (Cannot use with Rear Anti-Tip Tubes)

**Price****\$351****\$144****\$177****\$151**<sup>1</sup>Manufacturer Installation Required.

CTR12	CTR14	CTR16	CTR18
C905611	C905611	C905611	C905611
C905526	C905526	C905526	C905526
C905527	C905527	C905527	C905527
C510737	C510737	C510737	C510737

**Transit Required Items\***

- ☐
- Transit Anchor Brackets
- <sup>1</sup>
- 
- ☐
- Headrest Extension
- 
- ☐
- H-Harness with Padded Covers (3-Point Pelvic Belt Required)
- 
- ☐
- 3-Point Pelvic Belt

**Price****\$308****\$199****\$151****\$75****Foot Positioners - Select One Size Below:**

- ☐
- Small 7"-10"
- 
- ☐
- Medium 8"-11" (Default if no size is selected)
- 
- ☐
- Large 9"-12"

**\$240****\$240****\$240**<sup>1</sup>Manufacturer Installation Required.<sup>\*</sup>Crash Tested - WC19 compliant

CTR12	CTR14	CTR16	CTR18
C905531	C905531	C905531	C905531
C902863	C900484	C900138	C905571
C903988	C903988	C903988	C903988
C904017	C904017	C904017	C904017
C903415	C903415	C903415	C903415
C903416	C903416	C903416	C903416
C903417	C903417	C903417	C903417

**Transit Option**

- ☐
- Q'Straint Transit Lap Belt

**Price****\$285**

CTR12	CTR14	CTR16	CTR18
C903094	C903094	C903094	C903094

**Coaster Scout All Terrain Option****Coaster Scout Standard Features:**

- Height and Tension Adjustable Back Positioner
  - Height Adjustable, Flip-Up Armrest
  - Attendant Hand Brake
  - 8" x 2" Front Ribbed, 12.5" Rear Pneumatic Tires
- ☐
- 8" x 2" Front Ribbed, 12.5" Rear Solid Knobby Tires

**Price****\$468****\$351****\$88****\$245****\$459**

CTR12	CTR14	CTR16	CTR18
C905545	C905545	C905545	C905545

**Tire Options**

- ☐
- 7.5" x 2" Front Ribbed, 12.5" Rear Solid Knobby Tires
- 
- ☐
- 7.5" x 2" Front, 12.5" Rear Pneumatic Tires

**Price****\$459****\$459**

CTR12	CTR14	CTR16	CTR18
C905535	C905535	C905535	C905535
C905561	C905561	C905561	C905561

**Head Supports**

- ☐
- Padded Headwing
- 
- ☐
- Occi Headwings
- 
- ☐
- Headrest Extension

**Price****\$374****\$351****\$199**

CTR12	CTR14	CTR16	CTR18
C905635	C905636	C905637	C905637
C905271	C905272	C905272	C905272
C902863	C900484	C900138	C905571

**Trunk Positioning**

- ☐
- H-Harness with Padded Covers (Requires 3-Point Pelvic Belt)
- 
- ☐
- Full Torso Support Vest
- 
- ☐
- Reducer Seat Insert (Reduces Seat Width 2" & Seat Depth 1")
- <sup>\*</sup>
- 
- ☐
- Side Panels

**Price****\$151****\$406****\$589****\$325**<sup>\*</sup>Requires Side Panels

CTR12	CTR14	CTR16	CTR18
C903988	C903988	C903988	C903988
C905598	C905599	C905600	C905601
C904211	C904110	C904110	C904110
C905552	C905553	C905554	C905555

**Pelvic Positioning**

- ☐
- 3-Point Pelvic Belt

**Price****\$75**

CTR12	CTR14	CTR16	CTR18
C904017	C904017	C904017	C904017

**Foot & Leg Positioning**

- ☐
- Medial Thigh Support, Abductor
- 
- ☐
- Lateral Thigh Support, Adductor
- 
- ☐
- Calf Panel
- <sup>\*</sup>

**Price****\$249****\$249****\$95****Front Rigging Options: (80° Fixed Front Riggings Standard)**

- ☐
- Elevating Leg Rest (90° - 180°)
- 
- ☐
- Padded Elevating Legrest Cover

**\$384****\$144****Footplate Options:**

- ☐
- Individual Angle Adjustable Footplates (Elevating Leg Rest Required)
- ☐ Elevating Leg Rest (90° - 180°)
- ☐
- Angle Adjustable Footplates
- 
- ☐
- One Piece Angle Adjustable Flip-Up Footplate (-15° - +15°)
- 
- ☐
- Fully Padded Footbox (One Piece Footplate Required)
- 
- ☐
- Footplate Securement Strap

**\$245****\$384****\$245****\$170****\$500****\$48****Foot Positioners - Select One Size Below:**

- ☐
- Small 7" - 10"
- 
- ☐
- Medium 8" - 11"
- 
- ☐
- Large 9" - 12"

**\$240****\$240****\$240**<sup>\*</sup>Calf Panel cannot be used when footplate extension tube is rotated, and in the highest position.

CTR12	CTR14	CTR16	CTR18
C905576	C905577	C905578	C905579
C905580	C905581	C905582	C905583
C904580	C904580	C904580	C904580
C905541	C905542	C905543	C905544
C904619	C904619	C904619	C904619
C905624	C905625	C905626	C905627
C905541	C905542	C905543	C905544
C903445	C903268	C903268	C903269
C905557	C905558	C905559	
C904954	C904955	C904956	
C903584	C903584	C903584	C903584
C903415	C903415	C903415	C903415
C903416	C903416	C903416	C903416
C903417	C903417	C903417	C903417

**Medical Support Accessories**

- ☐ Medical Necessity Storage Basket  
☐ IV Pole (Height Adjustable and Collapsible)\*\*  
☐ Adj. Swing Away Hardware (Required for IV Pole)

\*\*Cannot be used with canopy.

\*Rear Anti-tip Tubes Required.

**Price**

**\$454**  
**\$294**  
**\$212**

CTR12	CTR14	CTR16	CTR18
C905587	C905588	C905588	C905589
C904566	C904566	C904566	C904566
C905612	C905612	C905612	C905612

**Accessories**

- ☐ Headrest Cover (Canopy)  
☐ Extended Headrest Cover (Canopy)  
☐ Upper Extremity Support Surface Tray (Clear Tray)\*\*  
☐ Adj. Removable Hardware (Required for Clear Tray)  
☐ Under Seat Storage Basket  
☐ Utility Bag  
☐ Rain Cover (Canopy Required)  
☐ Mosquito Net (Canopy Required)

**Price**

**\$351**  
**\$384**  
**\$171**  
**\$350**  
**\$176**  
**\$65**  
**\$232**  
**\$232**

CTR12	CTR14	CTR16	CTR18
C905017	C905018	C905019	C905020
C905602	C905603	C905604	C905605
C905547	C905548	C905549	C905550
C905551	C905551	C905551	C905551
C905584	C905585	C905585	C905586
C904531	C904531	C904531	C904531
C905590	C905591	C905592	C905593
C905594	C905595	C905596	C905597

\* Note: At the 30° seat angle setting with the tray, there may be interference with the patient's knees

^Cannot be used during Transit

°Height Adjustable Flip-Up Armrests and Hardware required.

The Coaster transit option includes four manufacturer installed wheelchair transportation anchors. The Coaster transit model has been crash tested and conforms to RESNA WC4:2017 section 19 and ISO 7176-19:2022 standards.

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

**Standard shipping rates apply • Expedited shipping available, please call 1-844-876-6245 or email [convoidsales.us@etac.com](mailto:convoidsales.us@etac.com)**

Prices in CAD & subject to change, ex works Torrance, California.

Effective: 1/5/2026

Convoid  
 2830 California Street  
 Torrance, CA 90503  
 1 844-US-MOBILITY (844-876-6245)  
[www.etac.us.com](http://www.etac.us.com)

OF0201-ENG-REV.00

 **Convoid**  
 by Etac