

Convoid Rodeo Tilt - WMK6CO060

Canadian Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____
 Contact: _____ ATP/RTS: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email:	
<u>Rodeo Tilt Standard Features:</u> 5° to 45° of Adjustable Tilt - Silver Powder Coated Frame - One Piece Height Adjustable Push Handle - One Piece Angle Adjustable Footplate 7.5 x 2" Front, 11.5" Rear Solid Tires - Quick Release Rolko Rear Tires - Single Foot Activated Wheel Lock - Zippered Pockets in Seat Bottom Upholstery 2-Point Pelvic Belt - Seat Depth Growth of 5" - Elevating Legrest - Tool-less Upholstery Removal <u>Standard Models:</u> ☐ Rodeo RD10 - C904340 ☐ Rodeo RD12 - C902689 ☐ Rodeo RD14 - C902691 ☐ Rodeo RD16 - C903134 <u>Price</u> \$3,798 \$3,798 \$3,798 \$3,798 All prices in CAD	

<u>Upholstery Colors</u>	<u>Price</u>
☐ Panther Black	N/C
☐ Royal Blue	N/C
☐ Sassy Purple	N/C
☐ Princess Pink	N/C
☐ Forest Green	N/C
☐ Cherry Red	N/C

RD10	RD12	RD14	RD16
C903854	C903854	C903854	C903854
C903850	C903850	C903850	C903850
C903851	C903851	C903851	C903851
C903759	C903759	C903759	C903759
C903857	C903857	C903857	C903857
C903855	C903855	C903855	C903855

<u>Frame Option</u>	<u>Price</u>
☐ Recline Option¹	\$710
☐ Height Adjustable Flip-Up Armrests	\$351
☐ Attendant Hand Brakes¹	\$363
☐ Caster Lock for 2" Wide Front Tires¹	\$144
☐ Curb Tipper (Cannot use with Rear Anti-Tip Tubes)	\$151
☐ Rear Anti-Tip Tubes (Cannot use with Curb Tipper)	\$177

RD10	RD12	RD14	RD16
C904342	C904343	C904344	C904345
C904346	C904346	C904346	C904346
		C904370	C904371
C904650	C904650	C904650	C904650
C905310	C905310	C905310	C905310
C904649	C904649	C904649	C904649

<u>Transit Required Items**</u>	<u>Price</u>
☐ Transit Anchor Brackets¹	\$308
☐ Headrest Extension	\$199
☐ H-Harness with Padded Covers (Requires 3-Point Pelvic Belt)	\$151
☐ 3-Point Pelvic Belt	\$75
Foot Positioners - Select Size Below:	
☐ Small 7" - 10"¹	\$240
☐ Medium 8" - 11"¹	\$240
☐ Large 9" - 12"¹	\$240

RD10	RD12	RD14	RD16
C904341	C903167	C903168	C903169
C903332	C902863	C900484	C900138
C903988	C903988	C903988	C903988
C904017	C904017	C904017	C904017
C903415	C903415	C903415	C903415
C903416	C903416	C903416	C903416
C903417	C903417	C903417	C903417

¹Manufacturer Installation Required.
*****Crash Tested - WC19 and ISO 7176-19 compliant.***

<u>Transit Options</u>	<u>Price</u>
☐ Q'Straint Transit Lap Belt	\$285

RD10	RD12	RD14	RD16
C903094	C903094	C903094	C903094

<u>Head Supports</u>	<u>Price</u>
☐ Padded Headwings	\$374
☐ Occi Headwings	\$351
☐ Headrest Extension	\$199

RD10	RD12	RD14	RD16
C905634	C905635	C905635	C905637
C905270	C905271	C905271	C905271
C903332	C902863	C900484	C900138

<u>Trunk Positioning</u>	<u>Price</u>
☐ H-Harness with Padded Covers (Requires 3-Point Pelvic Belt)	\$151
☐ Full Torso Support Vest	\$406
☐ Reducer Seat Insert (Reduces Seat Width 2" & Seat Depth 1")	\$589

RD10	RD12	RD14	RD16
C903988	C903988	C903988	C903988
C904364	C903261	C903260	C903259
C904211	C904211	C904110	C904110

Swing Away Laterals**Price****Right Adj Firm Lateral Trunk Sup - Contoured (Hrdw Required)**

- ☐ X Small= 2.75"x 3.75** \$181
- ☐ Small = 3"x 5" \$181
- ☐ Medium = 3.75"x 5.5" \$181
- ☐ Large = 4.75"x 6" \$181
- ☐ Right S/A Adj. Hardware for Lateral Trunk Support \$366

Left Adj Firm Lateral Trunk Sup - Contoured (Hrdw Required)

- ☐ X Small= 2.75"x 3.75** \$181
- ☐ Small = 3"x 5" \$181
- ☐ Medium = 3.75"x 5.5" \$181
- ☐ Large = 4.75"x 6" \$181
- ☐ Left S/A Adj. Hardware for Lateral Trunk Support \$366

RD10 requires X Small laterals.*Pelvic Positioning****Price**

- ☐ 3-Point Pelvic Belt \$75
- ☐ Support Cushion, General Use* \$370
- ☐ Position Cushion, Medial Thigh Support/Anti Thrust* \$717
- ☐ Align Cushion, Lateral Pelvis/Lateral Thigh Support* \$717
- ☐ Incontinence Cushion Liner (Cushion Required)* \$111

Manufacturer Installation Required.**Solid Seat Insert included with all Cushions.****Foot & Leg Positioning****Price**

- ☐ Medial Thigh Support, Abductor \$249
- ☐ Lateral Thigh Support, Adductor \$249
- ☐ Padded Elevating Legrest Cover \$144
- ☐ 2" Legrest Extension** \$161
- ☐ Calf Panel \$95
- ☐ Individual Angle Adjustable Footplates* \$530
- ☐ Full Padded Footbox (One-Piece Footplate Required) \$500

Foot Positioners - Select Size Below:

- ☐ Small 7" - 10" \$240
- ☐ Medium 8" - 11" \$240
- ☐ Large 9" - 12" \$240

Manufacturer installation required if model was made prior to June 2019.**2" Legrest Extension interferes with caster swivel when seat is in upright position.****Seat must be tilted back at least 5 degrees OR Legrests must be elevated to allow caster swivel.****Tire Options****Price**

- ☐ 12.5" Rear Solid Knobby Tires \$307
- ☐ 7.5" x 2" Front, 12.5" Rear Pneumatic Knobby Tires \$459

Medical Support Accessories**Price**

- ☐ Medical Necessity Storage Basket \$454
- ☐ IV Pole (height adjustable and collapsible)** \$294
- ☐ Adj. Swing Away Hardware (Required for IV Pole) \$212
- ☐ Oxygen Tank Bag \$269

Rear Anti-tip Tubes Required***Cannot be used with canopy****Accessories****Price**

- ☐ Headrest Cover (Canopy) \$351
- ☐ Extended Headrest Cover (Canopy) \$384
- ☐ Upper Extremity Support Surface Tray (Clear Tray)** \$171
- ☐ Adj. Removable Hardware (Required for Clear Tray) \$350
- ☐ Under Seat Storage Basket \$176
- ☐ Utility Bag \$65
- ☐ Solid Side Panels \$325
- ☐ Mesh Side Panels \$229
- ☐ Rain Cover (Canopy Required) \$232
- ☐ Rain Cover (Extended Canopy Required, RD14) \$232
- ☐ Mosquito Net (Canopy Required) \$232

Flip-Up Armrests and Hardware required.***Cannot be used during Transit.**

RD10	RD12	RD14	RD16
C904385	C904385	C904385	C904385
	C904071	C904071	C904071
	C904072	C904072	C904072
	C904073	C904073	C904073
C904362	C904362	C904362	C904362
C904465	C904465	C904465	C904465
	C904466	C904466	C904466
	C904467	C904467	C904467
	C904468	C904468	C904468
C904363	C904363	C904363	C904363

RD10	RD12	RD14	RD16
C904017	C904017	C904017	C904017
C904347	C904350	C904353	C904356
C904348	C904351	C904354	C904357
C904349	C904352	C904355	C904358
C902324	C902324	C902324	C902324

RD10	RD12	RD14	RD16
C904359	C902667	C902647	C903140
C904360	C902668	C902648	C903141
C904619	C904619	C904619	C904619
C904095	C904095	C904096	C904097
C904361	C904361	C904361	C904361
	C904957	C904958	C904959
C904953	C904954	C904955	C904956
C903415	C903415	C903415	C903415
C903416	C903416	C903416	C903416
C903417	C903417	C903417	C903417

RD10	RD12	RD14	RD16
C904117	C904117	C904117	C904117
C904118	C904118	C904118	C904118

RD10	RD12	RD14	RD16
C904476	C904476	C904477	C904477
C904566	C904566	C904566	C904566
C905612	C905612	C905612	C905612
C902870	C902870	C902870	C902870

RD10	RD12	RD14	RD16
C904366	C902679	C902659	C903152
C904728	C904729	C904730	C904731
C904367	C903405	C903401	C903402
C903369	C903369	C903369	C903369
C904529	C904529	C904530	C904530
C904531	C904531	C904531	C904531
C904388	C902672	C902652	C903145
C904389	C902680	C902681	C903146
C904939	C904939	C904939	C904940
		C904940	
C904655	C904655	C904655	C904656

The Rodeo Transit option includes four manufacturer installed black wheelchair transportation anchors. The Rodeo Transit model has been crash tested and conforms to RESNA WC19 and ISO 7176-19:2022 standards.

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Standard shipping rates apply • Expedited shipping available, please call 1-844-876-6245 or email convaidsales.us@etac.com

2830 California St., Torrance, CA 90503 • 1-844-876-6245, 310-618-0111 • Fax: 310-618-2166

Prices in CAD & subject to change, ex works Torrance, California.

Effective: 1/5/2026