

Convoid Rodeo Tilt w/out Seating

Canadian Price List / Order Form

Account #: _____ Purchase Order #: _____ Date: _____
 Contact: _____ ATP/RTS: _____ Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email:	

Rodeo Tilt w/out Seating Standard Features:

- 5° to 45° of Adjustable Tilt
- Silver Powder Coated Frame
- Adjustable One-Piece Push Handle
- Elevating Legrest / One-Piece, Angle Adjustable Footplate
- 7.5 x 2" Front, 11.5" Rear Solid Tires
- Quick Release Rolko Rear Wheels
- Single Foot Activated Wheel Lock
- Seat Depth Growth of 5"

Standard Models:

Price

- | | |
|--|----------------|
| <input type="checkbox"/> Rodeo w/out Seating RD10WOS - C904971 | \$3,568 |
| <input type="checkbox"/> Rodeo w/out Seating RD12WOS - C904972 | \$3,568 |
| <input type="checkbox"/> Rodeo w/out Seating RD14WOS - C904973 | \$3,568 |
| <input type="checkbox"/> Rodeo w/out Seating RD16WOS - C904974 | \$3,568 |

All prices in CAD

Seating Module

Price

- ☐ Depth Adjustable Solid Seat Pan¹
***Manufacturer Installation Required.**

\$232

RD10	RD12	RD14	RD16
C509549	C509550	C509551	C509552

Frame Options

Price

- ☐ Recline Option¹
☐ Height Adjustable Flip-Up Armrests
☐ Attendant Hand Brakes¹
☐ Caster Lock for 2" Wide Front Tires¹
☐ Curb Tipper (Cannot use with Rear Anti-Tip Tubes)
☐ Rear Anti-Tip Tubes (Cannot use with Curb Tipper)
***Manufacturer Installation Required.**

\$710

\$351

\$363

\$144

\$151

\$177

RD10	RD12	RD14	RD16
C904342	C904343	C904344	C904345
C904346	C904346	C904346	C904346
		C904370	C904371
C904650	C904650	C904650	C904650
C905310	C905310	C905310	C905310
C904649	C904649	C904649	C904649

Transit Required Items*

Price

- ☐ Transit Anchor Brackets¹^A
☐ 2-Point Pelvic Belt¹^A
☐ Using Other Manufacturer's Head Rest
***Manufacturer Installation Required.**
***Crash Tested - WC19 and ISO 7176-19 Compliant.**

\$308

\$69

\$0

RD10	RD12	RD14	RD16
C904341	C903167	C903168	C903169
C903931	C903931	C904559	C904559
C995681	C995681	C995681	C995681

Transit Options

Price

- ☐ Q'Straint Transit Lap Belt

\$285

RD10	RD12	RD14	RD16
C903094	C903094	C903094	C903094

Foot & Leg Positioning

Price

- ☐ Individual Angle Adjustable Footplates
☐ Padded Elevating Legrest Cover
☐ Full Padded Footbox (One-Piece Footplate Required)
☐ 2" Legrest Extension*

\$530

\$144

\$500

\$161

Foot Positioners - Select Size Below:

- | | |
|--|--------------|
| <input type="checkbox"/> Small 7" - 10" | \$240 |
| <input type="checkbox"/> Medium 8" - 11" | \$240 |
| <input type="checkbox"/> Large 9" - 12" | \$240 |

***2" Legrest Extension interferes with caster swivel when seat is in upright position.**
Seat must be tilted back at least 5 degrees OR Legrests must be elevated to allow caster swivel.

RD10	RD12	RD14	RD16
	C904957	C904958	C904959
C904619	C904619	C904619	C904619
C904953	C904954	C904955	C904956
C904095	C904095	C904096	C904097
C903415	C903415	C903415	C903415
C903416	C903416	C903416	C903416
C903417	C903417	C903417	C903417

Tire Options

- ☐ 12.5" Rear Solid Knobby Tires
☐ 7.5" x 2" Front, 12.5" Rear Pneumatic Knobby Tires

\$307

\$459

RD10	RD12	RD14	RD16
C904117	C904117	C904117	C904117
C904118	C904118	C904118	C904118

Medical Support Accessories

- ☐ Medical Necessity Storage Basket **\$454**
- ☐ IV Pole (Height Adjustable and Collapsible)^{2^^} **\$294**
- ☐ Adj. Swing Away Hardware (Required for IV Pole) **\$212**
- ☐ Oxygen Tank Bag **\$269**
- ²**Rear Anti-tip Tubes Required.**
- ^{^^}**Cannot be used with canopy**

Accessories

- ☐ Headrest Cover (Canopy) **\$351**
- ☐ Extended Headrest Cover (Canopy) **\$384**
- Canopy Colors (Choose One Below)
- ☐ Panther Black **N/C**
- ☐ Royal Blue **N/C**
- ☐ Sassy Purple **N/C**
- ☐ Princess Pink **N/C**
- ☐ Forest Green **N/C**
- ☐ Cherry Red **N/C**
- ☐ Upper Extremity Support Surface Tray (Clear Tray)^{^^} **\$171**
- ☐ Adj. Removable Hardware (Required for Clear Tray) **\$350**
- ☐ Solid Side Panels **\$325**
- ☐ Mesh Side Panels **\$229**
- ☐ Under Seat Storage Basket **\$176**
- ☐ Utility Bag **\$65**
- ☐ Rain Cover (Canopy Required) **\$232**
- ☐ Rain Cover (Extended Canopy Required, RD14) **\$232**
- ☐ Mosquito Net (Canopy Required) **\$232**
- [^]**Flip-Up Armrests and Hardware required.**
- ^{*}**Cannot be used during Transit.**

RD10	RD12	RD14	RD16
C904476	C904476	C904477	C904477
C904566	C904566	C904566	C904566
C905612	C905612	C905612	C905612
C902870	C902870	C902870	C902870

RD10	RD12	RD14	RD16
C904366	C902679	C902659	C903152
C904728	C904729	C904730	C904731
C903854	C903854	C903854	C903854
C903850	C903850	C903850	C903850
C903851	C903851	C903851	C903851
C903759	C903759	C903759	C903759
C903857	C903857	C903857	C903857
C903855	C903855	C903855	C903855
C904367	C903405	C903401	C903402
C903369	C903369	C903369	C903369
C904388	C902672	C902652	C903145
C904389	C902680	C902681	C903146
C904529	C904529	C904530	C904530
C904531	C904531	C904531	C904531
C904939	C904939	C904939	C904940
C904655	C904655	C904655	C904656

The Rodeo Transit option includes four manufacturer installed black wheelchair transportation anchors. The Rodeo Transit model has been crash tested and conforms to RESNA WC19 and ISO 7176-19:2022 standards.

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Standard shipping rates apply • Expedited shipping available, please call 1-844-876-6245 or email convaidsales.us@etac.com

Prices in CAD & subject to change, ex works Torrance, California.

Effective: 1/5/2026

Convaid
2830 California Street
Torrance, CA 90503
1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

OF1501-ENG-REV.00

by Etac